

# Emergency & OSHA Supplies

## Allegiance



### Safe Care Isolation Gowns

(Techstyles) Basic coverage limited use isolation gowns are available in either standard-weight non-woven fabric or economical lighter-weight fabric. These gowns have a waist belt. Yellow color. One size. Gown 46969-075 features 46.5" length and a tape-tab neck closure. Glued construction

BGH46969-075  
Lightweight gown, 100/cs



### Safe Care Polyethylene Fluid-Barrier Gowns

(Techstyles) These single-use halfback polyethylene gowns slip easily over the head. They are available in three styles: with thumb-loop wrist, with tape-tab wrist. The thumb loop provides secure holding inside gloves. While the tape tab allows for flexibility in sizing and gloving. All styles may be removed easily without violating good technique

BGH46969-098  
Yellow Gown with Tape-Tab wrist, large, 100/cs

## SAGE Products, Inc



### P2 Disposable Gowns

P2 Disposable Gowns offer comfort and polyethylene protection. One-piece torso/sleeve designs reduce the number of vulnerable seam areas. Open back panels offer cool comfort, and yoke-style necks eliminate cumbersome ties for easy donning and removal. Stirrup cuffs provide added protection

to the vulnerable wrist areas, and make donning gloves easy.

SAG8575 Medium, 20/pk, 5pk/cs  
SAG8576 Universal, 15/pk, 5pk/cs



### Eye Safe Faucet-Mounted Eyewashes

(Water Saver Faucet) Convert any faucet into an emergency eyewash station without special plumbing or interfering with normal faucet operation. Threaded inlet attaches to any standard laboratory faucet. Pull knob to active. Water pressure holds eyewash in operation. Angle of aerated outlets is adjustable to permit full coverage and avoid splashing the counter or floor, Intended as supplement to, but not replacement for, standard plumbed-in eyewash equipment. Eyewash S1396-10 has outlets 3" part and angled forward. Eyewash S1396-11 has outlets mounted on bends ending 5" apart and angled in toward face. Use only cold water when eyewash is in operation

BGHS1396-10

Eye Safe eyewash, 3”

## BANYAN



### **Banyan Stat Kit 600**

The Stat Kit 600 is your best choice if you already have oxygen, resuscitation and suction equipment in your office. You have 25 different drugs, plus the equipment to start treatment on virtually any emergency associated with a medical practice. Quarterly update notification reminds you to inspect the kit for expired items. Weights 16lbs  
BYN600 Stat Kit 600



### **Banyan Stat Kit 800**

Same contents as the Stat Kit 600, plus an oxygen system, disposable Ambu resuscitator and manual suction unit, the oxygen system is preset at 8LPM for approximately 30 minutes. Weight 28lbs

BYN800 Stat Kit 800

### **Banyan Stat Kit Plastic Seal**

BYN421 10/pk

### **Banyan Stat Kit Auto-Update Program**

Annual replacements for all items in the Stat Kit Which have expiration dates. This insures that the Stat Kit is fully equipped and ready for any medical emergency requiring your attention

## Allegiance



### **Resuscitation Devices**

Reliable and affordable, our latex-free disposable resuscitation devices answer your contamination concerns in today's healthcare environment. Our device provides low-profile textured bag, Preattached PEEP, manometer port, double swivel elbow and removable accumulator.

### **Manual Resuscitators without Peep Valve**

BGH2K8008 Pediatric, with mask, 40” oxygen reservoir tubing, 6/cs

BGH2K8004 Adult, with mask, oxygen reservoir bag, 6/cs  
BGH2K8005 Adult, with mask, 40” oxygen reservoir tubing, 6/cs

## AMBU



### **Resuscitator Mask**

The Ambu SPUR is fully disposable and features a see-through patient valve allowing for visual check of operation. The Ambu SPUR features a textured surface and comes with a unique support which facilitates a steady grip and uniform tidal volumes. The Ambu SPUR is 100% latex free

AMB420-011-000  
Adult Mask

## Hudson RCI



### **Adult Resuscitators**

The Hudson RCI LIFESAVER disposable manual resuscitator combines quality, features and performance.

Designed to meet ISO and ASTM resuscitation bag standards, Available with optional 60cm pressure relief pop-off valve 100% latex-free. Individually packaged in convenient, reseal able hanging bag, 6/cs

HUD5370 Adult Mask

**Johnson Johnson Medical**

**BARRIER Protective Goggles**

These goggles can be worn comfortably over eyeglasses. They features solid side shields per OSHA guidelines

JJ1701 1pr/bx, 30bx/cs

**BARRIER Protective Glasses**

Lightweight, comfortable and designed to give you maximum peripheral vision without distortion they features solid sides shields per OSHA guidelines

JJ1072 10 pr/bx, 3bx/cs

**MADA**



**Disposable Adult Manual Resuscitator**

Disposable manual resuscitator with oxygen bag reservoir, Designed for single patient use. Anti-locking one-way valve prevents lock-up under high flow conditions. 6/cs

MAD2200 Adult



**Disposable Manual Resuscitator**

Disposable manual resuscitator with flexible tube reservoir designed for single patient use. Anti-locking one-way valve prevents lock-up under high flow conditions. 6/cs

MAD2201 Adult  
MAD2203 Child

**Olympic Medical**



**Papoose Board**

Temporary fast, safe patient restraint during examination or treatment Canvas flap set allows any part of patient to be selectively exposed for treatment while patient is securely immobilized. Includes flap set, head strap, and two arm straps

OLY50500 Papoose Board, 2-6 years

**Quality America**



**OSHA Safety Program Manual**

Protect your staff and comply with OSHA requirements with our fill-in-the-blanks Exposure Control Plan. Two-inch binder contains everything needed for compliance with OSHA's Bloodborne Pathogens, Hazard Communication, and TB regulations. Includes required forms, copies of actual OSHA regulations, eyewash sign, OSHA posters, compliance checklist, sample test for training sessions, two-inch MSDS

binder, and much more.  
Includes PC disk for  
Microsoft word 6

QAI3001 OSHA  
Safety Program Manual

### **OSHA Yearly Retraining Video Learning System**

Have lunch with the folks  
at Dr. Vascular's office  
who just found out that  
the OSHA inspector is on  
his way. Luckily, a  
quality America  
consultant arrives to show  
them how to pass an  
OSHA inspection with  
flying colors. Comical  
vignettes throughout this  
interactive video prompt  
your staff to identify "What's  
Wrong with this  
Picture"? Painless full  
OSHA's yearly retraining  
requirements. Covers  
Bloodborne Pathogens,  
Hazard Communication  
and TB regulations

QAI5000 OSHA  
yearly retraining Video  
Learning System

QAI5001 OSHA  
Video Learning System  
workbook

### **OSHA Watch Newsletter**

Eight page, newsletter  
keeps you informed of  
current OSHA regulations  
and most recent changes.  
Latest breaking news and

practical tips to ensure  
medical practice  
compliance Includes  
helpful Q & A section.  
Six issues/year, 1year  
subscription

QAI4001 OSHA  
Watch Newsletter

### **Sharps Injury reduction Video Program**

Includes: OSHA Manual,  
Annual Training Video "I  
can believe they do  
that!", facilitator's guide,  
workbook, eyewash sign,  
MSDS binder, OSHA  
Watch newsletter (1 yr) &  
new sharps Injury  
Reduction Program  
"Point of No Return"  
Video, w/ guide

QAI2002 Sharps  
Injury Reduction Video  
Program

### **Sharp Plus – Complete Compliance Program**

Includes: 15 Minute  
"point of No Return!"  
video, w/ guide. The  
guide includes detailed  
instructions for  
performing sharps  
evaluation and forms to  
document these  
evaluations. Also  
includes Sharps injury log

QAI2001 Sharps  
Plus – Complete  
Compliance Program

### **Complete Compliance Package**

Simply the best Get in  
compliance and stay in  
compliance with our  
Manual, Video and  
Newsletter

QAI3003 Complete  
Compliance Package

## **BioMed**



### **NICSx Compliance Assistance Program**

NICSx web-based  
software publishes your  
Sharps Exposure Control  
Plan. Trains and certifies  
caregivers on plan,  
records safety product  
evaluations, and logs any  
needle stick incidents.  
Review data preventing  
future needle stick  
incidents

N005 NICSx for 1 to 5  
Employees

N010 NICSx for 6 to 10  
Employees

N10P NICSx for  
Additional 10 Employees

N100 NICSx for up to  
100 Employees

N025 NICSx for up to  
25 Employees

N050 NICSx for up to  
50 Employees

NMAN NICSx paper  
Manual

## Quality America

### OSHA Compliance Eyewash Station Sign

QAI3005 Eyewash  
Station Sign

## RUSCH



### Berman Airway Oral Tips

RUS122005 Lg.  
100mm, 10/bx  
RUS122004 Med.  
90mm, 10/bx  
RUS122003 Small,  
80mm, 10/bx

## Graham-Field

### Berman-Type Oral Airway Kits

Seven, piece emergency  
Kit contains one of each  
size Guedel-Type oral air  
ways 110mm. Disposable

GF13787 KT Snap-lock  
plastic box, 1/kit

### Berman-Type oral Airways

Smooth, Flexible,  
Medical-grade white  
polypropylene with,  
reinforced block to  
prevent collapse.  
Recommended for single  
use only. Disposable.

3786-12#2 60mm  
Child, 12/bx  
3786-12#1 50mm  
Infant, 12/bx  
3786-12#5 100mm  
large adult, 12/bx  
3786-24#5 100mm  
large adult, 24/bx  
3786-12#4 90mm  
medium adult, 12/bx  
3786-24#4 90mm  
medium adult, 24bx  
3786-12#3 80mm  
small adult, 12/bx

### Stocked First Aid Kits

25 Count Kits contain a  
balanced assortment of  
essential first aid supplies  
adequate for the number  
of people indicated for  
each kit. Construction of  
heavy gauge enameled  
steel. Weather, rust and  
dust resistant can be wall  
hung. Cover acts as  
additional shelf when  
open. Contains: 100 -3/4"  
x 3" plastic Strips, 1-1/2  
oz. Absorbent cotton, 1-  
2" Bandage, 4 oval Eye  
pads, 4 Antiseptic wipes,  
1-1/2" x 5 yds. Adhesive  
tape, 1-1" Bandage, 10-  
2"x3" Non-Adherent  
pads, 1 Triangular  
bandage, 1-1-1/2 oz. tube  
first aid cream, 2 -  
Ammonia inhalants, 1  
cold pack, 1 Scissors, 1  
First Aid Booklet &  
Content Card, Serves up  
to 25 people

GF1799-25 10-7/8" x  
7-1/2" x 2-5/8" 1ea

## Welch Allyn



### Disposable Laryngoscopes

These battery powered  
KleenSpec disposable  
laryngoscopes are  
designed for single use.  
The reusable light source  
inserts in the handle.  
After use, remove the  
light source and discard  
the blade and handle. A  
clean, new laryngoscope  
is always ready for instant  
use

WA56903 Mac 3-1/2,  
10/bx  
WA56901 Mac 1-1/2,  
10/bx  
WA60500 Illuminator  
for Disposable  
Laryngoscope

## Sharps Containers

## Bemis



### Wall Safe Sharps Container

Ideal for patient rooms, Bemis wall safe mounts unobtrusively at any point-of-use location. The rotating cylinder top permits one-handed disposal of sharps, dropping needles horizontally to maximize capacity.

BEM150020 5 qt

### BioMed



### Sharpx Needle Destruction Unit

Sharpx needle destruction unit saves \$25,000 per 100,000 needles, assisting in compliance with the new, revised OSHA regulations. Incinerate existing needles eliminating blood borne pathogens instead of purchasing expensive safety needles

8001 Sharpx Needle Destruction Unit

0024 Sharpx One Year service contract

0012 Sharpx Two Year service Contract

### Sharpx Needle Destruction Accessory Pack

Includes 10 wipes and 10 cups for use with the Sharpx Needle Destruction Unit

5009 Sharpx Accessory Pack

### Futura

The New Thinking Safety

### Pro-Tec Sharps Containers

For safe disposal of sharps, Pro-Tec containers are constructed of rigid, puncture-resistant polypropylene, supplied pre-assembled and ready to use with grip handles and locking lids for safe transport.

Containers are available for both biohazard us (red) and chemotherapy (yellow) materials

FTAC1GRED 1 gallon (150-180 syringe capacity), Red, 18/cs

FTAC2GRED 2 gallon (260-280 syringe capacity), Red, 12/cs

FTAC3GRED 3 gallon (460-500 syringe capacity), Red, 8/cs

### MAXXIM Medical



### Maxxim Sharps Collectors

Red sharps collector with clear locking lid

MXX8702 1 quart, 72/cs

MXX8705 16 quart, 10/cs

MXX8703 4 quart, 24/cs

MXX8704 8 quart, 24/cs

### Post Medical



### Automatic Needle Disconnect Sharps Collector

The Automatic Needle Disconnect (A.N.D.) phlebotomy container from Post medical is the ideal system for safe disposal of blood collection needles, lancets and butterflies. Molded of durable, puncture-resistant plastic, this unit's compact size fits easily into phlebotomy trays. Separate A.N.D. units are designed to work with BD

Vacutainer, Sherwood Monoject, Terumo, Nissho, and Greiner blood collection needles. Simply push holder downward with ONE HAND – the blood collection needle is automatically unwound and deposited safely into the container. 24 per case

PMIBE-100 200  
Needles Capacity

## BIOHAZARD BAGS

UNISOURCE

**Biohazard Bags**  
Red, 1.25 mil bags

UNJLHCP2423HPR 24”  
x 23”, 10 Gal, 200/cs  
UNJLHCP3043HPR 30”  
x 43”, 33 Gal, 60/cs  
UNJLHCP3339HPR 33”  
x 39”, 33 Gal, 60/cs

## SPILL KITS

DOW



### Bio Hazard Spill Kit with DISPATCH

Single use spill, kit with 8oz bottle disinfectant, 6kits/case. Pre-assembled kits with choice of

disinfectant (bleach-non-bleach)

DOWSK960 Spill Kit

**Bio Hazard Spill Kit with PRECISE QTB**  
Kit with precise QTB, non-bleach based disinfectant or kit with Dispatch, bleach based disinfectant. Convenient kit for containing, removing and disposing of spilled body fluids. Meet OSHA & CDC standards for blood borne pathogens

DOWSK008 Spill Kit

## Futura

The New Thinking Safety



### BASK Biohazard Spill Kits

The Bask safety systems, provides everything needed for safe, easy and effective clean up, handling, disposal and disinfection of biohazard spills. Available in multiple or single use kits. Components also sold separately

FTA8300 Multi-use  
kit, 4 kt/cs

FTAH-2500 Single-use  
Kit, 36 kt/cs

## Premium Plastics

### Emergency Response Kit

These convenient emergency response kits are great for quick biohazard us spill cleanup. Select from four kits featuring effective Red Z solidifiers specially designed for encapsulating spills

### Emergency Response Kit

Kit includes: two Vinyl gloves, one .75 oz. Red Z Packet, one shovel, one plastic bag (red), Bag Ties, one Sanizide Solution/Dry Towel, one Towel, one Antiseptic Hand Wipe, Instructions, and one Scraper

PP2038 Up to .47  
gallons, 72/cs

### Deluxe Emergency Response Kit

Deluxe Kit includes: two vinyl gloves, one 2 oz. Red Z Packet, One shovel, one plastic bag (red), one plastic (white), two bag ties, one Germicidal towelette, one towel, one Antiseptic Hand wipe, Instructions, one Splash Goggles with Mask, One disposable Gown and one scraper

PP2037 Up to 1.25  
Gallons, 72/cs

### Red-Z Spill Solidifier

Premium plastics offers, the most complete line of no biodegradable solidifiers for virtually every situation involving blood or other potentially infectious fluids. Red Z has an extremely powerful polyacrylate base for optional solidification

PP2030 8oz., 6/cs

## BATTERIES

### Duracell



#### Duracell Procell Professional Batteries

Procell alkaline batteries are manufactured to the same performance and quality specifications as DURACELL COPPER TOP

PC1300 Alkaline  
"D" cell, 12/ct, 72/cs

PC1400 Alkaline  
"C" cell 12/ct, 72/cs

PC1500 Alkaline  
"AA" cell 4/box 12/ct,  
72/cs

PC1604 Alkaline,  
"9 volt" cell 12/ct, 72/cs

PC2400 Alkaline  
"AAA" cell 4/box, 12/ct,  
72/cs

PC7K67 Alkaline,  
6.0v "Flat J" cell 6/ct,  
24/cs



DL2032B Lithium,  
3.0v .787" diameter,  
24/cs

DL2450B Lithium,  
3.0v .965" diameter,  
24/cs

#### Duracell Copper Top Batteries

MN1604 Alkaline,  
"9 volt" cell 12/ct, 72/cs

MN1500 Alkaline,  
"AA" cell 4/box 24/ct,  
72/cs

MN2400 Alkaline,  
"AAA" cell 12/ct, 72/cs

### Welch Allyn



#### Welch Allyn Replacement Batteries

Replacement batteries that will keep your Welch Allyn medical products working for years to come

WA72200 NiCad  
Rechargeable Battery  
(Black)

WA72300 Ni-Cad  
Rechargeable Battery  
(Orange) for 71000-A  
WA72600 2.5v  
Rechargeable Battery for  
pocket Scope  
WA72700 7.2v  
Rechargeable Battery for  
Original Micro Tympanometer  
WA72000 2.5v  
Nickel-Cadmium  
Rechargeable Battery  
(Red)  
WA72100 2.5v  
Rechargeable Battery  
(Green)  
WA72500 7.2v  
Rechargeable Battery (2)