

Furniture & Equipment

Bone Density Testing

HOLOGIC



QDR 4500 Elite

The QDR 4500 Elite combines the superiority of fan-beam technology with a unique, value-added package of clinical, patient reporting and support software. It delivers high-performance imaging and bone density measurement with unsurpassed precision and low patient dose. The Elite protects your investment from obsolescence with a modular design that easily accommodates upgrades and enhancements as the technology advances and your practice grows. In just 30 seconds or less, Elite's fan beam x-ray completes hip, PA spine and forearm scans with uncompromised precision and image quality and low patient dose. That's nearly five times faster

than the fastest pencil-beam system takes to produce similar, clinically meaningful scans, Faster scanning lets you accommodate more patients in your busy schedule and contributes to your potential for practice growth.

Equipped with a powerful computer system, large monitor and fast color printer, the Elite package contains a JAZ drive for data storage, CD-ROM, and a modem to connect HOLOGIC customer support to you for assistance with equipment operation and clinical procedures. Two reporting software programs – QDR One Page and optional Osteo Report – provide concise, attractive patient reports for referring physicians and to facilitate reimbursement

HLG49796 QDR
4500C Elite



QDR 4000 X-Ray Bone Densitometer

Introduced in 1987, Hologic QDR bone densitometers quickly became the standard for assessing bone mineral status. The QDR 4000

combines the reliability and economy of classic DXA technology with a unique package of value-added peripherals and software to deliver precise bone density measurements of the hip, spine and forearm. The QDR 4000 features high-quality imaging, low patient dose, unsurpassed precision, and a compact, low-profile design for easy patient positioning. Like all QDR instruments, the QDR 4000 utilizes the NHANES reference database and is the only Bone Densitometer in its class to meet national standardization criteria. Equipment with powerful computer system, large monitor and fast color printer, the QDR 4000 package contains a JAZ drive for data storage, CD-ROM, and a modem to connect you to HOLOGIC customer support for help with scan analysis. Two reporting software programs – QDR One Page and the optional Osteo Report – provide concise, attractive patient reports for referring physicians and to facilitate reimbursement

HLG45712 QDR4000



SAHARA Clinical Bone Sonometer

In response to a growing array of effective drug therapies now available for the treatment of osteoporosis, the evaluation of bone status has become an essential component of women's health. Ultrasound bone densitometry a safe, non-ionizing modality – provides precise quantitative assessment of skeletal status, information particularly useful for identifying patients at risk of developing osteoporosis and for assessing their risk of fracture. Now you can add bone testing capabilities to your practice or clinic with the waterless, portable, easy-to-use SAHARA Clinical Bone Sonometer. In less than a minute and with the push of a button, you can determine a patient's bone mineral density (BMD), based on an ultrasound measurement of the calcaneus (heel bone) – the preferred peripheral site proven in numerous prospective studies to predict fracture risk. With Sahara, there's finally a dry ultrasound modality for bone

assessment that is simple, convenient and practical enough for the office-based physician

HLGSAHARA-02
SARARA



Delphi Clinical Bone Densitometer

The Delphi raises the standard in the practice of clinical bone densitometry by integrating bone mineral density (BMD) with instant Vertebral Assessment (IVA). Delphi's IVA, possible only with fan beam technology, generates a high resolution image of the entire lateral spine in just 10 seconds, enabling physicians to visually assess vertebral status for a accurate determination of fracture risk than just BMD alone

HLGDELPHI Delphi

Bone Density Testing Accessories

HLG295-0480 Bone Density Couplant gel, 60 ml tube, 36/cs
HLGPAPER-SAH-4SA Exam Paper 500/bx
HLG010-0919 Bulk Exam Paper, 500/bx

HLG080-0630 Paper, Patient Report Form Pad, 100 sheets/bx, 10bx/cs
HLG295-0552 Thermal paper Roll, 5/bx

CARTS

Lakeside MFG. INC.



Modular Utility Cart

Made of high-density polyethylene, the shelves resist stains and odors. Attractive beige color blends harmoniously into any environment. Full-width handle for easy cart handling and maneuvering, the three 16" x 24" shelves have upright edges on three sides to keep contents in place. Open side for easy loading/unloading. 4" swivel casters. Carrying capacity is 200 pounds. Overall size is 16-3/4" x 32-1/2" x 34-1/2"

LAK2000 Utility Cart

Chairs & Stools

The brewer Company



Step Stool

Ideal for exam room and x-ray areas, Features skid resistant rubber mat and chrome-plated 18G steel construction

BRW11200-1 Step Stool, 1/bx
BRW11200-2 Step Stool, 2/bx



Step Stool with Handrail

Same as BRW11200-1 but with 50014 handrail. Handrail is 7/8" chrome-plated steel, 35" height from top of rail to floor

BRW11220 Step Stool with Handrail, 1/bx

MIDMARK



Midmark 193 Exam Stool

The Ritter 193 features an amply padded, wide seat cushion and five-caster, black composite base with fixed height rotation. The thick, but firm, foam padding enhances comfort, while the dual hooded silent casters allow the physician to sit on the edge of the seat...roll around the room...and never tip. Available in Sand Grey (shown)

RIT193001SD 193 Exam Stool

Dopplers

Huntleigh healthcare

Fetal Dopplex 1 Plus

The Fetal Dopplex 1+ fetal heart rate display Doppler offers all the traditional features of an obstetric Doppler, combined with an advanced heart rate display capability. Rugged and reliable, the Fetal Dopplex 1+ is ideal for use in almost every setting. Features include: Real time fetal heart rate display with autocorrelation signal processing, permanently attached 2MHz probe with integral probe storage and low battery indicator

HUNFD1+ Fetal Dopplex 1 plus

NICOLET Vascular



Pocket-Dop II – Vascular

Pocket-Dop II – Known worldwide as the leading Doppler. Four interchangeable probes, two obstetrical and two vascular, and a built-in speaker are features that makes this device versatile for many applications. Other features include belt clip, sturdy, protective, carrying case, recharger, and a five-year parts warranty. Vascular applications include helping to detect arterial disease and venous disease, assessing blood pressures, checking for impotence and locating vessels for IV s and any other application where detection and monitoring of blood flow is important

IMXP300 Pocket-Dop II

IMXC620 Pocket-Dop II Battery

Hewlett Packard

Doppler Monitoring Charts

HPM1910A Fetal
Orange Grid, 150mm x
49' FHR/UA, Scales: 30-
240/0-100, 40/cs

EENT INSTRUMENTS

Welch Allyn



AM 232 Manual audiometer with AC Power

Allows quantifying hearing loss related to otitis media, ototoxic drugs, presbycusis, or any of the many factors that affect hearing. Obtains the precise threshold of patient's hearing to better identify specific patterns of hearing loss. Enables testing in a full range of frequencies and intensity levels with steady, frequency-modulated, and pulsed stimulus modes. Test frequency values are easily visible at the dial. Displays HL levels in large, easy-to-see numbers on the LCD.

Silent-functioning front panel controls prevent patient cueing; removable cover functions as a screen between the operator and subject. Sleek, ergonomic design for easy room-to-room portability; takes up minimal counter space. Microprocessor based technology assures years of dependable service

WA23200 AM232
Manual Audio meter

Electrosurgical

Welch Allyn



Audio Path Screener with printer/charger Module

EOAE technology stimulates and detects sounds created by a biomechanical process originating from the outer hair cells within the cochlea. This eliminates the need for active patient response required in traditional audiology screening tests. Ensures early detection of hearing loss in infants, children, and disabled persons. Acquires test results in as little as 10 seconds. Portable and light, weight

for performing the test virtually anywhere in your office. Print and record test results from the screener's built-in printer in 8 seconds

WA29200 Audio
Path Screener



Audio Scope 3

Screen your patient's quickly, accurately and with more flexibility, Has 20, 250 and 40 dB HL levels in one unit. Lightweight, handheld design allows complete portability. Viewing capability eliminates inaccurate results due to occlusion by wax or collapsed ear canal. Rechargeable with charging transformer or convenient charging stand

WA92680 With set of
Audio Specs 110-130v
charging transformer
stand, recording forms

WA23300 Audio
Scope 3 with Audio
Specs, no transformer
WA71123 Charger
Stand

WA55230 Recording
Forms, 10/pkgs/bx



Micro Tym 2

Diagnosing chronic otitis media is often a difficult task. But the delivery of sound and air pressure to the ear, precisely measured by advanced microelectronic components, brings otitis media into sharp focus. Welch Allyn's Micro Tym 2 tympanometry instrument produces an accurate, objective measurement of middle ear status. It also helps detect perforated tympanic membrane, patent tympanostomy tube, ossicular disruption, tympanosclerosis, and cholesteatoma

WA93690 Micro Tym 2 set with Charging Stand and printer
WA23630 Set of Micro Tym Tips (four color-coded sizes)



Sure Sight Vision Screener

Finally, a child-friendly vision test that's automatic, objective, and takes just seconds. Ensures reliable early detection of refractive error, the primary vision disorder in children. Five-second, automatic test addresses compliance problems and makes testing easier. Lights and sounds engage the patient's attention. Automatically measures and indicates abnormal readings – no provider interpretation is required. Operates at a non-threatening 14-inch (35cm) distance. Requires minimal cooperation, making it ideal for use on young children, the disabled, and when there is a language barrier. Screens for common vision problems, including near-and farsightedness (myopia/hyperopia), astigmatism (asymmetrical focus), and anisometropia (unequal power between eyes)

WA14000 Sure Sight Vision Screener

EXAM LIGHTS

Ritter
By Midmark



Ritter 152 Halogen Exam Light

The Ritter 152 Exam uses a 35-watt quartz halogen bulb, which provides 1,500 foot-candles at 24" for an optimum amount of quality light. The 3,000K, color temperature provides cool white light for optimum color recognition, yet, remains cool to the touch. A unique heat sink effectively removes heat from the bulb and its socket to increase the life of the bulb for minimal downtime. A small head section allows for high quality lighting, that's not burdensome. With its flexible arm and gooseneck, the Ritter 152 is easy to adjust, and won't drift while you're performing your procedures, and is covered in durable plastic for easy cleaning. Cord at the bottom of the light prevents tangling

RIT152001 152 Exam Light

BURTON



Genie

The new Genie Operating Room Surgery Light is manufactured in the U.S. by Burton Medical and offers over 6500 foot candles at one meter, 4400oK color temperature, shadow reduction, non-drift positioning, focus ability and, best of all, an elegant, simple design. Plus, the Genie features a unique, easy bulb replacement and storage system for maximum reliability and service. (Available in floor stand, single or double ceiling models and combo with Genesis)

BUT0300122 Single Ceiling
BUT0300132 Double Ceiling
BUT0600192 Genie/Genesis Ceiling



Outpatient II

The Outpatient II, the undisputed industry leader in the minor surgical category, offers maximum intensity, drift-free arms, shadow-free operation and focusable flood-type light. With 7000 foot, candles at 24" and 3600oK, the Outpatient II is the perfect low maintenance, high output light for nearly every application in the hospital, surgical suite, or physician's office. (Available in single or double ceiling, single or double on Fastrac, wall mount and floor stand models)

BUT0224112 Floor stand
BUT0224122 Single Ceiling
BUT0224132 Double Ceiling
BUT0224142 Wall
BUT0204470 Single Outpatient on Fastrac
BUT0204472 Double Outpatient on Fastrac



Cool Spot II

The new, aesthetically appealing Cool Spot II offers incredibly pure,

white light along with other impressive features such as guaranteed no-drift arm, intensity control, and variable spot size. No light in its price category can match the Cool Spot it's range of 10,000 foot candles at 24" (and 4000 foot candles at one meter) as well as color temperature of up to 4700oK. (Available in single or double ceiling, single or double on Fastrac, wall mount and floor stand models)

BUT0244312 Floor stand
BUT0244322 Single Ceiling
BUT0244332 Double Ceiling
BUT0244342 Wall
BUT0204370 Single Cool Spot II in Fastrac
BUT0204372 Double Cool Spot II on Fastrac



Disposable Light Covers

Saves time by eliminating the need to autoclave handles and allows compliance with OSHA's standard for sterile procedures

BUT0008100PK Set of 25

Welch Allyn



Exam Light III

Fiber optic technology provides homogenous illumination with no heat. Welch Allyn's Exam light III takes flexible illumination to a new level. Exam Light III is versatile enough to accommodate your needs. It mounts easily to fit on a wall, stand, or examination table. The Exam Light III delivers bright, white, fiber optic illumination for better visibility and more accurate diagnoses

WA48700 Exam
Light III

Accessories

WA04200 Halogen
Lamp for Exam Light III
WA48859 Universal
Table Mount
WA48850 Mobile
Stand
WA52640 Disposable
Sheaths for Exam Light,
25/bx

Welch Allyn



LS-135 Examination Light

Long-lasting, 35-watt halogen lamp produces significantly whiter, brighter light than ordinary gooseneck units. Dual intensity control switch allows you to instantly adjust light output to match your exact requirements. Molded lexan light head does not conduct heat keeping the luminaire cool to the touch. Lamp has been offset at a 12 degree angle relative to the luminaire housing, allowing an unobstructed view of the field of interest. Compact base frees up valuable exam room floor space. Smooth surfaces are easy to clean. Encased, flexible 24" arm enables accurate one hand positioning you don't get with conventional gooseneck units. Multiple installation options – including floor stand, table mount, wall mount, and caster base – enable you to customize your setup to best meet your needs

WA44300 LS-135
Examination Light

Measuring Devices

Seca

Stadiometer Accu-Hite

Measuring Range 24 to 83" form its built-in leveling bubble to its handsome design. Accu-Hite has long been the Stadiometer of choice. The headpiece comes with a lock knob for measurement accuracy and safety

SEC216 Stadiometer



Growth Chart

Growth Charts for boys (blue) and for girls (pink) are the fast, convenient way to track length/height and weight percentiles. Plot new born up to 36 months on one side and 2 to 18 years on the other side

SEC405B Boy,
100/pad
SEC405G Girl,
100/pad

Pulse Oximeters

Spacelabs Medical Burdick Products



Pulse Oximeter

Spacelabs Burdick 100 pulse Oximeter offers accurate non-invasive blood-oxygen saturation analysis and heart rate monitoring in a portable, hand-held unit. Accurate SPO2, Pulse rate, and pulse strength measurements, for neonate to adult

BUR92521-1A Pulse Oximeter

Nellcor



NPB-40 handheld Pulse Oximeter

The NPB-40 Handheld Pulse Oximeter is compact and a affordable. It is simple to operate, provides spot-check

monitoring of SPO2 and pulse rate, and is compact enough to fit in your pocket. Uses convenient AA-sized batteries. Monitor includes: 4-AA Alkaline Batteries, one DuraSensor Finger Clip Sensor (DS-100A), and Operator's Manual

PURNPB-40 Handheld Pulse Oximeter



Oxisensor II Adhesive Sensors

Nellcor Oxisensor II Adhesive Sensors maintain strong signals from weak pulses right from the start and provide a comfortable, secure fit for long-term monitoring

- PURD25 Oxisensor II Adult Sensor (18" cable)
- PURD20 Oxisensor II Pediatric Sensor
- PURI20 Oxisensor II Infant sensor
- PURN25 Oxisensor II Neonatal/Adult sensor

Scales



Physician Scale

350 lb capacity in 4 oz. increments in easy to read die cast beam. Platform is 10.5" x 14.5" with removable poly alloy waterfall cover. Dual reading height rod reads 30x78 and 75cm x 78cm. Steel base and torque lever system. Electro powder baked epoxy finish available in white, champagne and gray

DET439 Physician Scale

Scale Accessory

Counterweight, which increase scale capacity by 100lbs

DET3P7030X 100lb. Counterweight



Pediatric Scale

130lb capacity in 1 oz. increments, Washable styrene tray measures 22"x14-3/4"x4-1/2". Dia-cast, easy to read beam.

Heavy-duty rust resistant understructure 33lbs

DET450 Pediatric Scale

HEINE



Electronic Infant Scale

The professional Electronic Infant Scale designed to meet the highest demands in accuracy and reliability. Durable and Light weight, portable and long lasting. Detachable Tray is ideal for weighing of toddlers (up to 44lbs) and cleaning. Battery Operated (9V) or optional mains power supply. Heine's exclusive "Quattro Tronic" technology ensures precise readings, no matter where the patient is positioned on the weighing platform. "Weight Lock" function provides fast, accurate readings even if patient is moving

HEIS-831.61 Electronic Infant Scale

SECA



Seca 700 Scale

The new simplicity: the seca 700 mechanical weighing station saves time and unnecessary effort when weighing and measuring. Body height and weight are measured precisely in one go. Its stability and sturdiness as well as a big, non-slip platform are in an added plus for nursing staff and patients. Measuring rod included

SEC700 700 Scale

Sigmoidoscopes

Welch Allyn



KleenSpec Disposable Sigmoidoscope with Obturator

Individually wrapped for safety and convenience. Highly visible graduations, Disposable Sigmoidoscopes are individually wrapped for safety and convenience with highly visible graduations.

WA53130 KleenSpec Disposable Sigmoidoscope with Obturator

Spirometers

Spacelabs Medical Burdick Products.



Spirotouch Spirometer

Now all of your spirometry needs can be met at the touch of a finger – Spirotouch. Whether your application demands a full featured portable or a complete office based system. Spirotouch adapts to the physician's office, bedside or occupational healthcare setting. It boasts new features such as a dynamic touch screen, context sensitive help screens and automatic operation. Compact and portable, Spirotouch's rechargeable battery allows more than eight hours of continuous use -Large bright graphic display for easy viewing under any lighting condition

- Real-time graphs provide instantaneous diagnostic information
- Simple, intuitive operation
- Built-in printer for spirometry on the go
- Choice of animated patient incentives
- Large memory for storage of up to 145 tests

BUR100616 Spirotouch
 BUR100612 Spirometer
 Filter Disk 25/package
 BUR100287 Spirometer
 Mouthpiece and
 Pulmonary Filter PFT
 Refill Kit, 48 of each/kit

Puritan Bennett



Renaissance II Spirometry System

Sophisticated, full-featured, portable Spirometer, Real time graphic display. 1,000 patient memory prints 8-1/2" x 11" color reports, Operates from three different power sources: AA, custom-rechargeable or Ni-cad

P-000700-00 Spirometer,
 10 FSII flow sensors
 P-000710-00 Base Station
 P-000715-00 Starter Kit –
 Spirometer, base, 10 flow

sensors, AC adapter,
 syringe adapter, 2 nose
 clips

P-000720-00 Kit
 w/printer – Spirometer,
 base, 10 flow sensors,
 printer, AC adapter,
 syringe adapter, 2 nose
 clips

P-000725-00 Kit
 w/Syringe – Spirometer,
 base, 10 flow sensors,
 calibration syringe, AC
 adapter, syringe adapter,
 2 nose clips

P-000-730-00 Complete
 Kit – Spirometer, base, 10
 flow sensors, printer,
 calibration syringe, AC
 adapter, syringe adapter,
 2 nose clips

P-495219-00 Battery
 Pack

Spirometry Accessories

PURP-000300-00 VS300
 3 liter, calibration Syringe

PURP-006548-00 Nose
 Clips, 25/bx

PURP-000220-00 Basic
 Starter kit: Includes

PB100 Spirometer,
 PB110 Base station,

PURPB130 Memory
 Card, FS200 sensors

PURP-00221-00 Basic
 Starter Kit, w/ addition of
 VS300 calibration
 Syringe

PURP-000201-00 FS200
 Disposable Expiratory
 Flow sensors, 50/bx

PURP-000200-00 FS200
 Disposable Expiratory
 Flow Sensors, 250/bx

PURP-000250-00
 Disposable Bi-directional
 Flow Sensors, 25/box
 PURP-000251-00
 BD250 Disposable Bi-
 directions Flow sensors,
 100/bx

PURP-000250-00
 Disposable Bi-directional
 Flow Sensors, 25/box
 PURP-000251-00
 BD250 Disposable Bi-
 directions Flow sensors,
 100/bx

Sterilization

Ritter By Midmark



Ritter M7 Speedclave Self-Contained Steam Sterilizer

Sterilize wrapped or unwrapped instruments, linen packs, and glassware. Three-prong plug-in operation for 115V, No plumbing, M-7 has its own reservoir. Low water cutout switch prevents "dry-burning" or operating when not enough water is present. Reset buttons. Cycle time regulates steam time from 0 to 30 minutes. Automatic shut-off when cycle is complete, Safety door protects you. Preset the temperature at a desired setting between 240" and 270" F. Three corrosion-resistant stainless steel trays are standard. Chamber dimension is 7-1/2" x

15". Dimensions are 19"L x 14"W x 13"H

MIDM7001 Speed
clave



Ritter M9 Ultra Clave Automatic Sterilizer

Press the Start/Standby button and the M9 automatically fills with the correct amount of water. Select the cycle desired. Cycle time from start to finish with unwrapped instruments is 22 minutes cold and 12 minutes hot. The door opens automatically to dissipate steam and dry instruments completely. Printer accessory provides a permanent record of cycles, times, temperatures and pressures. Two (each) large and small trays come standard. Chamber dimensions are 9" x 15". Counter area is 16-1/2" x 14"

MIDM9001 Ultra
Clave



Ritter M11 Ultra Clave Automatic Sterilizer

Also the one with the door that opens automatically! The M11, like the M9, features a door that opens automatically making sterilization a set-it-and-forget-it function. It's the newest, largest addition to the Ritter line of automatic sterilizers and in this case... biggest also means best

MIDM11001 Ultra Clave

Autoclave Cleaner Speed Clean

MID002-0396-01 16oz
12/cs

TABLES

Ritter By Midmark



Ritter 104 Examination Table

Patient positioning is easy with the Ritter 104 Examination Table. It's spring-assisted backrest means no more pulling or lifting patient's into position. Hide-away stirrups pull out and lock

in numerous positions. The 104 can be customized with numerous positions. The 104's versatility is further enhanced with an electrical receptacle, non-skid locking footstep and five large drawers for storage. Select from among numerous accessories like pelvic tilt, stainless steel irrigation pan or knee crutches to customize the 104. Steel construction makes the 104 durable for years of service. Available in sand Grey

RIT104035SDSD 104
Exam Table, no Heat
RIT104037SDSD 104
Exam Table, Heat



405 Power Examination Table

The 405 Power Examination Table raises effortlessly from flat to a 75o chair position. Hand-held Commander controls table height and back angle. A foot control accessory allows height and back positioning at the touch of a toe. Vinyl-coated stirrups lock in one of three lateral positions and an infinite

number of lengths. Paper roll holder is concealed in the rear of the back section. Pelvic tile (standard) provides 50° of angle to the seat section. Optional, factory-installed drawer heater available, a large, autoclave polystyrene treatment pan pulls out from under the extra wide foot extension. ISO 9001 certified. 24-1/2" – 36-3/4"H x 54-15/16" – 73-7/8"L x 27"W

MID405 Exam Table



411 Universal Power Table

The 411, adjusts height tilt, back and leg section angles simultaneously, in one-third the time required by sequentially positioning power tables. At the touch of a button, the 411 moves quickly to one of eight precise positions you pre-program into the on-board computer. The hand-held control plugs into the command console or head end of the pedestal base. A

convenient Auto Return button lowers table to the low height position. Several accessories may be added to the table, making in the perfect table for examination, procedures and minor surgery in any size practice

MID411 Power Table

Ultrasound

Siemens



SONOLINE Adara

Sonoline Adara platform reflects what in an ultrasound service. It is an affordable, light-weight (155 pound or 70kg), highly mobile system with premium grayscale imaging. The primary applications of the system include: general abdominal imaging, OB/GYN, urology, orthopedics, emergency medicine, and ear, nose and throat imaging. This mobile grayscale ultrasound has user-focused design which gives you better

exam results with less effort and enhanced image quality

SMN5927301 Sonoline Adara with Vidprint



SONOLINE Sienna

Sonoline Sienna is a next-generation digital ultrasound system that incorporates the latest advances in broadband transducer technology, digitally controlled beamformer, and signal processing design to provide outstanding spatial and contrast. User-defined Quick Sets instantly configure the system to the specific preferences of each user providing consistent, high-quality images and saving valuable time

SMN5684738 Sonoline Sienna, 115V

SMN7294171 Sonoline Sienna LC, 115V



SONOLINE Prima

Sonoline Prima is an innovative ultrasound system that integrates the latest technological advancements in a compact and mobile design. The system is enhanced with a panoramic digital scan convert, dual echo processing, dynamic receive focusing and dynamic aperture control for improved resolution and signal-to-noise ratio. SONOLINE Prima with an integrated video printer is just 19" wide and can be moved conveniently and quickly anywhere you need to your exams

SMN4900606
SONOLINE Prima w/
B/W Printer

Biosound Esaote



Megas Ultrasound System

192 channels, portable 48lb mainframe, digital beamformer, B-mode, M-mode, color and power mapping, spectral PW & CW Doppler. Multi-specialty applications w/phased sector, linear, and curved array probes

BIOMEGAS Megas
Ultrasound System

Megas Probes

BIOPA122E Probe
Megas 5.0/7.5mhz Pa Mf
BIOPA220E Probe
Megas 2.5/3.5mhz Pa Mf



AU3 System

Ideal for the large group, mobile service or hospital, The AU3 system features a powerful, futuristic extended spectrum Platform (ESP) that supports multiple advanced electronic array imaging technologies for unprecedented image

quality. Doppler and CFM sensitivity and clinical versatility, Offering wide-band/high-density linear, curved and annular phased array imaging, The AU3 was developed to address numerous cardiovascular and general imaging applications of hospitals or group practice. It is ideal for the shared services environment

BIO9807050021
AU3 System



Caris Plus Ultrasound System

64 channels, portable 35lb, mainframe, digital beamformer, B-mode, M-mode, color and power mapping, spectral PW & CW Doppler, Multi-specialty applications w/phased sector, linear, and curved array probes

BIO Caris plus
Ultrasound System

Waste Cans

DETECTO



Step-on Waste Can

Built for durability by
Detecto in the USA.

Features, includes folded
and hemmed edges,
heavy-duty foot pedal and
linkage, self-closing lid
with full gasket, internal
bail mechanism.

Available in non-rusting
stainless steel, Each step-
on can includes two bi-
lingual “Bio-Hazard” and
“Infectious Waste” self-
adhesive labels for OSHA
requirements

DETC-24 Stainless
Steel, 24 qt. (6 gal)