


## Nursing

### Nonsterile Gloves

#### Allegiance



#### Flexam Nonsterile Latex Exam Gloves

Our ambidextrous, lightweight gloves give you the protection you need while keeping you comfortable. Featuring a natural latex formulation, Flexam latex exam gloves are powdered for easy donning and meet ASTM D3578 requirements. Average length is 9.1". 100 gloves per box, 20 boxes (2000 gloves) per case

#### With Beaded Cuff

BGH8826 Small  
BGH8827 Medium  
BGH8828 Large

#### With Fitted Cuff

BGH8856 Small  
BGH8857 Medium  
BGH8858 Large



#### Flexam Powder-Free Latex Exam Gloves

Made of natural rubber latex for durability and comfort, our Nonsterile Flexam powder-free gloves offer excellent barrier protection and reduce the risk of skin irritations. The ambidextrous gloves fit securely while providing tactile sensitivity, durability and comfort. Multiple finishes help meet your tactile needs. Gloves contain 50 micrograms or less of total water-extractable protein per gram. They have a 9.1" average length. Caution: safe use of this glove by or on latex-sensitized individuals has not been established. 100 gloves per box 10 boxes (1000 gloves) per case

#### Textured Micro-Roughened Finish Gloves with Polyurethane

##### Coating


BGH8836 Small  
BGH8837 Medium  
BGH8838 Large

#### Positive Touch maximum Grip Gloves with polyurethane Coating

BGH8841 Small  
BGH8842 Medium  
BGH8843 Large

#### Sure-Grip Medium Grip finish Gloves with Beaded Cuff

BGH8876 Small  
BGH8877 Medium  
BGH8878 Large



#### Synthetic Examination Gloves

Contain no natural rubber latex proteins or chemical accelerators. These ambidextrous gloves offer excellent tactile sensitivity and provide maximum dexterity. Powdered for smooth easy donning. They meet or exceed ASTM barrier properly limits and pass ASTM Bacteriophage Penetration Test F1671. Available sterile or nonsterile gloves have a 9" average length and 3.3mil average finger thickness

#### Sterile pairs in peel-open pouch with inner wrap, 200pr/cs

BGH2D7012LI  
Medium  
BGH2D7013LI  
Large

#### Nonsterile singles in dispenser box, 100/bx, 1000/cs

BGH88731 Large


BGH88471 X-large

**Powder-Free Synthetic Examination Gloves**

Contain no natural rubber latex proteins. Available in four sizes, our ambidextrous powder-free synthetic exam gloves feature seamless vinyl construction, increased product life and superior resistance to tears. They have excellent tactile sensitivity and provide maximum dexterity. Gloves are Nonsterile and meet or exceed ASTM barrier properly limits. They have a 9" average length and 6.9-mil average finger thickness. Packaged as singles in a dispenser box. 100 gloves per box, 10 boxes (1000 gloves) per case

BGH8886I Small  
BGH8887I Medium  
BGH8888I Large

**Ansell**



**Micro-Touch Powdered Latex Exam Gloves**  
Premium quality white natural rubber latex examination gloves, powdered with

absorbable dusting powder, U.S.P. Superior comfort and smooth finish Non sterile with fitted cuff. 100/bx, 10bx/cs

ANS5721 Small  
ANS5741 Medium  
ANS5761 large



**DermaClean X-AM Powder-Free Exam Gloves**

Cream natural rubber latex examination gloves. Helps eliminate the risks and costs associated with powder contamination and interference. Textured wet-grip surface and beaded cuff. 100/bx, 10bx/cs

ANS3181 Small  
ANS3182 Medium  
ANS3183 Large



**Perry X-AM Powdered Exam Gloves**  
Premium quality white natural rubber latex examination gloves, powdered with absorbable dusting

powder, U.S.P. satin smooth finish. Non-sterile with beaded cuff. 100/box, 10bx/cs

ANS5760802 Medium  
ANS5760803 large

**Cypress Medical Products**

Cypress Medical products, LP offer a full line of quality examination gloves to fit all health professional needs. Cypress remains committed to the health care professional by providing the highest quality and most widely used examination glove types on the market, while continually researching new materials to bring the health care professional the very best in fit, comfort and durability in barrier protection. Premium or economical, latex or synthetic, smooth or textured-Cypress Medical offers an examination glove to fit any health professional's needs



**Synthesis Synthetic Vinyl, Powdered**

Latex-free, 100 gloves per box

PYP25-02 Small  
PYP25-04 Medium  
PYP25-06 Large  
PYP25-08 X-large



**Synthesis Synthetic Vinyl, Powder-Free**

Latex-free. 100 gloves per box

CYP25-92 Small  
CYP25-94 Medium  
CYP25-96 Large  
CYP25-98 X-Large

**Tillotson Healthcare Corporation**

**Preserve With Aloe Vera Gloves**

Powder free natural rubber latex gloves with aloe Vera added to the inner polymer lining. 100/bx, 10bx/cs

HP02-4400 X-Small  
HPI02-4401 Small  
HPI02-4402 Unisize  
HPI02-4403 Large  
HPI02-4404 X-Large (90/bx)

**Ultra Preserve With Aloe Vera Gloves**

Made from Allotex enzyme treated latex. Enhanced comfort and fit with aloe Vera added to

the inner lining 100/bx, 10bx/cs

HPI02-5350 Small  
HPI02-5351 Unisize  
HPI02-5352 large



**AccuTouch Latex Gloves**

Powder free natural rubber latex gloves. Non-chlorinated, polymer lined and micro textured. 100/bx, 10bx/cs

HPI02-4301 Small  
HPI02-4302 Medium  
HPI02-4303 large  
HPI02-4304 X-Large (90/bx)

**Next generation Vinyl Gloves**

Stretch Poly Vinyl Chloride (PVC) with no donning powder added. Does not contain natural rubber latex. 100/bx, 10bx/cs

HPI02-7500 Small  
HPI02-7501 Medium  
HPI02-7502 Large  
HPI02-7503 X-Large



**Formula One Non-Sterile Gloves**

Soft and Supple natural rubber latex gloves, 100/bx, 20bx/cs

HPI02-5001 Unisize  
HPI02-5200 Small  
HPI02-5002 Large

Sterile Pairs, 50/bx, 8bx/cs

HPI02-5005 Unisize  
HPI02-5008 Small  
HPI02-5006 large

Sterile Singles, 100/bx, 8bx/cs

HPI02-5003 Unisize  
HPI02-5007 Small  
HPI02-5004 Large



**Powder Free Plus Non Sterile Gloves**

Premium soft and supple non-chlorinated powder free protection 100/bx, 20bx/cs

HPI02-4500 Small  
HPI02-4501 Medium  
HPI02-4502 Large

**Ultra Care 2 Non-Sterile exam Gloves**

Starch free non-sterile exam gloves made form Allotex, an enzyme treated natural rubber latex, 100/bx, 10bx/cs

HPI02-5321 Unisize  
HPI02-5320 Small  
HPI02-5323 Large



**Ultra Care PolyTex Blend Non-Sterile Gloves**

One-of-a-Kind PolyTex Blend Starch Free exam Gloves 100/bx, 20bxs/cs

- HPI02-5301 Unisize
- HPI02-5302 Large
- HPI02-5300 Small



**Entrust Powder free Exam Gloves**

Powder free exam gloves with reliable grip control. 100/bx, 20bxs/cs

- HPI02-3000 Small
- HPI02-3001 Medium
- HPI02-3002 Large



**Dental Care Gloves**

Quick donning starch free PolyTex lined dental exam gloves 100/bx, 10bxs/cs

- HPI02-5401 Unisize
- HPI02-5400 Small
- HPI02-5402 Large

**Non-latex Powder Free Nitrile Gloves**

For sensitive hands  
100/bx, 10bxs/cs

- HPI02-6400 Small
- HPI02-6401 Medium
- HPI02-6402 Large
- HPI02-6403 X-large (90/bx)

**Multi Care Synthetic Gloves**

Soft and supple synthetic PVC protection 100/bx, 10bxs/cs

**Powder free**

- HPI02-7000 Small
- HPI02-7001 Medium
- HPI02-7002 large

**Lightly Powdered**

- HPI02-2000 Small
- HPI02-2001 Medium
- HPI02-2002 Large

**MAXXIM Medical**



**Maxxim medical Examination Gloves**

All of Maxxim's non-latex gloves have proven viral barrier protection against blood borne pathogens. Unique, proprietary formulations offer increased touch sensitivity, excellent strength and are available in a variety of sizes. Exam gloves include SensiCare (non-latex

medical exam gloves), SensiCare Powder Free (non-latex powder free medical exam gloves), and TruTouch (vinyl medical exam gloves)

**SensiCare Non-Latex Exam Gloves**

- MXX484001 Small
- MXX484002 Medium
- MXX484003 Large
- MXX484004 X-large

**SensiCare PF Non-Latex Exam Gloves**

- MXX484301 Small
- MXX484302 Medium
- MXX484303 Large
- MXX484304 X-Large

**Tru-Touch Clear Vinyl Exam Gloves**

- MXX480201 Small
- MXX480203 Medium
- MXX480205 Large
- MXX480207 X-Large

**SAFESKIN**



**Safeskin PFE Powder Free latex Exam Glove**

Powder-free eliminates powder associated transfer of latex protein and chemical allergens Fully textured finish; exceptional wet and dry grip

Protein content: This latex glove contains 50 micrograms or less of total water extractable protein per gram

Lower levels of chemicals: Reduces the risk of irritant or allergic contact dermatitis (Type IV)

Tapes and labels: Do not stick or tear gloves

Beaded cuff: For added Strength and easier donning

**100/bx, 10bxs/cs**

KC110 X-small


KC220 Small

KC330 Medium

KC440 Large

**90/bx, 10bxs/cs**

KC550 X-Large



**Safeskin PFE-XTRA Powder Free Latex Glove**

Powder-Free: Eliminates powder associated transfer of chemical allergens

Extra length & thickness: 11" in length and 9 mils thick in fingertips

Meets OSHA, ASHP and ONS guidelines for safe handling of

chemotherapeutic drugs Fully textured finish: Exceptional wet and dry grip

Protein content: This latex glove contains 50 micrograms or less of

total eater extractable protein per gram

Lower levels of chemicals: Reduces the risk of irritant or allergic contact dermatitis (Type IV)

Beaded cuff: added strength and easier donning

Tapes and labels: Do not stick or tear gloves

**50/bx, 10bxs/cs**

KC50500 X-small

KC50501 Small

KC50502 Medium

KC50503 Large

KC50504 X-large



**Safeskin Satin Plus Powder Free Latex Exam Glove**

Powder-Free: Eliminates powder associated transfer of latex protein and chemical allergens

Smooth Finish: Exceptional wet and dry grip

Protein Content: This latex glove contains 50 micrograms or less of total water extractable protein per gram

Lower Levels of Chemicals: reduces the risk of irritant or allergic contact dermatitis (Type IV)

Tapes and labels: Do not stick or tear gloves

Beaded cuff: For added strength and easier donning

**100/bx, 10bxs/cs**

KCSP2110 X-small

KCSP2220 Small

KCSP2330 Medium

KCSP2440 Large

**90/bx, 10bxs/cs**

KCSP2550 X-Large



**Safeskin Powdered Exam Glove**

Powdered: Easy to don and comfortable to wear

Lower Levels of Chemicals: reduces the risk of irritant or allergic contact dermatitis (Type IV)

Smooth Finish:

Exceptional wet and dry grip

Beaded cuff: For added strength and easier donning

**100/bx, 10bxs/cs**

KC100 X-small

KC200 Small

KC300 Medium

KC400 Large

**90/bx, 10bxs/cs**

KC500 X-large



**Safeskin Powder Free Purple Nitrile Exam Glove**

Powder-Free: Eliminates powder associated transfer of chemical allergens

Natural rubber latex free  
Alternative: Made from 100% Nitrile, a synthetic co-polymer, eliminates the risk of an immediate type. Hyper sensitivity (Type I) allergic reaction to latex protein

Lower levels of chemicals: Reduces the risk of irritant or allergic contact dermatitis (Type IV)

Tapes and labels: Do not stick or tear gloves

Textured Fingertips:

Enhanced wet and dry grip

Beaded Cuff for added strength and easier donning

**100/bx, 10bx/cs**

KCN55080 X-small

KCN55081 Small

KCN55082 Medium

KCN55083 Large

**90/bx, 10bx/cs**

KCN55084 X-large

**Safeskin Powder Free Purple Nitrile-Xtra Exam Glove**

Same as above, but 12" in length

KCN50602 Medium



**Safeskin Blue Nitrile Glove**

Textured Fingertips:  
Enhanced wet and dry grip

Powder-Free: Eliminates powder associated transfer of chemical allergens

Natural rubber latex free  
Alternative: Made from 100% Nitrile, a synthetic co-polymer, eliminates the risk of an immediate type. Hyper sensitivity (Type I) allergic reaction to latex protein

Lower levels of chemicals: Reduces the risk of irritant or allergic contact dermatitis (Type IV)

Beaded Cuff for added strength and easier donning

9-1/2" Length

**100/bx, 10bx/cs**

KC53100 X-small


KC53101 Small

KC53102 Medium

KC53103 Large

**90/bx, 10bx/cs**

KC53104 X-large



**Safeskin Synthetic**

Powder-Free: Eliminates powder associated transfer of chemical allergens

Natural rubber latex free

Alternative: Made from 100% Nitrile, a synthetic co-polymer, eliminates the risk of an immediate type. Hyper sensitivity (Type I) allergic reaction to latex protein

Lower levels of chemicals: Reduces the risk of irritant or allergic contact dermatitis (Type IV)

Beaded Cuff for added strength and easier donning

10" Length

**100/bx, 10bx/cs**

KC55030 X-small

KC55031 Small

KC55032 Medium

KC55033 Large

**90/bx, 10bx/cs**

KC55034 X-large

**Sempermed**



**Polymed Powder-free latex Examination Gloves**

Look what separates POLYMED from other gloves. It's called fusion bonding. Our special manufacturing process actually infuses synthetic

polymers into the latex. This creates a multi-layer glove that combines the elasticity of latex and the strength of copolymers protecting against cracking and peeling in normal usage. Without fusion bonding, copolymers can fracture, creating particulate residue, and possible delamination 100/bx, 10bx/cs

SATPM101 XS  
 SATPM102 S  
 SATPM103 M  
 SATPM104 L  
 SATPM105 XL



**Satari Low Powder Latex Examination Gloves**

Providing optimal protection, Sempermed's Satari Latex Powdered gloves are easy to don because of the special light coating of U.S.P. absorbable dusting powder. And because the Satari latex powdered glove is PH balanced, there is less irritation and greater comfort. Satari gloves provide the ideal balance between tactile sensitivity and protection, economy and value. 100/bx, 10bx/cs

SATP102 S

SATP103 M  
 SATP104 L



**Sempermed Supreme Powder-Free latex Surgical Gloves**

Designed with the surgeon in mind. Surgeon's hands are the most delicate instruments in the surgical suite. They need to be prepared. And protect that's why we developed Sempermed Supreme surgical gloves. These unique gloves are anatomically de-signed to reduce hand fatigue and increase comfort. They are powder-free to eliminate aerosolization and skin irritation. A special multi-layer copolymer composition provides optimal protection during surgery. 40/bx, 6bx/cs

SATSPF650 6-1/2  
 SATSPF700 7  
 SATSPF750 7-1/2



**SemperCARE Vinyl Powder-Free Examination Gloves**

The look and feel of a latex glove minus the

latex Strength durability sensitivity fit the SemperCare Vinyl glove is designed as the ultimate solution for those who are at risk of becoming sensitive to natural rubber latex. Its cream color and new resin is an upgrade from regular synthetic vinyl gloves, giving the SemperCare Vinyl glove the look and feel of latex without the latex allergens. It offers exceptional strength and protection, yet maintains maximum tactile sensitivity. What's more, this composition contributes to a better fit and feel, while eliminating odor and skin irritation 100/bx, 10bx/cs

SATSCVNP102 S  
 SATSCVNP103 M  
 SATSCVNP104 L



**Sempermed Synthetic Surgical Gloves**

Eliminating latex protein without sacrificing strength, fit, feel or function Synthetic non-latex examination gloves by Sempermed provide the solution for professionals and their patients who are at risk becoming sensitive to

natural rubber latex, while maximizing tactile sensitivity Strong and durable. The minimal chemical residual in Sempermed's synthetic Glove eliminates odor and skin irritation

**Powdered, 100/bx, 10bx/cs**

- SATEVP102 S
- SATEVP103 M
- SATEVP104 L

**Non-powdered, 100/bx, 10bx/cs**

- SATEVNP101 XS
- SATEVNP102 S
- SATEVNP103 M
- SATEVNP104 L
- SATEVNP105 XL



**SemperCARE Powder-Free Nitrile Examination Gloves**

SemperGuard Industrial Nitrile powder-Free Gloves are the ideal solution for individuals who are at risk of becoming sensitive to natural rubber latex. Nitrile gloves are made from a synthetic polymer exhibiting rubber-like characteristics when vulcanized. They are designed to provide the fit, feel, tensile and tear strengths and other properties of natural rubber latex gloves without the protein allergens. Nitrile gloves

effectively dissipate electrostatic charges and have a higher solvent resistance than natural rubber latex gloves, 100/bx, and 10bx/cs

- SATINIPFT102 S
- SATINIPFT103 M
- SATINIPFT104 L

**Sterile Gloves**

**Allegiance**



**Flexam Sterile Latex Examination Gloves**

Our ambidextrous latex exam gloves offer long wear with excellent fit and dexterity. Powdered for easy donning the gloves has a 9.1" average length and 7.5mil average finger thickness Available in peel-pouch single or pairs

**Pairs with Inner Wrap, 50pr/bx, 200pr/cs**  
 BGH2D7028 Small  
 BGH2D7029 Medium  
 BGH2D7030 Large  
**Singles in Individual Peel Pouch with on Inner Wrap, 100/bx, 1000/cs**

- BGH8861A Medium
- BGH8862A Large



**Sterile Procedure Gloves**

Approximately 3" longer than exam gloves for extra protection Hand-specific for added comfort. Our procedure gloves are ideal for sterile exam room procedures and are powdered for smooth, easy donning. They offer excellent tactile sensitivity and are a cost-effective alternative to regular surgical gloves for noninvasive procedures. The patented molded design helps reduce palm strain, and the anatomically designed fingers allow and excellent range of movement. Gloves feature an anti lip finish for a firm grip. A mechanically locking cuff guards against roll down. Gloves have a 12" average length and 9.1-mil average finger thickness. Packaged in peel-open pouches for sterility. 50 pairs gloves/box, four boxes (200 pairs) per case



BGH2D7122T  
6-1/2  
BGH2D7123T 7  
BGH2D7124T  
7-1/2  
BGH2D7125T 8



**Triflex Latex Surgical Gloves**

Our Triflex surgical gloves offer excellent tactile sensitivity and an anti slip finish for a firm grip. They feature a patented natural anatomical design that reduces thumb and palm strain for greater comfort. These gloves meet or exceed ASTM thickness requirements and have a mechanically locking cuff to prevent roll down. Gloves have a 12" average length and 9.1-mil average finger thickness. Supplied in easy-to-store Allegiance pack, 40 pairs of gloves per box, five boxes (200 pairs) per case

BGH2D7250 5-1/2  
BGH2D7251 6  
BGH2D7252 6-1/2  
BGH2D7253 7  
BGH2D7254 7-1/2  
BGH2D7255 8


BGH2D7256 8-1/2  
BGH2D7257 9



**Ultraderm Latex Surgical Gloves**

These Gloves feature an anti slip finish that provides a firm grip when handling instruments, and the tapered fingers offer an excellent fit. The mechanically locking cuff fits snugly to help prevent roll down. Cuff is tapered with reinforced ribs to distribute pressure evenly and to help keep the glove from slipping. Gloves contain reduced concentrations of chemical accelerators that are known to be potential sensitizers. Gloves have a 12" average length and 8.7-mil average finger thickness. Packaged in easy-to-store Allegiance pack. Caution: Safe use of this glove by or on latex-sensitized individuals has not been established. 40 pairs of gloves per box, five boxes (200 pairs) per case


BGH2D7264 7-1/2



**Perry X-AM Sterile Exam Gloves**

Premium quality white natural rubber latex examination gloves  
Lightly powdered, sterile with based cuff, 50/bx, and 4bxs/cs

ANS5761402 Medium, pairs  
ANS5761102 Medium, Singles  
ANS5761103 Large, Singles



**No Powder Sterile Surgical Gloves**

Highly cost-effective, cream natural rubber latex surgical gloves  
Powder-free to help eliminate powder-related complications  
Textured wet-grip surface with unbeaded reinforced cuff  
50pr/bx, 4bx/cs

ANS8603 6.5  
ANS8604 7  
ANS8605 7.5  
ANS8606 8

**Ansell**



**Perry natural Powder-Free Sterile Surgical Gloves**

Cream natural rubber latex surgical gloves Powder-free with beaded cuff 50pr/bx, 4bx/cs

- ANS5717002 6.5
- ANS5717003 7
- ANS5717004 7.5
- ANS5717005 8



**Micro-Touch Sterile Surgical Gloves**

White natural rubber latex surgical gloves Powdered with absorbable dusting powder, U.S.P. Anti-slip bisque finish improves grip. Inverted, beaded cuff multiple quality checks ensure high quality and user confidence, 50 pr/bx, and 4bx/cs

- ANS5860 6
- ANS5865 6.5
- ANS5870 7
- ANS5875 7.5
- ANS5880 8
- ANS5885 8.5
- ANS5890 9



**NEUTRALON Powdered Brown Sterile Surgical Gloves**

Brown natural rubber latex surgical gloves Powdered with absorbable dusting powder, U.S.P. Micro-thin inner synthetic membrane provides an extra layer on the inside of the glove. Bisque finish for enhanced grip. Inverted, beaded cuff helps prevent roll down. 50 pr/bx, 4bx/cs

- ANS5375 7



**The original Perry Style 42 sterile Surgical Gloves**

White natural rubber latex surgical gloves. Satin smooth finish for maximum tactile sensitivity. Powdered with absorbable dusting powder, U.S.P. Beaded cuff 60ug/g or less of water soluble protein 50pr/bx, 4bx/cs  
ANS5711101 6  
ANS5711102 6.5  
ANS5711103 7  
ANS5711104 7.5

- ANS5711105 8
- ANS5711106 8.5
- ANS5711107 9



**Sensi-Touch Sterile Surgical Gloves**

White natural rubber latex surgical gloves Powdered with absorbable dusting powder, U.S.P anti-slip bisque finish improves instrument handling Regular, reinforced cuff 50pr/bx, 4bx/cs

- ANS7825 7.5

**REGENT**



**Biogel Powder-Free latex Surgical Gloves**

Sterile Powder-free latex surgical gloves with Biogel coating for easy donning with damp or dry hands 100% air inflated. Curved finger design for added comfort. Contain 50 micrograms or less of water-extractable protein per gram. Tested with selected cytotoxic drugs commonly used in

chemotherapy. (1) Low dermatitis potential. (2) 50 pairs/bx, 4bx/cs

30465	6.5
30470	7
30475	7.5
30480	8



### **Biogel M – Micro-Textured for Improved Grip**

Sterile powder-free latex surgical gloves with enhanced grip. Biogel coating for easy donning with damp or dry hands 100% air inflated. Straight finger design with additional width in the palms Contain 50 micrograms or less of water-extractable protein per gram Tested with selected cytotoxic drugs commonly used in chemotherapy. (1) low dermatitis potential. (2) 50 pairs/bx, 4bx/cs

30575	7.5
-------	-----

**MAXXIM**  
Medical



### **Maxxim Surgical gloves**

Maxxim Medical offers a complete line of surgical gloves uniquely formulated to offer a tiered reduction of exposure to natural rubber latex (NRL) protein content. Healthcare workers can select the glove features needed without compromising the benefits associated with the reduction or elimination of exposure to higher concentrations of NRL protein

Eudermic Latex Surgical Glove  
MXX376189 Size 7.5

Tradition Latex Surgical Glove  
MXX376199 Size 7.5

## **SAFESKIN**

### **Safeskin Satin Plus Powder Free Sterile Single/Pair**

Powder-free: Eliminates powder associated transfer of chemical allergens  
Sterilized: In accordance with AAMI Guidelines  
Smooth finish:  
Exceptional wet and dry grip  
Protein content: This latex glove contains 50 micrograms or less of total water extractable protein per gram

Lower levels of chemicals: Reduces the risk of irritant or allergic contact dermatitis (Type IV)  
Beaded cuff: Added strength and easier donning  
Tapes and labels: Do not stick or tear gloves

### **Sterile Singles. 100 singles/box, 4boxes/case, 400 singles/case**

KC50301	Small
KC50302	Medium
KC50303	Large

**Sterile Pairs, 50 pairs/box, 4boxes/cs**

KC50401	Small
KC50402	Medium
KC50403	Large

### **Safeskin PFS Powder Free Latex Surgical Glove**

Powder-free: Eliminates powder associated transfer of chemical allergens  
Protein content: This latex glove contains 50 micrograms or less of total water extractable protein per gram  
Lower levels of chemicals: Reduces the risk of irritant or allergic contact dermatitis (Type IV)  
Bisque finish: provides exceptional wet or dry grip  
Donning agent: Facilitates wet and dry donning with exceptional mildness to the skin

50 pairs/box, 4boxes/case,  
200pairs/case

KC54365 6-1/2  
KC54370 7  
KC54375 7-1/2  
KC54380 8

## Utility Gloves

### Ansell



#### Latex/Nitrile Utility Gloves

Latex-Nitrile blend non-sterile utility gloves  
Protects against acids, oils, disinfectants, detergents and abrasives  
Flock-lined for comfort and easy don ability.  
Embossed sure-grip surface and unbeaded cuff. 12pr/box, 4bx/cs

ANS5789016 Medium

## Allegiance

#### Wall-Mount Glove Dispenser

Universal wall bracket  
glove dispenser

BGH8000R

## Masks, Caps & Shoe Covers

### 3M Health Care



#### 3M 1860 Health Care Particulate Respirator Type N95

NIOSH approved as a Type N95 respirator. It meets CDC guidelines for TB exposure control and is designed specifically for use in a health care setting

MMM1860 Regular,  
20/box, 6bx/cs  
MMM1860S Small,  
20/box, 6bx/cs



#### 3M Fluid Resistant Tie- On Surgical Masks

Made of a comfortable, lightweight network of synthetic non-irritating fibers. Sonic stitching for greater strength and comfort Additional mask folds for assured comfort, 100% fiberglass free. High filtration and asepsis

MMM1818 50/box,  
12bx/cs  
MMM1818FS With anti-fog/antireflective coated face shield, 25/box, 4bx/cs  
MMM1838FS Duck bill, with anti-fog/antireflective coated face shield, 25/box, and 4bx/cs



#### 3m Aseptex Surgical Mask

Non-fiberglass, high-efficiency, disposable molded mask

MMM1800+50/box,  
12bx/cs



#### 3M Earloop Fluid Resistant Surgical Masks

A single-use Earloop mask available with or without face shield

MMM1820 50/box,  
12bx/cs



**3M Filtron High Performance (Duck Bill) Surgical Masks**

Features a fluid resistant construction that is efficient against wearer-generated particles. The unique filter combination allows for easy breathing and reduces fogging of glasses or goggles, while the distinctive shape helps keep the mask away from nose and mouth

MMM1838 50/bx, 12bxs/cs  
MMM1838FS With anti-fog-antireflective coated face shield, 25/bx, 4bxs/cs

**Busse Hospital Disposable**



**Standard Style Shoe Cover**

Features quality and durability at an economical price

BUS340  
150pairs/carton (3 dispensers of 50)

**Kimberly-Clark Corporation**



**Kimberly-Clark Headwear Products**

**Bouffant Cap, Blue**  
KC69221 regular, 300/cs  
KC69231 Large, 300/cs



**Surgical Cap, Blue**  
KC69240 Universal, 300/cs  
**Protective Surgical Cap, Blue**  
KC69520 Universal, 300/cs



**Kimberly-Clark Tecnol Face mask Products**

**Standard Face Masks**  
KC47095 Procedure mask, Pink, 50/bx, 10bxs/cs  
KC47117 High Filtration ISO Mask II, 50/bx, 10bxs/cs  
KC47200 Cone Classic Surgical Mask, 50/bx, 6bxs/cs

KC48201 Classic Surgical mask, 50/bx, 6bxs/cs

**Fluidshield Face Masks Without Face Shields**

KC47107 Procedure Mask, 40/bx, 10bxs/cs

**Fluidshield Face Masks With Face Shields**

KC47137 procedure mask w/ Splashguard Visor, 25/bx, 4bxs/cs

KC47147 Procedure mask w/ Wraparound Splashguard Visor, 25/bx, 4bxs/cs

KC48247 Surgical Mask w/ Wraparound Splashguard Visor, 25/bx, 4bxs/cs

KC48207 Fog-Free Surgical Mask, w/ orange Fluidshield, 50/bx, 6bxs/cs

KC48237 Fog-Free Surgical mask, w/ Wraparound Splashguard Visor, 25/bx, 4bxs/cs

**Specialty Face Masks**

KC48100 Surgical mask, Blue, The Lite One, 50/bx, 6bxs/cs

KC47500 Surgical Mask, Blue, Soft Touch II, 50/bx 6bx/cs

KC49214 Fog-Free Surgical Mask, 50/bx, 6bxs/cs

KC49215 Fog-Free Surgical mask w/ Derma-Touch Tape Green, 50/bx, 6bxs/cs

KC47625 Lazer Surgical Mask Silver, 35/bx, 6bxs/cs

KC47650 Lazer Fog-Free Surgical Mask, Silver, 35/bx, 6bxs/cs  
 KC48220 Surgical mask, Blue Duckbill, 50/bx, 6bxs/cs  
 KC49230 Anti-Fog Surgical Mask Blue, 50/bx, 6bxs/cs

**Face Mask Accessories**


KC41204 Guardall Shield Full Length Visor, 40/cs



**Ulti-Mask**

Hypoallergenic mask made from premium “linen-like” material containing no dyes to help to reduce facial irritation Thin foam band acts as a cushion between the nose and the metal nosepiece and helps reduce eyeglass fogging

KC49310 50/bx, 6bxs/cs



**Cone Classic Surgical Mask**

An enclosed nosepiece wire eliminates the chance of metal falling


onto the surgical site, while a soft knitted headband offers maximum comfort for the wearer with a secure fit. Mask provides excellent filtration at over 95% efficiency

KC47200 50/bx, 300/cs

**Nurse’s Surgical Mask**

Nurse’s mask is created to provide a comfortable, proportionately correct fit for the smaller face. The three-fold design allows necessary fullness

KC48210 50/bx, 300/cs



**So Soft Surgical Mask**

Designed especially for sensitive skin, this surgical mask features special inner facing with a smooth satin finish. The pure white materials do not contain any color dyes, reducing the risk of irritation

KC48390 50/bx, 300/cs



**Procedure Mask**

Easy to don and comfortable to wear, the Procedure Mask is ideal for non-surgical procedures. Mask features a surgical-grade filter medium and has convenient ear loops. 50 per box, 500 per case

KC47080 Blue



**Kimberly-Clark Footwear Products**

**Basic Shoe Cover**

Coverage: Shoe only, materials, protective 3-layer fabric, exposure anticipated, no potential fluid contact

KC69561 Universal, 300/cs

**Heavy Duty Shoe Cover**

Coverage: shoe only, materials, extra thick protective 3-layer fabric for extended wear, exposure anticipated, light fluid contamination

KC69345 X-large, 24/cs

KC69551 Universal,  
300/cs

**X-TRA Traction Shoe  
Cover**

Coverage: shoe only,  
materials: protective 3-  
layer fabric, six adhesive  
traction strips, exposure  
anticipated: light fluid  
contamination

**X-TRA Traction Shoe  
Cover**

KC69352 Universal,  
300/cs

**Basic X-TRA Traction  
Shoe Cover**

KC69252 Universal,  
300/cs

KC69254 X-large,  
80/bx

**HOT/COLD  
PACKS**

**3M Health Care**



**3M ColdHot Pack**

A soft, gel-filled pillow  
that can be used on any  
part of the body to  
provide hot or cold pain  
relief The pack can be  
stored in a freezer until  
ready to use

MMM1570 2 Pack and  
2 covers 4" x 10"

**Allegiance**

**Kwik-Heat Single-Use  
Instant hot Pack**

Kwik-Heat hot pack  
reaches 110oF (43o)  
almost instantly when  
you squeeze it. Each hot  
pack is 6.25" x 8' 12  
packs per case

BGH211 Instant hot  
pack



**Kwik-Kold instant Cold  
packs**

Our single-use instant  
cold packs offer a  
convenient and cost-  
effective method of  
delivering localized cold  
therapy Ready when you  
need them; the cold packs  
require no special  
preparation or  
refrigeration. Store them  
at room temperature and  
activate them in an instant  
with a firm squeeze and a  
gentle shake.

Approximate peak  
temperature and  
therapeutic time duration  
33oF/0.5oC, 30 minute.  
Choose from three sizes

BGH103 Kit, 4' x  
5.75", individually boxed,  
50/cs

BGH102 Junior,  
5.5" x 6", Bulk, 16/cs

BGH107 Regular,  
6.25" x 8", Bulk, 16/cs

BGH101 Regular,  
6.25" x 8", 16/cs

**Reusable Hot/Cold Gel  
Packs**

Our reusable hot/cold gel  
packs are filled with a  
versatile gel product that  
stays flexible during use.  
Packs are designed for  
both hot and cold therapy  
needs. 24 packs per case

BGH70304A Gel pack  
kit, Medium, includes 2  
gel packs and 2  
disposable sleeves 4-1/2 x  
10-1/2

**Chattanooga  
Group, Inc**



**Moist heat steam Packs**

Simple to use-place pack  
in heating unit or in a hot  
pan of water for home  
treatment Remove and  
wrap with a terry cover or  
regular towel, and apply  
to treatment area. Store  
packs in a plastic bag,  
place in your freezer

CHT1006 Standard,  
10"x12"

CHT1002 neck  
contour, 24" long



### Reusable Cold packs

These reusable cold packs eliminate the need for messy ice packs and soggy towels. Simply chill in refrigerator freezer compartment and apply to treatment area. Fit any part of the body Non-toxic silicate gel filter. Flexible vinyl covering is durable and won't become brittle, crack or leak

CHT1500 Standard, 11"x14"  
CHT1508 Neck Contour, 23" long

## Graham-Field

### Disposable instant Cold Packs

Standard, 6'x9"

GF9100 24/cs

### Reusable Hot/Cold Gel Packs

Insulated on one side. Microwavable. May be used hot or cold. Stays flexible

GF9121 Standard 4-1/2" x 10-1/2", 24/cs

## Kimberly-Clark Corporation



### Kimberly-Clark Cold Therapy products

#### Ice Filled Packs

KC33101 Bilateral Facial 12/bx, 2bxs/cs  
KC33150 EyeCare (EENT), 15/bx, 2bxs/cs  
KC33600 Secure-All, Large, 15/bx, 2bxs/cs  
KC33495 Stay-Dry, Small, 5"x12", 25/bx, 2bxs/cs  
KC33620 Soft'n Cold, Large, 6-1/2" x 14", 25/bx, 2bxs/cs



#### Instant Packs

COLD RUSH Instant Cold Pack offers unique dual compartment design to separate the reactive ingredients; packs provide immediate cold therapy. Latex-free

KC59648 small, 4-1/2"x10-1/2", 10/bx, 4bx/cs

## CPS

Consolidated products and Services, Inc

### Instant Cold Compress

Single-use Instant Cold Packs provide convenient and cost effective pain relief, reducing swelling and inflammation. They are quickly activated and deliver a consistent, long-lasting therapeutic temperature

CPT42991 5.5"x8", standard, 24/cs

## NEEDLES & SYRINGES

### Becton Dickinson



Regular Bevel

### Precisionglide Needles

An advanced low-angle bevel design, a unique polishing process and a low-friction lubricant have made BD Precisionglide Needles the world leader in sharpness. Now, a new sharpening method electronically hones the needle surface for a finer, smoother cutting edge. The BD exclusive four step process results in unequalled sharpness and penetration ease that you and your patients will appreciate



**Regular bevel & Specialty use, 100/bx, 1000/cs**

BD305122 25Gx5/8"  
 BD305143 23Gx3/4"  
 BD305145 23Gx1"  
 BD305155 22Gx1"  
 BD305159 22Gx1-1/2"  
 BD305156 18Gx1/2"  
 BD305196 18Gx1-1/2"  
 BD305109 27Gx1/2"  
 BD305111 28Gx1/2"  
 BD305110 26Gx3/8"  
 BD305115 26Gx5/8"  
 Sub-Q  
 BD305127 25Gx1-1/2"  
 BD305124 25Gx7/8"  
 BD305165 21Gx1"  
 BD305167 21Gx1-1/2"  
 BD305129 21Gx2"  
 BD305175 20Gx1"  
 BD305176 20Gx1-1/2"  
 BD305195 18Gx1"  
 BD305125 25Gx1"  
 BD305177 21Gx1"  
 BD305185 18Gx1-1/2" TW  
 BD305186 19Gx1"  
 BD305187 19Gx1-1/2"  
 BD305189 19Gx1-1/2"  
 BD305193 23Gx1"  
 TW IM  
 BD305194 23Gx1-1/2" Sub Q  
 BD305106 30Gx1/2"  
 BD305128 30Gx1"  
 BD305136 27Gx1-1/4"  
 BD305197 16Gx1"

BD305198 16Gx1-1/2"  
 BD305199 18Gx1-1/2"  
 BD305200 19Gx1-1/2" filter  
 BD305201 18Gx1-1/2"



**Syringes**

Sterile, single use syringes with Tip Shields  
 BD309586 3cc, slip Tip, 100/bx, 800/cs  
 BD309585 3cc, Luer-Lok Tip, 100/bx, 800/cs  
 BD301603 3cc, Slip Tip, 100/bx, 400/cs  
 BD309603 5cc, Luer-Lok tip, 100/bx, 400/cs  
 BD301604 10cc, Luer-Lok Tip, 100/bx, 400/cs  
 BD309604 10cc, Luer-Lok Tip, 100/bx, 400/cs  
 BD301627 60cc, slip Tip, 30/bx, 120/cs  
 BD309653 60cc, Luer-Lok tip, 40/bx, 120/cs  
 BD309661 20cc, Luer-Lok Tip, 40/bx, 120/cs



**Syringe/Needle Combo**  
 Luer-Lok Tip. 100/bx, 800/cs

**3cc**  
 BD309570 25Gx5/8"  
 BD309571 23Gx1"  
 BD309572 22Gx1"  
 BD309573 22Gx1-1/4"  
 BD309574 22Gx1-1/2"  
 BD309575 21Gx1  
 BD309577 21Gx1-1/2"  
 BD309578 20Gx1"  
 BD309579 20Gx1-1/2"  
 BD309581 25Gx1"  
 BD309580 18Gx1-1/2"  
 BD309582 25Gx1-1/2"  
 BD309587 26Gx5/8"  
 BD309588 23Gx1"  
 BD309589 23Gx1-1/2"  
**1cc**  
 BD309628 Luer-Lok Tip, 100/bx, 800/cs  
**5cc**  
 100/bx, 400/cs  
 BD309630 22Gx1"  
 BD309631 22Gx1-1/2"  
 BD309632 21Gx1"  
 BD309633 21Gx1-1/2"  
 BD309634 20Gx1"  
 BD309635 20Gx1-1/2"  
**10cc**  
 100/bx, 400/cs  
 BD309640 22Gx1"  
 BD309642 21Gx1"  
 BD309643 21Gx1-1/2"

BD309644 20Gx1"  
BD309645 20Gx1-  
1/2"

**Allergist Trays**

Permanent needle, 25  
syringe per tray. 1000/cs  
BD305535 1/2cc  
27Gx1/2"  
BD305536 1/2cc  
27Gx1/2' intradermal  
bevel-testing  
BD305537 1cc  
26Gx1/2" regular bevel-  
testing  
BD305538 1cc  
26Gx1/2" intradermal  
bevel-testing  
BD305539 1cc  
26Gx3/8" intradermal  
bevel-testing  
BD305540 1cc  
27Gx1/2" intradermal  
bevel-testing  
BD305541 1cc  
27Gx3/8" intradermal  
bevel-testing  
BD305542 1cc  
27Gx3/8" regular bevel-  
testing  
**No Tray**  
100/bx, 5bx/cs  
BD305500 1cc  
28Gx1/2"

**Insulin Syringes**

With permanently  
attached needle  
**1cc**  
BD309309 28Gx1/2"  
1cc, 100/bx, 5bx/cs  
BD309331 27Gx1-  
1/2" 1cc, 40/bx, 4bx/cs  
BD329410 28Gx1/2"  
U100, 100/bxs, 5bx/cs  
BD329411 29Gx1/2'  
Ultra-Fine 100/bx, 5bx/cs

BD329412 27Gx5/8"  
U100, 100/bx, 5bx/cs  
BD329420 28G U100,  
100/bx, 5bx/cs  
**1/2cc**  
BD329461 28Gx1/2"  
1/2cc U100, 100/bx,  
5bx/cs  
BD329465 28Gx1/2"  
1/2cc U100, 100/bx,  
5bx/cs  
BD329466 29Gx1/2"  
Ultra-Fine 1/2cc U100,  
100/bx, 5bx/cs



**TB Syringes & Needles**

1cc Tuberculin syringe  
with detachable needle  
with slip tip  
BD309624 21Gx1",  
100/bx  
BD309626 25Gx5/8",  
100/bx  
BD309625 26Gx3/8"  
intradermal bevel, 100/bx  
BD309623 26Gx3/8",  
with tip shield, 100/bx  
BD309602 With tip  
shield, 100/bx

**Spinal Needles**

Sterile, single use  
BD405074 22Gx2-  
1/2", 25/bx  
BD405148 23Gx5",  
10/bx  
BD405149 22Gx7",  
10/bx  
BD405161 22Gx1-  
1/2", 25/bx  
BD405170 25Gx3",  
25/bx

BD405171 22Gx3",  
25/bx  
BD405180 25Gx3-  
1/2", 25/bx  
BD405181 22Gx3-  
1/2", 25/bx  
BD405182 20Gx3-  
1/2", 25/bx  
BD405184 18Gx3-  
1/2", 25/bx  
BD405211 20Gx6",  
10/bx



**Vacutainer Needles**

Multi-sample thin wall  
needles are individually  
packaged  
BDV367210 22Gx1",  
100/bx, 400/cs  
BDV367211 22Gx1-  
1/2", 100/bx, 400/cs  
BDV367212 21Gx1",  
100/bx, 400/cs  
BDV367213 21Gx1-  
1/2", 100/bx, 400/cs  
BDV367214 20Gx1",  
100/bx, 400/cs  
BDV367215 20Gx1-  
1/2", 100/bx, 400/cs



**ECLIPSE Safety  
Shielding Blood  
Collection Needle**

Integrated, single-use  
blood collection needle

for fast and easy needle protection

**Advanced Protection Technology**

- Single handed activation
- Immediate protection
- No change in technique
- Unique double locking system
- Audible click indicator

**Ease of use**

- Integrated safety needle shield design
- Intuitive in use
- Single-use
- Bevel up orientation

BDV368607 21G, 1-1/4", 48/bx, 480/cs  
BDV368608 22G, 1-1/4", 48/bx, 480/cs



**BD Safety-Lok Syringe**

BD329464 1ml Insulin Syringe with Permanently Attached 29gx 1/2in. ULTRA-FINE Needle 9U-100) 100/bx, 500/cs

BD305553 1ml Tuberculin Syringe with Permanently Attached 27gx1/2in. Needle, 100/bx, 500/cs

BD305554 1ml Tuberculin Syringe with permanently Attached 25gx5/8in. Needle, 100/bx, 500/cs

BD309606 3ml Syringe only 100/bx, 800/cs  
BD309592 3ml syringe with 25gx 5/8in. Needle, 100/bx, 800/cs  
BD309564 3ml Syringe with 23g x 1in. Needle, 100/bx, 800/cs  
BD309593 3ml Syringe with 22g x 1-1/2in. Needle, 100/bx, 800/cs  
BD309596 3ml Syringe with 22g x 1in. Needle, 100/bx, 800/cs  
BD309595 3ml Syringe with 21g x 1-1/2in. Needle, 100/bx, 800cs  
BD305558 5ml Syringe only 50/bx, 400/cs  
BD305559 10ml Syringe only 50/bx, 400/cs  
BD305561 5ml Syringe with 21g x 1-1/2in. Needle, 50/bx, 400/cs  
BD305564 10ml Syringe with 21g x 1-1/2in. Needle, 50/bx, 400/cs



**BD Medsaver Syringe**

You can make a significant impact on reducing your facility's drug costs – typically up to 15% - by using the BD

Medsaver syringe. The Medsaver syringe is designed with virtually no dead space, so there is minimal residual medication in the syringe. The result less waste, lower costs 100/bx, 500/cs

BD305591 25Gx 1 3ml Syringe with needle  
BD305592 25gx 5/8-3ml Syringe with needle  
BD305595 23Gx 1 3ml Syringe with needle  
BD305605 25G x 5/8 1ml Syringe with needle



**BD SafetyGlide Shielding Hypodermic Needles & Syringes**

Activation-assist technology allows for fast and easy, single-handed needle shielding. It can help prevent needle shielding. It can help prevent needle stick injuries by allowing easy and secure shielding of the needle tip immediately after injection. BD SafetyGlide Syringe a line extension of the BD SafetyGlide Shielding Hypodermic Needles, is similar to the larger BD SafetyGlide sizes, BD SafetyGlide Syringe offers single-

handed immediate activation with minimal change in technique

BD305930 1ml Insulin Syringe with 28g 1/2in Needle, 100/bx, 400/cs  
BD305932 1/2ml Insulin Syringe with 28g 1/2in Needle, 100/bx, 400/cs  
BD305935 3/10ml Insulin Syringe with 29g 1/2in Needle, 100/bx, 400/cs  
BD305945 1ml Tuberculin Syringe with 27g 1/2in Needle, 100/bx, 400/cs  
BD305946 1ml Tuberculin Syringe with 26g 3/8in Needle, (IB) 100/bx, 400/cs  
BD305901 25g 5/8in Subcutaneous Needle, 50/bx, 500/cs  
BD305903 1ml Syringe with 25g 5/8in Sub cut, Needle, 50/bx, 400/cs  
BD305904 3ml Syringe with 25g 5/8in Sub cut Needle, 50/bx, 400/cs  
BD305916 25g 1in IM Needle, 50/bx, 500/cs  
BD305902 23g 1in IM Needle, 50/bx, 500/cs  
BD305900 22g 1-1/2in. IM Needle, 50/bx, 500/cs  
BD305915 21g 1in IM Needle, 50/bx, 500/cs

BD305917 21g 1-1/2in. IM Needle, 50/bx, 500/cs  
BD305905 3ml Syringe with 23g 1in IM Needle, 50/bx, 400/cs  
BD305906 3ml Syringe with 22g 1-1/2in IM Needle, 50/bx, 400/cs  
BD305907 5ml Syringe with 22g 1-1/2in IM Needle, 50/bx, 400/cs  
BD305908 10ml Syringe with 22g 1-1/2in IM Needle, 50/bx, 400/cs  
BD63190 1ml Allergist Tray with 27g 1/2in Needle, (RB) 25/Tray 1000/cs  
BD66235 1ml Allergist Tray with 26g 3/8in Needle, (IB) 25/Tray 1000/cs

## OMRON



### Chrome-Plated self-filling Syringe Utility Kit

Includes: 1 double valve assembly with spare parts, 26" tubing with sinker, protection shield and 4 tips-pomeroy, fistula, catheter, kraemer

OMR16-100 Syringe Utility Kit

## EXEL



### Hypodermic Needles

Latex-free, non-toxic, non-pyrogenic, each needle is individually packaged, sterilized and ready for use. Needles and boxes are color-coded for easy identification. Many sizes available

EXL26426 27Gx1-1/2  
EXL26437 30G x 1/2

## Allegiance



### Tru-Cut Soft Tissue Biopsy Needles

Our Tru-Cut biopsy needles are designed for manual capture of high quality tissue samples with minimal trauma to the patient. They have a full 20mm specimen notch, centimeter depth markings that assist in placement and removable stylet for multiple sampling capability A scalpel-sharp surgical steel cutting edge retrieves cleanly cut, high-quality specimens.

Thin-wall needle yields larger core. Individually packaged. 10 needles per case

BGH2N2702X  
14G x 4-1/2"

### 60cc Irrigating Syringes

Our irrigating syringes include a protector cap. Individually packaged sterile. 50 syringes per case

BGH3T7600 Bulb syringe

### 60cc Irrigation Syringe

(Kendall) This bulb syringe is sterile and individually packaged. 50 syringes per case

BGHKB7525 60cc irrigation syringe



### Jamshidi Bone Marrow Needles

The tapered distal tip of our disposable Jamshidi bone marrow needles is designed for successful specimen collection. 10 needles per case

**Needles with standard handle**  
BGHDJ4011X  
11G x 4"

**Futura**  
The New Thinking Safety



### Pro-Ject Blood Needle Holder

Pro-ject needle holders provide a safe and economical mechanism for the disposal of contaminated blood collection needles, Easy one-handed procedure. No need to recap or unscrew from holder

FTAPJ200-10 10/box  
FTAOJ200-25 25/box  
FTAPJ200-250 250/bx  
FTAPJ200-10H For use w/ Hemogard 10/bx  
FTAPJ200-25H For use w/ Hemogard, 25/box  
FTAPJ200-250H For use w/ Hemogard, 250/bx

**Tyco**  
Healthcare

**Kendall**



### Monoject Hypodermic Needles with Polypropylene Hub

Ultra-sharp, tri-beveled, anti-coring, stainless steel needle Translucent, color-

coded luer lock hub. Color-coded cap on rigid pack Autoclavable in rigid pack, for use in O.R/procedure packs sterile

### Regular (A), 100/bx, 1000/cs

KEN8881250032 30x3/4", Light Blue  
KEN8881250313 25x5-8", Red  
KEN8881250305 25x1", Red  
KEN8881250545 25x1-1/2", Red  
KEN8881250255 23x1", Orange  
KEN8881250222 22x1-1/2", Blue  
KEN8881250172 21x1", Lavender  
KEN8881250123 20x1", Pink  
KEN8881250040 18x1" Green  
KEN8881250016 18x1-1/2", Green  
**Intermediate (IN), 100/bx, 1000/cs**  
Ken8881250362 27x1/2", Yellow



### Monoject hypodermic Needles w/ Aluminum Hub

Ultra-sharp, tri-beveled, anti-coring, stainless steel needle Luer lock hub Color-coded needles sheath. White cap identifies short bevel. Autoclavable in rigid

pack for use in  
O.R./procedure packs  
sterile

**Regular (A), 100/bx,  
1000/cs**

KEN8881200508  
27x1-1/4", Yellow  
KEN8881200433  
25x1-1/4"

**Special Sizes (A), 100/bx,  
1000/cs**

KEN8881200441

25x2", Red  
KEN8881200037

16x1-1/2", Gray

**Special Sizes (IN), 100/bx,  
1000/cs**

KEN8881200011

14x1-1/2", Tan



**Monoject 3cc Syringes**

Bold graduations: 0.1cc  
Ultra-sharp, tri-beveled,  
anti-coring needle. Luer  
lock syringe tip on  
combination units. Color-  
coded translucent needle  
hub and rigid pack cap.  
Oriented, compact  
package. Autoclavable in  
rigid pack for use in  
O.R./procedure packs.  
Sterile

**100/bx, 1000/cs**

KEN8881513033

20x1", Pink

KEN8881513132

21x1" lavender

KEN8881513157

21x1-1/2", Lavender

KEN8881513231

22x1", Blue

KEN8881513256

22x1-1/2", Blue

KEN8881513330

23x1", Orange

KEN8881513512

25x5/8", Red

KEN8881513538

25x1", Red

KEN8881513744

27x1-1/4" Yellow

**Syringe Only, 100/bx,  
1000/cs**

KEN8881513918

Regular Luer Tip

KEN8881513934

Luer Lock Tip



**Monoject 6cc Syringes**

Bold graduations: 0.2cc  
Ultra-sharp, tri-beveled,  
anti-coring needle Luer  
lock syringe tip on  
combination units. Color-  
coded translucent needle  
hub and rigid pack cap.  
Oriented, compact  
package. Autoclavable in  
rigid pack for use in  
O.R./procedure packs  
sterile

**50/bx, 500/cs**

KEN8881516135

21x1", Lavender

KEN8881516259

22x1-11/2", Blue

**Syringe Only, 50/bx,  
500/cs**

KEN8881516911

Reg, Luer Tip

KEN8881516937

Luer Lock Tip



**Monoject 20cc Syringes**

Bold graduations: 1cc  
autoclavable in rigid pack  
for use in O.R./procedure  
packs. Oriented, compact  
package sterile. 50/bx,  
300/cs

KEN8881520673

Regular Luer

KEN8881520657

Luer Lock



**Monoject 35cc Syringe**

Bold dual graduations:  
35cc in 1cc increments  
and 1oz in 1/4oz  
increments Autoclavable  
in rigid pack for use in  
O.R./procedure packs  
oriented compact package  
Sterile

**30/bx, 180/cs**

KEN8881535796

Regular Luer

KEN8881535770

Catheter (Irrigation)

KEN8881535762

Luer Lock



### **Monoject 60cc Syringes**

Bold dual graduations:  
60cc in 5cc increments  
and 2oz in 1/4oz  
increments. Autoclavable  
in rigid pack for use in  
O.R/procedure packs  
oriented, compact  
package Sterile

**20/bx, 100/cs**  
KEN881560224  
Regular Luer  
KEN881560141  
Catheter (Irrigation)  
KEN881560125  
Luer Lock  
KEN881560265  
Toomey Type (Irrigation)



### **Monoject 3cc Safety Syringes**

Safety shield with  
patented two-step locking  
mechanism True transport  
position restrictive end  
collar. Bold graduations:  
0.1cc Ultra-sharp, tri-  
beveled, anti-coring  
needle Luer lock syringe  
tip. Translucent color-  
coded needle hub and  
rigid pack cap oriented,  
compact package sterile  
**100/bx, 500/cs**  
KEN881533510  
25x5/8", Red

## **Terumo**



### **Hypodermic Syringes**

Clear syringe barrel Bold,  
precise scale markings  
Translucent needle hub  
for quick observation of  
flashback. Packaging is  
designed for safety and  
cost efficiency. Plastic  
blister package is tamper  
evident, for assurance of  
sterility. Peel-packs  
provide waste disposal  
cost savings versus the  
bulky plastic tube  
containers of some other  
syringes

### **3cc Syringes**

Graduation is 0.1cc (3cc  
total) packed 100/bx,  
10boxes/case

**Without Needles**  
TERSS-03S Luer Slip  
Tip

**Luer Lock Tip with  
needles, Ultra Thin Wall**  
TERSS-03L2025  
20G x 1"

TERSS-03L2038

20G x 1-1/2"

TERSS-03L2138

21G x 1-1/2"

TERSS-03L2338

23G x 1-1/2"

TERSS-03L2516

25G x 5/8"

TERSS-03L2525

25G x 1"

**Luer Slip Tip with  
Needles, Regular Wall**

TERSS-03S2025V

20G x 1"

TERSS-03S2219V

22G x 3/4"

TERSS-03S2225V

22G x 1"

TERSS-03S2516V

25G x 5/8"

### **5cc Syringes**

Graduation is 0.2c (6cc  
total) packed 100/box,  
6boxes/case

**Without Needles**  
TERSS-05S Luer Slip  
Tip

**Luer Lock Tip with  
needles, Ultra Thin Wall**

TERSS-05L2025

20G x 1"

TERSS-05L2038

20G x 1-1/2"

TERSS-05L2125

21G x 1"

TERSS-05L2138

21G x 1-1/2"

TERSS-05L2225

22G x 1"

TERSS-05L2238

22G x 1-1/2"

### **10cc Syringes**

Graduation is 0.2cc (12cc  
total). Packed 100/box,  
6boxes/case

**Without Needles**  
TERSS-10S  
Luer Slip Tip

**Luer Lock Tip with  
needles, Ultra Thin Wall**

TERSS-10L2025

20G x 1"

TERSS-10L2038

20G x 1-1/2"

TERSS-10L2125

21G x 1"

TERSS-10L2138

21G x 1-1/2"

TERSS-10L2238

22G x 1-1/2"

### 20cc syringes

Graduation is 1.0cc (25cc total) packed 25/box, 8boxes/case

#### Without Needles

TERSS-20S

Luer Slip Tip

### 30cc Syringes

Graduation is 1.0cc (35cc total) packed 25/box, 8boxes/case

#### Without Needles

TERSS-30L

Luer Lock Tip

TERSS-30S

Luer Slip Tip

### 60cc Syringes

Graduation is 1.0cc (60cc total) packed 25/box, 4boxes/case

#### Without Needles

TERSS-60L

Luer lock Tip

TERSS-60ES

Ecc. Luer Tip

TERSS-60C

2oz cath Tip

### Allergy Syringes

TERUMO Allergy

Syringes are packaged in trays of 25 pieces and have low dead space for maximum cost savings the barrel is clear and easy to read Sharp needles minimize patient discomfort. Graduation is 0.01cc (1cc total)

TERSS-01A2713

1cc single Syringe 27G x 1/2"

### Allergy Syringe trays

1cc, packed 25 syringes/tray, 40 trays/case. Graduation is 0.01cc (1cc total)

TERSS01A2609T

26G x 3/8" Reg Bevel

TERSS01A2713T

27G x 1/2" Reg Bevel

TERSS10A2709IDT

27G x 3/8" Intradermal

Bevel needle for skin testing

TERSS10A2713IDT

27G x 1/2" Intradermal

Bevel needle for Skin Testing



### Hypodermic Needles

Terumo needles are known worldwide for extreme sharpness. They are ground and honed with a unique double-bevel design to provide a comfortable injection. Thin wall and Ultra-thin wall needles provide more comfort with an increased flow rate. With translucent hub, regular bevel packed 100/bx, 10bx/cs

TERNN-1825N

18G x 1", TW

TERNN-1925RN

19G x 1", TW

TERNN-1938RN

19G x 1-1/2", TW

TERNN-2025RN

20G x 1", UTW

TERNN-2038RN

20G x 1-1/2", UTW

TERNN-2125

21G x 1", UTW

TERNN-2138RN

21G x 1-1/2", UTW

TERNN-2225RN

22G x 1", UTW

TERNN-2338RN

23G x 1-1/2", UTW

TERNN-2613RN

26G x 1/2", RW, Low

Dead Space Hub

TERNN2713R

27G x 1/2", RW

TERNN3013R

30G x 1/2"

### Spinal Needles

Terumo spinal needles have smooth cutting edges, unbendable, unbreakable needle, and a finger-fitting grip-hub for quick entry. They are blister-packed 20 per box in peel-apart packages

TERSN-1890

18gx3-1/2"

TERSN-2090

20G x 3-1/2"

TERSN-2290

22G x 3-1/2"

TERSN-2590

25G x 3-1/2"



### Insulin Syringes

Low Dead Space design, for use with U-100 insulin only, with attached needle. Based on extreme needle sharpness



and smooth plunger movement, TERUMO insulin syringes are designed for maximum comfort. Clear, bold scale markings assure dosage control and accuracy. Latex-free. Packaged 100/box, 5 boxes/case

TERSS05M2713  
27G x 1/2", 1/2cc  
TERSS05M2813  
28G x 1/2", 1/2cc  
TERSS05M2913  
29G x 1/2", 1/2cc  
TERSS10M2713  
27G x 1/2", 1cc  
TERSS10M2813  
28G x 1/2", 1cc  
TERSS10M2913  
29G x 1/2", 1cc

### **Tuberculin Syringes**

TERUMO TB Syringes have clear, bold scale markings, which assure dosage control and accuracy. The package is tamper evident, for greater protection of the sterile contents. Latex-free. 100/bx, 10 boxes/case  
TERSS-01T  
1cc, without Needle  
TERSS01A2713T  
1cc, with 27G x 1/2"  
Regular bevel

### **SurGuard Needle Protection System**

Unique needle locking mechanism to help reduce the risk of needle stick injuries One-handed

activation meets key criteria for safety and ease of use. Available with or without attached needle. Compatible with luer lock and luer slip syringes that conform to ISO 594, ISO/ANSI fittings Cost-effective, easy to use positive lock design Latex-free. 100/box, 8 boxes/case

TERSG1825 SurGuard Sheath with 18G x 1" Needle  
TERSG1838 SurGuard Sheath with 18G x 1-1/2" Needle  
TERSG1925 SurGuard Sheath with 19G x 1" Needle  
TERSG1938 SurGuard Sheath 19G x 1-1/2" Needle  
TERSG2025 SurGuard Sheath with 20G x 1" Needle  
TERSG2038 SurGuard Sheath with 20G x 1-1/2" Needle  
TERSG2125 SurGuard Sheath with 21G x 1" Needle  
TERSG2138 SurGuard Sheath with 21G x 1-1/2" Needle  
TERSG2225 SurGuard Sheath with 22G x 1" Needle  
TERSG2238 SurGuard Sheath with 22G x 1-1/2" Needle  
TERSG2325 SurGuard Sheath with 23G x 1" Needle

TERSG2516 SurGuard Sheath with 25G x 5/8" Needle  
TERSG2525 SurGuard Sheath with 25G x 1" Needle  
TER129331 SurGuard Short Sheath used with 5/8", 3/4", 7/8" and 1" needles  
TER129349 SurGuard Long Sheath used with 1" - 1 1/2" needles

## **BARD**



### **Bulb Syringe**

Bulb syringe has a pliable bulb for responsive, controllable suction. 50cc barrel has 1cc graduations for accuracy. Finger flanges provide stable grip. Single use, sterile

DAV0035280 50cc Bulb irrigation Syringe, 50/cs



### **Toomey Syringe**

70cc piston syringe barrel fits all standard cystoscopes or resectoscopes. Double-flanged plunger tip ensures adequate vacuum. Thumb-ring plunger and smooth action permit easy one hand use.

Catheter adapter tip allows irrigation of Foley catheters. Luer adapter tip allows inflation of Foley Balloon. Single use, sterile

DAV0038480 Toomey Syringe with Catheter Tip and Luer Adapter, 50/cs



### Ear/Ulcer syringes

Available in reusable and single-use designs, each in 1,2, and 3oz. sizes Flat bottoms allow ear syringe to stand up without tipping. Slimmer tips on single-use bulbs, broader tips on reusable bulbs. 2 oz. Single use size has added suction power for efficient aspiration for newborns. Sterile

0035810 1oz,

Green, 50/cs

0035820 2oz, blue, 50/cs

0035830 3oz, green, 50/cs

## Premium Plastics



### Bulb Syringe w/ Tip Protector

Direct the flow of liquids quickly and easily with promotes visualization of fluids and measuring ease. A pliable, easy-to-grip bulb simplifies fluid retrieval/administration. Sterile bulb syringes w/ tip protector are graduated in 5cc (ml) and 5 ounce increments up to 60cc's Peel pouch

PP4090 50/cs



### Ear/Ulcer Syringe

Soft pliable material makes these syringes safe and easy to use. Whether irrigating or suctioning, the smooth tapered tip facilitates fluid administration/retrieval and reduces trauma. The ribbed bulb provides a nonslip grip. Ear/ulcer syringes are available in 2 or 3 oz sizes. Peel pouch

PP4172 2oz, sterile, 50/cs

PP4173 3oz, sterile, 50/cs

## Miltex



### Metal Ear Syringe

Metal Ear Syringe, with one short bulbous, one tapered tip, with shield, 4 ounce (12 liter)

MIL19-385 Metal Ear Syringe

### Vantage Metal Ear Syringe

MILV919-383

3-ounce capacity

MILV919-385

4-ounce capacity

## Miscellaneous

### Graham-Field

#### Fade-A-Dyne

Clear, water-soluble liquid safely and effectively removes Betadine and other iodine-based stains from tile, vinyl, upholstery, fabric, carpet, bedding and clothing. Ideal for hospitals, nursing homes, outpatient clinics, kidney dialysis, medical and veterinary offices.

Principle component: 1% Butyl Cello-solve Isopropanol. Non-flammable

GF3379 16oz, 12/cs

#### Latex Finger Cots

Not for exam use

GF3908L Large,  
144/bx  
GF3908M Medium,  
144/bx  
GF3908S Small,  
144/bx

## Dynarex Corporation



### Adhesive Tape Remover Pads

Aids in the removal of adhesive tapes Leaves the skin free of tape residue. Non-irritating, reduces patient discomfort

DYN1505 100/bx,  
10bx/cs



### Ammonia Ampules

Liquid products packaged in economical ampule form for significant savings.

- Precise measure for uniformity
- Hermetically sealed
- Unit dose
- Increased storage shelf life
- Effective delivery system

DYN1401 Ammonia  
Inhalants, 0.33cc, 10/bx



### Alcohol Swabsticks

Each swabstick is saturated with 70% alcohol. Convenient and clean application

DYNP-S84925 3s,  
25pk/bx, 10bx/cs

### Tincture of Benzoin Swabsticks

Contains compound Benzoin Tincture U.S.P. containing 77% alcohol, assures good adhesion for tapes. Protects against skin breakdown

P-S42450 1s, 50  
pk/bx, 10bx/cs



### PVP Prep Pads

Ideal for antiseptic skin preparation venipuncture, IV starts, renal dialysis, prep-op prepping, and other minor invasive procedures. Each pad is saturated with a 10% povidone-iodine solution

DYN1108 Medium,  
100/bx, 10bx/cs

### Acetone Alcohol Prep Pads

Alcohol prep pad with 10% acetone for easy skin defatting prior to injection or venipuncture

DYNP-B35900  
100/bx, 10bx/cs

## Triad Medical Products

### PVP Prep Pads

Non-irritating, non-staining, fast acting, anti-microbial solution containing 10% weight by volume povidone iodine Acts as an antiseptic germicide. Ideal for surgical site preparation

TRI10-3201 Prep Pads

### Nail Polish Remover Pads

Non-woven pads saturated with three different solvents, an emollient, a conditioning agent and a pleasant fragrance, individually packaged in sealed foil. Enhances patient monitoring by permitting observation of nail beds to detect circulation problems and loss of blood oxygen

TRI10-6200 100/bx,  
10bx/cs

### **Adhesive tape Remover Pads**

Non-woven pads saturated with an adhesive tape remover liquid, individually packaged in a hermetically sealed foiled laminated pouch. Effectively and painlessly removes adhesive tape residue without the need for harsh chemicals. Utilized throughout the hospital, in therapy, emergency rooms, orthopedic minor surgery and outpatient clinics. 100/bx, 10 boxes/case

TRI10-6300 Tape  
Remover Pads

### **Alcohol Prep Pads & Swabsticks**

Primarily used as a topical antiseptic, anti-infective prior to administering injections. Each soft, absorbent, non-woven pad or rayon-tipped swabstick is impregnated with 70% Isopropyl Alcohol, USP our sterile alcohol pads are gamma sterilized to meet AAMA/ANSI/ISO 11137 standards. Swabsticks are packaged three swabsticks per pouch, 25 pouches/box. Stick length: 4 inches

### **Alcohol Prep Pads**

TRI10-3001B  
Medium, 200/bx, 10bx/cs  
TRI10-3002  
Large 100/bx, 10bx/cs  
**Alcohol Swabsticks**  
TRI10-4300  
3/pouch, 25 pouches/bx,  
10bx/cs

### **PVP-I Prep Swabsticks**

Each rayon-tipped swabstick is saturated with a 10% PVP-I, USP solution equivalent to 1% available iodine. Our swabsticks not only kill microbes on contact, they are sturdy tools for removing small particles from wounds, such as glass fragments and gravel. Packaged in singles and triples for your convenience. Swabstick size: 4 inches

TRI10-4402 Triples, 25  
pouches/bx, 10bx/cs

### **Lemon Glycerin Swabsticks**

Lemon Glycerin swabsticks rely on tangy citrus flavor to stimulate the patient's taste buds and excite the salivary glands. And for added effectiveness, we've supplemented these swabsticks with a second water-producing catalyst-glycerin-to stimulate the salivary glands. Effective for patients with parched lips; patients on oxygen or fluid restrictions;

comatose patients on IV's; and part of routine oral care in elderly patients. Conveniently packaged in singles and triples. Swabstick size: 4 inches

TRI10-4001 Singles,  
50/bx, 10bx/cs  
TRI10-4002 Triples,  
25/bx, 10bx/cs

### **PVP-I Prep Solution**

Contains 10% weight by volume Povidone iodine, USP solution equivalent to 1% available iodine

TRI10-8204 4oz Bottle,  
72/cs

**Tyco**  
Healthcare

**Kendall**




### **T.E.D. Knee Length Anti-Embolism Stockings**

Designed to reduce venous stasis below the knee when thigh length styles are medically contraindicated. Correct graduated compression helps ensure optimum blood flow velocity. Inlay circumferential knit helps assure a proper pressure pattern and helps to

prevent slippage. Contoured seamless construction eliminates harsh ridges, which might cause constriction. Specially designed knee band minimizes slippage without constricting circulation. Wide range of sizes facilitates proper fit for correct compression pattern. Inspection toe opening provides easy access for medical examination

KEN7115 Medium, Regular, White, 12 Prs/ctn  
 KEN7203 Large, Regular, Blue, 12prs/ctn  
 KEN7604 X-Large, Regular, Green, 12prs/ctn




**T.E.D. Thigh Length Anti-Embolism Stockings**

Designed to reduce venous stasis in the entire leg. Correct graduated compression from ankle to upper thigh helps ensure optimum blood flow velocity. Inlay circumferential knit helps assure a proper pressure pattern and helps to prevent slippage. Contoured seamless construction promotes comfortable fit without the constriction caused by harsh seams and ridges. Patented interrupted band and functional gusset

minimize stocking slippage without constricting circulation. Pressure break protects critical popliteal area and promotes free leg movement. Distinctive stitch pattern in knee area assists in correct positioning. Defined heel pocket facilitates application and assures correct stocking placement. Inspection toe opening provides easy access for medical examination

KEN3416 Medium, Regular, White, 6prs/ctn  
 KEN3728 Large, Regular, Blue, 6prs/ctn  
 KEN3856 Large, Long, Blue, 6prs/ctn

**SAGE Products, Inc**




**Toothette**  
 Designed for patients who can perform all or most of their oral care, but who may also require some nursing assistance. Individually wrapped. Features soft, secure foam heads, with distinct ridges that:

- Gently lift and remove debris and mucus from the oral cavity
- Stimulate oral tissue
- Help clean between teeth
- Mini-flavored dentifrice or untreated

SAG5602 250ea/ic, 4ic/cs

**Mentor**




**EvaCare Pessary Products**

Designed to help manage uniquely feminine conditions using minimally invasive techniques and products, EvaCare pessaries are 100% silicone. They are used in the management of vaginal prolapse of the urogenital organs

MNTR2.50S #3, w/ Support  
 MNTR2.75S #4, w/ Support  
 MNTR3.00S #5, w/ Support  
 MNTR3.25S #6, w/ Support

**BARD**



## Feeding Tubes

Smooth plastic two-eyed tube used for nasogastric feeding of neonates, infants and adults as well as diagnostic and therapeutic aspiration. Transparent tube for easy observation of fluids. Flexible connector mates with standard luer-tip syringes. Attached cap protects lumen from contamination. X-ray opaque stripe. Single-use, sterile.

**15 in.w/1cm depth markings, 200/cs**

DAVO036400

5 Fr. Premature Infant

DAVO036410

8 Fr. Infant

## ENEMAS

### Allegiance

#### Flexible enema Administration kit

Our flexible enema administration kit contains a 1500cc plastic bag with a preattached tube, side clamp, liquid castile soap and a water-resistant underpad. It also includes a lubricant package, which you can use to lubricate the tube. Kits are individually packed. Caution: federal law (USA) restricts this device to sale by or on the order of a physician. 50 kits per case

BGH6010B Flexible  
Enema Administration kit

## Fleet

Pharmaceuticals

#### Fleet Adult Enema

A ready-to-use saline laxative with pre lubricated Comfortip. For use as a bowel cleanser or the relief of occasional constipation in adults and children 12 years and order

FLE201 4.5oz,  
48/cs

#### Fleet Mineral Oil Enema

A sodium-free stool softener and lubricant laxative with pre-lubricated Comfortip. For relief of fecal impaction or occasional constipation in adults and children 2 years and older

FLE301 4.5oz,  
48/cs

#### Fleet prep Kit 1 (The Dry 1)

In a clinical study comparing Fleet Prep Kit 1 vs. PEG lavage in CT Colonography, Fleet Prep Kit 1 was proven effective in leaving less residual fluid in the colon than PEG lavage. A clean and dry colon is crucial

for accurate results in radiological exams

FLE3-01320-21101  
48/cs

## Chester Labs



#### Gent-L-Tip Enemas

The Gent-L-Tip Phosphate Enema produced by Chester Labs in a registered process that is not commonly used with this type of product. After years of research we have been able to create a one-piece tube that allows our tip to be 7/8" longer than the leading national brand. This process makes a superior product; the extra benefit is the cost savings our special tube provides. Softer tubes with a lubricated tip mean that our enema really is a Gent-L-Tip

CHE101503 4.5fl. Oz.  
Phosphate Enema,  
individually Boxed, 48/cs

## RAZORS

### Personna Medical

## Razor Blades

### Prep

PRS74-0002

Prep Razors, Blade Only,  
and 100/cs

## Allegiance



### Shave Prep Razors

Our disposable razors are easy to hold and give a close shave. 100 razors/box

BGH4415B Double-sided prep razor w/ comb, 500/cs

BGH4417GN Gillette Good News twin blade razor, 100/cs

## Johnson & Johnson Medical



### PREPTIC Shave Prep Razors

PREPTIC Shave Razors assure a smooth effective shave with a precise stainless steel Persona razor blade. "Sure-Grip" razor handle feels like a reusable razor

JJ9888 Double Edge,  
Disposable razors

## Premium Plastics

### Prep Razor w/ Stainless Steel Blade

Gent-L-Kare razors are designed to provide close, smooth prep shaves. Ideal for surgical skin preps, these sharp double-edged razors help reduce skin irritation, Sonically welded blade guards stabilize/secure blades, and the nonslip grip allows for maximum control and guidance. While shaving the surgical area, weighted razors provide additional support and stability. Gent-L-Kare razors can be purchased separately or in convenient 2-compartment trays

PP4776-02 Weighted, 100/cs

PP4777-02 Un-weighted, 100/cs

## BLOOD PRESSURE MONITORING

### Baumanometer



### Baumanometer Wall Unit 33

Simple in design, the Baumanometer wall unit 33 is convenient to use, functions on the immutable law of gravity and is guaranteed to remain scientifically accurate. A 300mm Hg swivel-type instrument for wall mounting which includes a Calibrated V-Lok Inflation System, Wall Basket or Wall Clip and 8 feet of Extendex Tubing Available with nylon coated or chrome plated steel basket or stainless steel clip

BAU0850 adult Calibrated V-Lok, Wall Basket (Nylon Coated Steel)

BAU0820 Adult Calibrated V-Lok, Basket, Wall Clip (St. Steel)



### **Baumanometer Standby Sphygmomanometer**

The Baumanometer Standby is a precisely calibrated, clinical sphygmomanometer that functions on the immutable law of gravity. It is simple in design, convenient to use and is guaranteed to remain scientifically accurate. The Standby model has been providing accurate readings in hospitals, nursing homes, clinics and physicians offices since its introduction in 1938

BAU0250 Adult Calibrated V-Lok



### **Baumanometer desk Sphygmomanometer**

The Baumanometer Desk Model features a large bore, Mylar Clad glass cartridge tube, and our Calibrated V-Lok Inflation System and diecast aluminum alloy case. It is a precisely calibrated clinical sphygmomanometer, designed for medical practitioners who require a strong, lightweight,

portable instrument that offers scientific accuracy and years of trouble-free service. Like all Baumanometer instruments, it is backed by our manufacturer's guarantee. The case is made of lightweight aluminum alloy. The instrument color is Neutral Gray textured baked enamel. The aluminum scale has large white epoxy numerals on a black anodized background

BAU0320 300mm Hg, Adult Calibrated V-Lok

### **Replacement & Repair Parts**

#### **Bulbs and Valves**

This pressure source for all blood pressure instruments consists of the Baumanometer latex and Latex-free Bulb and Air-Flo Control valve

BAU1890 Latex Bulb and Air-Flo Control valve

BAU1891 Latex Bulb only

BAU1892 Air-Flo Control valve only

BAU1893 Large size latex Bulb only

BAU1894 Large Size Latex Bulb and Air-Flo Control Valve

#### **Inflation Bags**

Dipped Latex, seamless construction

BAU1840 Latex Bag (12cm x 23cm), Adult  
BAU1845 Latex bag (15cm x 33cm), Large Arm

#### **Connectors**

BAU2920 Tubing Connectors, Set

BAU2921 Tubing Connectors, Male Only

BAU2922 Tubing Connectors, Female Only

#### **Latex Tubing**

All Baumanometer latex Tubing is ozone-resistant and is supplied as permanently coiled, retractile Extendex tubing 4/32 I.D. x 9/32 O.D. In ordering tubing for various models of the Baumanometer follow this guide: Standby, 4ft Extendex, Panel Wall, Wall Unit 33: 8ft.

Extendex

BAU2910 Extendex Tubing (8 ft) with Connectors

BAU2919 Extendex Tubing (94ft) with Connectors

#### **Inflation system wall Basket**

Nylon coated steel storage basket, 6-1/2"W x 7"H x 3-3/4"D

BAU2420 Storage Basket

#### **Calibrated V-Lok Cuff And Inflation Bag**

The Calibrated system immediately and clearly indicates whether the



correct size cuff is being used. Genuine Velcro hook and loop fasteners are used to provide the optimal number of open-close cycles. The center of the bladder is clearly marked on the inner and outer surface of the cuff for correct application on either right or left limb. The cuffs are made of non-absorbent, stain resistant fabric. Sizes are specified by name; made to fit dimensions of corresponding name-size Latex inflation bags for limb circumferences shown

**Cuff and inflation Bag**

BAU1880 Adult  
BAU1869 large Arm

**Cuff Only**

BAU1870 Adult  
925cm to 35cm)  
BAU1868 Large Arm  
(33cm to 47cm)

**Calibrated V-Lok Inflation systems**

The Calibrated V-Lok System immediately and clearly indicates whether the correct cuff is being used. Genuine Velcro hook and loop fasteners are used to provide the optimal number of open-close cycles. The center of the bladder is clearly marked on the inner and outer surface for correct application on either the right or left limb. The cuff material is a tightly

woven, urethane coated, heat-set Dacron Polyester fabric. It's stain-resistant and dimensionally stable

BAU1820 Adult  
calibrated V-Lok  
BAU1821  
Child/Small adult  
Calibrated V-Lok  
BAU1822 infant  
Calibrated V-Lok  
BAU1824 Thigh  
Calibrated V-Lok  
BAU1825 Large Arm  
Calibrated V-Lok

**OMRON**



**IntelliSense Digital Blood Pressure Monitor**

Designed for professional use in a clinical setting. IntelliSense technology offers automated cuff inflation level settings and electronic deflation control. Average mode measures up to three readings and averages the total for a more accurate diagnosis, operates quietly

OMRHEM-907  
IntelliSense Digital Blood Pressure Monitor



**Exactus II Series Latex Free**

**Sphygmomanometer**

Bold Gradation with contrasting faceplate and nylon cuff. Black, nylon cuff for easy cleaning is offered in different sizes and provides accurate blood pressure reading. Deluxe, no stop pin, 300mmHg gauge and air release valve. LIFETIME gauge calibration warranty +/- 3mmHg

OMR11-0601 Adult  
OMR11-601C Child  
OMR11-601X Large Adult



**Exactus Series Cotton Cuff**

**Sphygmomanometer**

Long life cotton cuff and gauge with bold, black gradations on white faceplate. No stop pin, 300mmHg gauge. Deluxe air release valve, and vinyl carrying case is great for easy storage. LIFETIME gauge calibration warranty +/-3 mmHg

OMR115MO Large  
Adult

### **Blood Pressure Gauge**

For use with Standard  
Aneroid Units  
OMR148

## **ADC**



### **Diagnostix Deluxe Aneroid Sphyg**

Offering more features  
than any other Sphyg in  
its class, the 720 provides  
the finest value. Precision  
crafted 300mmHg no-pin  
stop manometer with  
classic black baked  
enamel finish withstands  
the toughest use. Comes  
complete with leatherette  
zippered carrying case

AD720            Adult  
AD720X        Large  
Adult

## **Heine**



### **Memotronic PC Digital Wrist Sphygmomanometer**

The easiest way to  
measure pulse and blood

pressure Automatic single  
button operation Inflation  
and deflation is fully  
automatic. Large LCD  
display

HEIM-00.09.150  
Memotronic PC with case

## **Graham-Field**

### **Patricia Sphygmomanometer**

Black enamel manometer,  
which reads pressures  
from 20mmHg to  
300mmHg, is color  
coordinated to the black  
cotton cuff. The cuff has  
artery indicator label and  
gauge holder for  
maximum convenience.  
Two tube latex bladder  
and standard chrome-  
plated air release and end  
valves. Comes complete  
with zippered leatherette  
Case and manual  
Lifetime Calibration  
Warranty

GF180X        Large  
Adult, 1ea

### **Labstar Sphygmomanometer**

This unit offers a  
calibrated nylon  
Accumax cuff with color-  
coordinated manometer,  
which measures pressure  
from 20mmHg to  
300mmHg. Deluxe  
inflation system has a two  
tube, heavy-duty latex  
bladder and a large,  
knurled air release valve

for precision deflation.  
Deluxe end valve with  
filter screen reduces dust  
build-up. Comes  
complete with zippered  
leatherette carrying case  
and instruction manual.  
Lifetime Calibration  
Warranty

GF202C        Blue,  
Child, 1ea  
GF202X        Blue,  
Large Adult, 1ea

## **Welch Allyn**



### **Tycos Classic Style Hand Aneroid**

This easy-to-use, classic  
sphygmomanometer fits  
comfortably in either  
hand with full reading  
position flexibility.  
Weights only 14-1/2  
ounces, has a feather-  
touch pressure control  
which releases air  
smoothly and evenly and  
luer-lock connectors,  
which make cuff  
changing faster. Comes  
complete with a zippered  
leatherette travel case for  
protection

WA5098-02    W/Adult  
Cuff

WA5098-03 Gauge,  
Bulb Only  
WA5098-05 W/Child  
Cuff  
WA5098-06 W/Infant  
Cuff  
WA5098-11 W/Adult  
Cuff (Less case)



**Tycos TR-1 Hand  
Aneroid**

Certified accurate to more or less 3mmHg. Trigger deflation valve provides smooth one-finger control. One-piece cuff is more durable and easier to use. Thermally treated diaphragm maintains calibration. Double-shell ABS housing for durability. Ergonomic design works well in right or left hand Luer lock for easy cuff-size changes. Comes with polyester zipper case. Latex free for safety

WA5098-27 W/Adult  
Cuff  
WA5098-28 W/Large  
Adult Cuff



**Tycos Economy Pocket  
Aneroids**

Tycos quality and accuracy in an easily affordable system certified accurate to more or less 3mmHg. One-piece cuff is more durable and easier to use. Thermally treated diaphragm maintains calibration. Comes with vinyl zipper case. Latex free for safety

WA7050-14 W/Adult  
Cuff  
WA7050-23 W/Child  
Cuff  
WA7050-24 W/Large  
Adult Cuff  
WA7050-15 Pocket  
Gauge Only



**Tycos Classic Pocket  
Aneroids**

The original classic pocket aneroids, cuff-mounted and time tested. Certified accurate to more or less 3mmHg. Ruby jeweled movement for long life. Laser-engraved

dial face for accuracy. One-piece cuff is more durable and easier to use. Thermally treated diaphragm maintains calibration. Zinc die-cast housing for mechanical stability Recessed dial face reduces parallax. Polyester zipper case Latex free for safety

WA5090-02 W/Adult  
Cuff and Case  
WA5090-03 Gauge  
Only  
WA5090-06 W/Child  
Cuff and Case  
WA5090-14 W/Adult  
Cuff, without Case  
WA5090-37 W/Thigh  
Cuff, without Case  
WA5090-41 W/Large  
Adult Cuff and Case



**Tycos Mercurial  
Sphygmomanometer**

Certified accurate to more or less 3mmHg. Break-resistant plastic tube is recessed for added safety. 45o side-to-side swivel for clear viewing (wall model only) One-piece blood pressure cuff is

more durable and easier to use. Large wire storage basket for cuffs (wall model only) Four-foot (desk and mobile models) or eight-foot (wall model) coiled tubing provides great range. Luer lock for easy cuff-size changes Latex free for safety

**Tycos Wall Mercurial Sphygmomanometer**

WA5097-26 W/Adult Cuff, charcoal gray scale  
WA5097-11 W/Adult Cuff, orange scale

**Tycos Mobile Mercurial Sphygmomanometer**

WA5097-29 W/Adult Cuff



**Tycos 509-Series Aneroids**

Certified accurate to more or less 3mmHg. Ruby jeweled movement for long life. Large, easy-to-read dial face (6.25"/15cm). Adjustable-height stand and 360o swivel for clear viewing. Large wire

storage basket for cuffs One-piece cuff is more durable and easier to use. Four-foot coiled tubing provides great range. Thermally treated diaphragm maintains calibration. Recessed dial face reduces parallax. Luer lock for easy cuff-size changes Latex free for safety

**Tycos 509-Series Mobile Aneroid**

WA5091-41 W/Charcoal Dial Face and Durable One-Piece Adult Cuff

**Tycos 509-Series Wall Aneroid**

WA5091-38 W/Charcoal Dial Face and Durable One-Piece Adult Cuff



**Tycos Wall Aneroid With Durable One-Piece Adult Cuff**

Certified accurate to more or less 3mmHg. Ruby jeweled movement for long life. Laser-engraved dial face for accuracy. Large, easy-to-read dial face (6.25"/15cm). 40o side-to-side swivel for clear viewing convenient built-in cuff storage basket. One-piece cuff is more durable and easier

to use. Eight-foot coiled tubing provides great range. Thermally treated diaphragm maintains calibration. Recessed dial face reduces parallax. Luer lock for easy cuff-size changes Latex free for safety

WA7670-01 W/Adult Cuff

**Two-Piece Blood Pressure Cuffs and Bladders**

Available in latex-free neoprene, our traditional style blood pressure cuffs and bladders also come in various sizes and combinations

**One-Tube Cuffs and Neoprene bladders**

Used for hand aneroid sphygmomanometer

WA5082-42 Child  
WA5082-43 Adult  
WA5082-44 Large Adult

**Two-Tube Cuffs and Neoprene Bladders**

Used for pocket aneroid, wall aneroid, and mercury sphygmomanometers

WA5082-08 Infant  
WA5082-24 Child  
WA5082-25 Adult  
WA5082-26 large Adult



### Blood Pressure Inflation Systems

Our latex-free blood pressure inflation systems come complete and ready to use. Comes with a durable hook-and-loop closure. Made of heavyweight nylon. Easy to clean

#### Cuff with bladder, Bulb, and valve, 6/pk

WA5082-21 Child  
 WA5082-22 Adult  
 WA5082-23 Large Adult  
 WA5082-45 Thigh



### Blood Pressure Cuffs

Heavyweight nylon cuff is economical. Easy to clean and water-resistant. Can be used with aneroid and mercury sphygmomanometers

WA5082-02 Child (gray), 12/pk  
 WA5082-01 Adult (Gray), 12/pk  
 WA5082-16 Large adult, 6/pk  
 WA5082-11 Thigh (Gray), 6/pk

### Blood Pressure Accessories

These individually wrapped blood pressure accessories are the ultimate in convenience, economy, and safety

#### Latex-free Bladders

WA5089-02H Child Two-Tube Blood Pressure Bladder, 12/pk  
 WA5089-01H Adult Two-Tube Blood Pressure Bladder, 12/pk  
 WA5089-22H Large Adult Two-Tube Blood Pressure Bladder, 6/pk

#### Black Neoprene Adult Inflation Bladder

WA5089-03 Black Neoprene Adult Inflation Bladder

#### Coil Tubing, 12/pk

WA5089-13H Coil Tubing

#### Bulbs and Valves, 12/pk

WA5088-01H Inflation Bulb and Air-Release Valve

#### Bulbs, 12/pk

WA5086-01H Inflation Bulb



#### Air Release Valves, 12/pk

WA5087-01H Blood Pressure Air-Release Valve

#### Miscellaneous

WA5098-07 Luer Lok Set

WA5091-13 Basket for Inflation System, Small  
 WA5091-47 Basket for Inflation System, Large  
 WA5089-11 Metal Twist Fit Connector  
 WA5090-45 Blood Pressure Crystal, Shatterproof

## IV Solutions

### BBraun

#### EXCEL

Large volume parenteral solutions, available in 250ml, 500ml, and one liter sizes. The EXCEL container is biologically inert, non-PVC, and non-DEHP.

#### 5% Dextrose injection USP

MCGL5100 1000ml, 12/cs  
 MCGL5101 500ml, 24/cs  
 MCGL5102 250ml, 24/cs

#### 0.9% Sodium Chloride Injection USP

MCGL8000 1000ml, 12/cs  
 MCGL8001 500ml, 24/cs  
 MCGL8002 250ml, 24/cs

#### 0.45% Sodium Chloride Injection, USP

MCGL8020 1000ml, 12/cs

#### 5% Dextrose and 0.9% Sodium Chloride Injection USP

MCGL6100 1000ml,  
12/cs  
**5% Dextrose and .045% Sodium Chloride Injection USP**  
MCGL6120 1000ml,  
12/cs  
**Lactated Ringer's Injection USP**  
MCGL7500 1000ml,  
12/cs  
MCGL7501 500ml,  
24/cs  
MCGL7502 250ml,  
24/cs  
**5% Dextrose in Lactated Ringer's Injection**  
MCGL7510 1000ml,  
12/cs  
MCGL7511 500ml,  
24/cs  
**Sterile Water for Injection USP**  
MCGL8500 1000ml,  
12/cs  
MCGL8501-01  
500ml, 24/cs  
MCGL8502 250ml,  
24/cs  
**0.15% Potassium Chloride in 5% Dextrose and 0.45% Sodium Chloride Injection USP (120mEq K+/liter)**  
MCGL6350 1000ml,  
12/cs

### **PAB**

Partial additive bags, available in 25ml, 50ml, and 100ml fill sizes. The PAB container is biologically inert, non-PVC, and non-DEHP. Also compatible with a wide range of drugs for pharmacy admixture

Catalog order numbers vary

### **5% Dextrose Injection USP**

MCGS5104-5264  
100/150ml, 64/cs  
MCGS5104-5384  
50/100ml, 84/cs  
MCGS5104-5410  
25/100ml, 116/cs

### **0.9% Sodium Chloride Injection USP**

MCGS8004-5264  
100/150ml, 64/cs  
MCGS8004-5384  
50/100ml, 84/cs

### **Glass IV Containers**

#### **5% Dextrose Injection USP**

MCGS1101-SS  
500ml, 12/cs

#### **Sterile Water for Injection USP**

MCGS9201-SS  
500ml, 12/cs  
MCGS9201-SS  
500ml, 12/cs

#### **Evacuated Containers**

MCGS9900-10  
1000ml, 6/cs  
MCGS9901-10  
500ml, 12/cs

### **Chemo Dispensing Pins**

BBraun chemo dispensing pins are highly effective in capturing the aerosolization of cytotoxic medications during the reconstitution and dispensing process

BMICDP2000

W/Mini-Spike Connector  
0.2 micron hydrophobic air-venting filter, 50/cs

### **Irrigation and Urology Solutions and Sets**

Dual-use plastic Irrigation Containers (PIC) used for both pouring and hanging. Available with a wide variety of irrigation and urology sets. Catalog order numbers vary

### **Solutions for Irrigation in Plastic Containers**

#### **0.9% Sodium Chloride Irrigation USP**

MCGR5205-01  
2000ml, 8/cs  
MCGR5200-01  
1000ml, 16/cs  
MCGR5201-01  
500ml, 16/cs  
**Sterile Water for Irrigation USP**  
MCGR5007 4000ml,  
4/cs

#### **MCGR5005-01**

2000ml, 8/cs  
MCGR5000-01  
1000ml, 16/cs  
MCGR5001-01  
500ml, 16/cs

#### **Lactated Ringer's Irrigation**

MCGR5410-01  
1000ml, 16/cs

### **Solutions for Urologic Irrigation in Plastic Containers**

#### **0.25% Acetic Acid Irrigation USP**

MCGR6600-01  
1000ml, 16/cs  
MCGR6601-01

500ml, 16/cs

**Irrigation and Urology  
Sets for BBraun PIC  
Irrigation Containers**

Continuous Overhead  
Irrigation Sets with  
UniSpike Connection

**Irrigation Set with  
UniSpike Connector,  
0.188" ID tubing, drip  
chamber, roller clamp  
with 4-1/2" distal latex  
free connector**

MCGV4508-20

Length 19", 20/cs

**Post-operative Irrigation  
Set with UniSpike, 0.188"  
ID tubing, Roberts clamp,  
4-1/2" distal latex-free  
connector**

MCGV4500 Length  
80"

**Two UniSpike  
Connectors, manifold,  
0.188" ID tubing, two drip  
chamber, roller clamp.  
Roberts clamps, 4-1/2"  
distal latex free connector  
(Consult Customer  
Support for Availability)**  
MCGV4538 Length  
92"

## BAXTER

**Evacuated Container**

Glass evacuated  
container, 6/cs

BGT1A8503 500ml

BGT1A8504 1000ml



**I.V. Solutions**

Dextrose 5% and 0.2%  
Sodium Chloride  
Injection, USP 500ml  
Viaflex  
BGT2B1093Q  
24/cs

Dextrose 5% and 0.45%  
Sodium Chloride  
Injection, USP 1000ml  
Viaflex  
BGT2B1074 12/cs

Dextrose 5% and .045%  
Sodium Chloride  
Injection, USP 250ml  
Viaflex  
BGT2B1072Q  
36/cs

Dextrose 5% and 0.45%  
Sodium Chloride  
Injection, USP 500ml  
Viaflex  
BGT2B1073Q  
24/cs

Dextrose 5% and 0.9%  
Sodium Chloride  
Injection, USP 1000ml  
Viaflex  
BGT2B1064  
12/cs

Dextrose 5% and 0.9%  
Sodium Chloride  
Injection USP 250ml  
Viaflex  
BGT2B1062Q  
36/cs

Dextrose 5% and 0.9%  
Sodium Chloride  
Injection USP 500ml  
Viaflex  
BGT2B1063Q  
24/cs

Dextrose 5% and  
Lactated Ringer's  
Injection, USP 1000ml  
Viaflex  
BGT2B2074  
12/cs

Dextrose 5% and  
Lactated Ringer's  
Injection, USP 500ml  
Viaflex  
BGT2B2073Q  
24/cs

Dextrose 5% and  
Ringer's Injection, USP  
1000ml Viaflex  
BGT2B2064  
12/cs

Dextrose 5% Injection,  
USP 1000ml Viaflex  
BGT2B0064  
12/cs

Dextrose 5% Injection,  
USP 100ml Viaflex  
BGT2B0082  
Mini-Bag Quad Pack  
4/pkg, 96/cs  
BGT2B0087

Mini-bag Single Pack 96/cs	Sodium Chloride 0.9% Injection, USP 100ml Viaflex BGT2B1302 Quad Pack, 4/pk, 96/cs BGT2B1307 Single Pack, 96/cs	BGT2B7124 12/cs
Dextrose 5% Injection, USP 150ml Viaflex BGT2B0061 36/cs		Sodium Chloride 0.9% Irrigation, USP 1500ml Pour Bottle in Plastic Containers BGT2F7125 9/cs
Dextrose 5% Injection, USP 250ml Viaflex BGT2B0062Q 36/cs	Sodium Chloride 0.9% Injection, USP 150ml Viaflex BGT2B1321 36/cs	Sodium Chloride 0.9% Irrigation, USP 250ml Pour Bottle I Plastic Containers BGT2F7122 24/cs
Dextrose 5% Injection, USP 500ml Viaflex BGT2B0063Q 24/cs	Sodium Chloride 0.9% Injection, USP 250ml Viaflex BGT2B1322Q 36/cs	Sodium Chloride 0.9% Irrigation, USP 3000ml in Arthromatic Plastic Containers BGT2B7477 4/cs
Dextrose 5% Injection, USP 50ml Viaflex BGT2B0088 Mini-Bag, Multi-Pack, 16/pk, 96/cs BGT2B0081 Mini-Bag Quad Pack, 4/pk, 96/cs BGT2B0086 Mini-Bag Single Pack, 96/cs	Sodium Chloride 0.9% Injection, USP 500ml Viaflex BGT2B1323Q 24/cs	Sodium Chloride 0.9% Irrigation, USP 3000ml Viaflex Bag in Uromatic Plastic Container BGT2B7127 4/cs
Sodium Chloride 0.45% Injection, USP 1000ml Viaflex BGT2B1314 12/cs	Sodium Chloride 0.9% Injection, USP 50ml Viaflex BGT2B1308 Multi- Pack, 16/pk, 96/cs BGT2B1301 Quad Pack, 4/pk, 96/cs BGT2B1306 Single Pack, 96/cs	Sodium Chloride 0.9% Irrigation, USP 500ml Pour Bottle in Plastic Containers BGT2F7123 18/cs
Sodium Chloride 0.45% Injection, USP 500ml Viaflex BGT2B1313Q 24/cs		Water Sterile Injection USP 1000ml Viaflex BGT2B0304 12/cs
Sodium Chloride 0.9% Injection, USP 1000ml Viaflex BGT2B1324 12/cs	Sodium Chloride 0.9% Irrigation, USP 1000ml Pour Bottle in Plastic Containers BGT2F7124 12/cs	Water Sterile Irrigation, USP 1000ml Pour Bottle in Plastic Containers BGT2F7114 12/cs
Sodium Chloride 0.9% Injection, USP 100ml Viaflex Multi-pack 16/pk BGT2B1309 96/cs	Sodium Chloride 0.9% Irrigation, USP 1000ml Viaflex Bag in Uromatic Plastic Container	Water Sterile Irrigation, USP 1000ml Viaflex bag in Uromatic Plastic Container BGT2B7114 12/cs



Water Sterile Irrigation,  
USP 1500ml Pour Bottle  
in Plastic Containers  
BGT2F7115 9/cs

Water Sterile Irrigation,  
USP 2000ml Viaflex Bag  
in Uromatic Plastic  
Container  
BGT2B7116 6/cs

Water Sterile Irrigation,  
USP 250ml Pour Bottle  
in Plastic Containers  
BGT2F7112 24/cs

Water Sterile Irrigation,  
USP 3000ml Viaflex bag  
in Uromatic Plastic  
Container  
BGT2B7117 4/cs

Water Sterile Irrigation,  
USP 500ml Pour Bottle  
in Plastic Containers  
BGT2F7113 18/cs

## **IV ADMIN. SETS**

### **BBraun**

#### **Specialty I.V. Sets**

BBraun offers a full line of specially administration and extension sets including, blood sets, 0.2 and 1.2-micron filter sets, and Rate Flow Regulator gravity flow control sets etc.

**Basic Administration Sets,  
15 drops/ml Controll  
clamp, slide clamp, two-  
piece make luer lock**  
MCGV1400 Priming  
volume: 12,1, length: 67"

MCGV1495 Priming  
volume: 16ml, Length:  
89"

**Basic Administration Sets  
with One Injection Site 15  
drops/ml Controll clamp,  
slide clamp, injection site  
28" above distal end, male  
luer lock**

MCGV1415 Priming  
Volume: 16ml, Length:  
93"

**Primary ADDitlv  
Administration sets with  
Check Valve, 15 drops/ml  
Check valve, Controll  
clamp, and male luer lock**  
MCGV1485 Injection  
sites 6" and 80" above  
distal end. Priming  
volume: 19,1, Length  
108"

MCBV1447 Injection  
sites 6", 28" and 90"  
above distal end. Slide  
clamp, Priming volume  
20ml, and length: 115"

#### **Peripheral Extension Sets, Latex-Free**

**Standard Bore  
IV Administration  
Extension Set, 2 injection  
sites, slide clamp, female  
adapter, male luer slip.  
Priming Vol. 6ml**  
MCGV5404 length 33",  
50/cs  
**IV Administration  
Extension Set, without Y-  
site, female adapter, male**

**luer slip. Priming Vol.  
3ml.**

MCGV5400 Length  
21", 50/cs

**IV Administration  
Extension Set, without Y-  
site, female adapter, male  
luer slip. Priming Vol.  
5ml.**

MCGV5401 length 33",  
50/cs

#### **Small Bore**

**ET06 Extension Set, small  
bore, luer lock connections  
with slide clamp, priming  
Vol. 0.2ml**

MCG471960 Length 6",  
100/cs

#### **WHIN Winged Huber Infusion Needle Sets**

Whether you are delivering chemotherapy, antibiotic therapy, or parenteral nutrition, BBraun has the WHIN Sets to meet your needs. Our sets offer Needle-free valves, latex-free y sites, and DEHP-free tubing

**Winged Infusion Set  
Huber Needle 20G x 1"  
90o bend, 6" tubing, Y-  
site, on/off clamp**

MCG471720 20/cs

**Winged Infusion Set  
Huber Needle 20G x 1"  
90o bend, 6" tubing, on/off  
clamp**

MCG471732 20/cs

**Winged Infusion Set  
Huber Needle 20G x 3/4"  
90o bend, 6" tubing, on/off  
clamp**

MCG471731 20/cs

### Basic IV Sets

BBraun offers a full line of basic administration sets. BBraun's administration sets feature latex free injection sites for patient safety, and are available with needle free injection sited for added clinical safety. DEHP free tubing is available to enhance compatibility when using Paclitaxel or Lipid suspensions

#### **Basic administration sets, 15 drop/ml, roller clamp, slide clamp, two-piece male luer lock**

MCGV1400 Priming volume: 14ml Length: 68"

MCGV1400-20

Priming volume: 14ml Length: 68", Flashtube Connector

#### **Basic administration sets, 15 drop/ml, unvented, Controll clamp, slide clamp, one injection site, two-piece male lock**

MCGV1416 Priming volume: 14ml, Length: 73"

#### **Basic administration set, 15 drops/ml, unvented, Controll clamp, injection site, two-piece male luer lock**

MCGV1417 Priming volume: 15ml, Length: 75"

#### **Basic IV administration set, 15 drops/ml, vented, Controll clamp, slide clamp, one injection site, two-piece male luer lock**

MCGV1418 priming Volume: 15ml, Length: 76"

**Basic administration set, 15 drops/ml, non-vented Controll clamp, slide clamp, injection site, two-piece male luer lock, Sigma compatible tubing**  
MCGV1482 Priming volume: 19ml, Length: 106"

### Standard Bore Extension Sets

#### **Female adapter, male luer slip**

MCGV5400 Priming volume: 3ml, Length: 21"

MCGV5401 Priming volume: 5ml, Length: 33"

#### **Female adapter, injection sites 6" and 24" above distal end, slide clamp, male luer slip**

MCGV5404 Priming volume: 6ml, Length: 33"

## BAXTER

### I.V. Administration Sets

I.V. Administration T-Connector Extension Set, 1 Injection Site Microbore Tubing Interlink System  
BGT2N3328 200/cs

I.V. Administration Extension Set, 1 Injection Site Minivolume "T" Connector 7"  
BGT2C5681 48/cs

I.V. Administration catheter Extension Set, 1 Injection Site Microbore

Tubing 6" Interlink System  
BGT2N3370 200/cs

I.V. Administration Catheter Extension set, 1 Injection Site Microbore Tubing 6" Interlink System  
BGT2N3374 200/cs

I.V. Administration Extension Set, 2 Injection Sites 21" Interlink  
BGT2C6606 48/cs

I.V. administration Extension Set, 2 Injection Sites 43"  
BGT2C5612 48/cs

I.V. Administration Extension Set, 2 Injection Sites 43"  
BGT2C5630 48/cs

I.V. Administration Extension Set, No Injection Site 20"  
BGT2C5625 48/cs

I.V. Administration Extension Set, No Injection Site 22"  
BGT2C6224 50/cs

I.V. Administration Extension Set, No Injection Site 21"  
BGT2C6226 50/cs

I.V. Administration Extension Set, No Injection Site 30"  
BGT2C5645 48/cs

I.V. Administration  
Extension Set, No  
Injection Site 35"  
BGT2C6225 50/cs

I.V. administration  
Extension Set, No  
injection Site 34"  
BGT2C5627 48/cs

I.V. administration  
Extension Set, No  
Injection site 7"  
BGT2C5641 48/cs

I.V. Administration  
Extension Set, No  
Injection Site Large Bore  
Tubing 20"  
BGT2C5628 48/cs

I.V. Administration  
Catheter Extension Set,  
No Injection Site  
Microbore Tubing 6' IV  
Interlink System  
BGTN1194 50/cs

I.V. Administration Filter  
Extension Set, 1 Injection  
Site 0.22 Micron High  
Pressure 16"  
BGT2C5671 48/cs

I.V. Administration Set  
Primary, 1 Injection Site  
Non-Vented Pin 68"  
BGT2C5417 48/cs

I.V. Administration Set  
primary, 1 Injection Site  
non-Vented Pin 68"  
BGT2C5418 48/cs

I.V. Administration Set  
Primary, 1 Injection Site  
Non-Vented Pin 68"  
BGT2C5431 48/cs

I.V. Administration Set  
Primary, 1 Injection Site  
Non-Vented Pin 76"  
Interlink System  
BGT2C6401 48/cs

I.V. Administration Set  
Primary, 1 Injection Site  
Duo-Vented Pin 68"  
BGT2C5419 48/cs

I.V. Administration Set  
primary, 1 Injection Site  
Duo-Vented Pin 92"  
BGT2C5435 48/cs

I.V. Administration Set  
Filter primary, 2 Injection  
Sites, 0.22 Micron High  
Pressure, 105" Continu-  
Flo, 48/cs  
BGT2C6571 48/cs

**I.V. Administration Sets**  
I.V. Administration set  
primary, 2 Injection Sites  
Non-Vented Pin 97"  
BGT2C5427 48/cs

I.V. Administration Set  
primary, 2 Injection Sites  
Non-Vented Pin 97"  
BGT2C5439 48/cs

I.V. Administration Set  
Primary, 2 Injection Sites  
Non-Vented Pin 101"  
Interlink System  
BGT2C6425 48/cs

I.V. Administration Set  
Primary, 2 Injection Sites  
Non-Vented Pin 81"  
Continu-Flo  
BGT2C5519 48/cs

I.V. Administration Set  
Primary, 2 Injection Sites  
Non-Vented Pin 85"  
BGT2C5423 48/cs

I.V. Administration Set  
primary, 2 Injection Sites  
Non-Vented Pin 85"  
BGT2C5424 48/cs

I.V. Administration Set  
Primary, 2 Injection Sites  
Non-Vented Pin 89"  
Continu-Flo  
BGT2C5531 48/cs

I.V. Administration Set  
primary, 2 Injection Sites  
Non-Vented Pin 89"  
Continu-Flo Interlink  
BGT2C6519 48/cs

I.V. Administration Set  
Primary, 2 injection Sites  
Duo-Vented Pin 105"  
Continu-Flo  
BGT2C5543 48/cs

I.V. Administration Set  
Primary, 3 injection Sites  
Non Vented Pin 105"  
Continu-Flo  
BGT2C5545 48/cs

I.V. Administration Set  
Primary, 3 Injection Sites  
Non-Vented Pin 109"  
Continu-Flo  
BGT2C5527 48/cs

I.V. Administration Set Primary, 3 Injection sites Non-Vented Pin 110" Continu-Flo Interlink BGT2C6537 48/cs	Male Luer Slip Large Bore BGT2C6202 50/cs	BGT2N3395 200/cs Injection Site Long w/ Luer Lock BGT2N1198 200/cs
I.V. Administration Set primary, 3 Injection Sites Non-Vented Pin 93" Continu-Flo BGT2C5521 48/cs	I.V. Administration Stopcock 3-Way Adapter  Rotating Male Luer Lock Large Bore BGT2C6201 50/cs	Injection Site Luer Lock Interlink System BGT2N3399 200/cs  Injection Site Luer Lock Interlink System BGT2N3379 800/cs
I.V. Administration Set Secondary 18G Needle Lock (Detached) Non- Vented Pin 37" BGT2C7431 48/cs	<b>I.V. Administration Sets</b> I.V. Administration Stopcock 3-Way Adapter Rotating Male Luer Lock Standard Bore BGT2C6240 50/cs	Injection Site Standard w/ Luer Lock BGT2N1199 200/cs
I.V. Administration Set Secondary lever Lock Cannula (Detached) vented Pin 37" Interlink System BGT2C7451 48/cs	I.V. Administration Stopcock 3-way with Extension Tubing Luer Lock 1.2ml 8" BGT2C5602 48/cs	Irrigation Set Cysto 78" Overall Length BGT2C4040 48/cs  Lactated Ringer's Injection, USP 1000ml Viaflex BGT2B2324 12/cs
I.V. Administration Sets Burette 150ml Add-On Set Valve less with Slide Clamp Buretrol BGT2C7505 48/cs	I.V. Administration Stopcock 3-Way with Extension Tubing Luer Lock 3.0ml 21" BGT2C5603 48/cs	Lactated Ringer's Injection, USP 250ml Viaflex BGT2B2322Q 36/cs
I.V. Administration Sets Burette 150ml Ball Valve, 2 injection Sites 84" Buretrol BGT2C7546 48/cs	I.V. Administration Stopcock 3-Way with Extension Tubing Male Luer Slip (2 Female Luer Slips) 2.3ml 20" BGT2C6229 50/cs	Lactated Ringer's Injection, USP 500ml Viaflex BGT2B2323Q 24/cs
I.V. Administration Sets Specialty Nitroglycerin Vented Pin, 1 Injection Site 105" BGT2C7552 48/cs	I.V. Administration Stopcock 3-Way with Extension Tubing male Luer Slip (2 Female Luer Slips) 3.6ml 33" BGT2C6230 50/cs	Lactated Ringer's Irrigation, 1000ml Pour Bottle in Plastic Containers BGT2F7154 12/cs
I.V. Administration Stopcock 3-Way Adapter	Universal Vial Adapter Interlink System	Lactated Ringer's Irrigation, 3000ml in

Arthromatic Plastic Containers  
BGT2B7487 4/cs

Lactated Ringer's Irrigation, 5000ml in Arthromatic Plastic Containers  
BGT2B7489 2/cs

Procedure Tray Caudral and Epidural  
BGT1T1071 10/cs

## IV THERAPY

### BAXTER

#### I.V. Catheters

I.V. Catheter, Quick-Cath. 200/cs  
BGT2N1114 18G x 1-1/4" Radiopaque  
BGT2N1115 20G x 1-1/4" Radiopaque  
BGT2N1116 22G x 1" Radiopaque  
BGT2N1117R 24G x 5/8" Radiopaque

### Becton Dickinson



#### BD Posiflow IV Access Products

BD Posiflow eliminates blood drawback into the catheter upon luer disconnect. Compatible for use with lipids and chemotherapy agents  
Intuitive, easy-to-use and

patient-friendly design  
Low-priming volume, dead space. Straight, two-way fluid path  
Flow rate capability to meet clinical needs. Flat swabable surface. Latex free. Contains no metallic materials

BD385302 6" Luer Slip Extension Set, 50/bx, 200/cs  
BD385303 6" Bi Luer-Lok Extension Set, 50/bx, 200/cs  
BD385304 6" Bi Luer Slip Extension Set, 50/bx, 200/cs  
BD385305 8" "Y" Extension Set, 25/bx, 200/cs  
BD385306 34" Dual "Y" Extension Set, 100/cs  
BD385307 100" Dual "Y" Drip Set, 25/bx, 100/cs  
BD385308 Multi Dose Vial Adapter, 25/bx, 100/cs  
BD385309 Port Converter, 50/bx, 200/cs  
BD385323 100" "Y" Drip Set, 25/bx, 100/cs  
BD385324 6" "Y" Extension Set, 25/bx, 200/cs



#### BD InterLink Blunt Plastic Cannula

Designed for improved performance, the BD Blunt Plastic Cannula enhances the Interlink IV Access System. The BD Blunt Plastic Cannula is a truly universal cannula, and can be used with all types of split septum IV injection sites and vials designed for needleless access. It provides simplified access to the interlink injection site with smoother insertion, reduced penetration force, and greater fluid dispersion

BD303367 Vial access Cannula 100/bx, 1000/cs  
BD303381 InterLink Blood Collection Assembly 25/bx, 200/cs  
BD303400 3cc Syringe with Syringe Cannula 100/bx, 800/cs  
BD303401 3cc Syringe with Vial Access cannula 100/bx, 800/cs  
BD303402 5cc Syringe with Syringe Cannula 100/bx, 400/cs  
BD303403 5cc Syringe with Vial Access Cannula 100/bx, 400/cs  
BD303404 10cc Syringe with Syringe Cannula 100/bx, 400/cs  
BD303405 10cc Syringe with Vial Access Cannula 100/bx, 400/cs



**BD Saf-T-Intima  
Integrated Safety IV  
Catheters**

It's the Complete Safety Catheter. BD Saf-T-Intima integrated safety IV catheters combine advanced needle-stick safety with an innovative closed catheter system design. The result is a virtually bloodless insertion and greater safety for you. The BD Saf-T-Intima integrated safety IV catheter with "Y" adapter allows for pre-priming and flushing during venipuncture.

Proven to optimize patient outcomes

**BD Saf-T-Intima  
Integrated safety IV  
Catheters, 25/bx, 200/cs**

DES383312 24ga x  
0.75in (0.7 x 19mm)

Yellow

DES383322 22ga x  
0.75in (0.9 x 19mm) Blue

DES383335 20ga x  
1.00in (1.1 x 25mm) Pink

**BD Saf-T-Intima  
Integrated Safety IV  
Catheters with Y Adapter,  
25/bx, 200/cs**

DES383313 24ga x  
0.75in (0.7 x 19mm)

Yellow

DES383323 22ga x  
0.75in (0.9 x 19mm) Blue

DES383336 20ga x  
1.00in (1.1 x 25mm) Pink

DES383346 18ga x  
1.00in (1.3 x 25mm)  
Green



**BD Insyte Autoguard  
Shielded IV Catheters**

BD Insyte Autoguard Shielded IV Catheters mean added safety for you. They're the safety-engineered version of the popular BD Insyte IV Catheter. The Patented BD Autoguard push-button shielding technology is proven to effectively safeguard healthcare workers from accidental needle stick injuries

**BD Insyte Autoguard  
Shielded IV catheters,  
50/bx, 200/cs**

DES381412 24ga x  
0.75in (0.7 x 19mm)

Yellow

DES381423 22ga x  
1.00in (0.9 x 25mm) Blue

DES381433 20ga x  
1.00in (1.1 x 25mm) Pink

DES381434 20ga x  
1.16in (1.1 x 30mm) Pink

DES381437 20ga x  
1.88in (1.1 x 48mm) Pink

DES381444 18ga x  
1.16in (1.3 x 30mm)

Green

DES381447 18ga x  
1.88in (1.3 x 48mm)

Green

DES381454 16ga x  
1.16in (1.7 x 30mm) Gray

DES381457 16ga x  
1.77in (1.7 x 45mm) Gray

DES381467 14ga x  
1.75in (2.1 x 45mm)  
Orange

**BD Insyte Autoguard  
Winged Shielded IV  
Catheters, 50/bx, 200/cs**

DES381512 24ga x  
0.75in (0.7 x 19mm)

Yellow

DES381523 22ga x  
1.00in (0.9 x 15mm) Blue

DES381533 20ga x  
1.00in (1.1 x 25mm) Pink

DES381534 20ga x  
1.16in (1.1 x 30mm) Pink

DES381537 20ga x  
1.88in (1.1 x 48mm) Pink

DES381544 18ga x  
1.16in (1.3 x 30mm)

Green


DES381547 18ga x  
1.88in (1.3 x 48mm)

Green

DES381554 16ga x  
1.16in (1.7 x 30mm) Gray

DES381557 16ga x  
1.77in (1.7 x 45mm) Gray

**EXEL**



**Huber Infusion Set**

Latex Free sterile, 90-degree 316L stainless steel needle, non-coring, "Y" injection site available, 7" Microbore tubing, rotating wing Non-DEHP tubing

available Color-coded for easy identification

EXL26940R 19G x 1  
EXL26944R 20G x 1  
EXL26945R 20G x 3/4  
EXL26948R 22G x 1

## TERUMO




### **SURFLO I.V. Catheters**

An ultra sharp, double-bevel, introducing needle plus an exclusive silicone lubricant (on needle and catheter) allow for easier penetration and smoother travel through tissue. Tri-layer catheter offers less resistance plus Radiopacity. Bevel-up indicator assures proper needle bevel alignment. Flashback chamber permits fast and easy confirmation of venous access. Large easy to grip chamber thin wall

**16G Needles**  
TERSR-OX1451CA  
14G x 2"  
**18G Needles**  
TERSR-OX1632CA  
16G x 1-1/4"  
**20G Needles**  
TERSR-OX1832CA  
18G X 1-1/4"  
**22G Needles**  
TERSR-OX2025CA  
20G x 1"  
TERSR-OX2032CA

20G x 1-1/4"  
**24G Needles**  
TERSR-OX2225CA  
22g X 1"  
**27G Needles**  
TERSR-OX2419CA



### **SURFLO Winged Infusion Sets**

A precisely honed needle with ultra thin wall design assures ease of penetration and continued access with a single insertion. An exclusive, silicone lubricant provides smoother needle passage through tissue and with the needles bonded to the wings at a slight downward angle, there is less risk of vessel damage or occlusion of the lumen. Soft and flexible, the wings readily conform to the body contour

**Regular, 12" Tubing,**  
**50/bx, 10bxs/cs**  
TERSV-18BLK  
18Gx3/4", TW  
TERSV-19BLK  
19Gx3/4", TW  
TERSV-21BLK  
21Gx3/4", UTW  
TERSV-22BLK  
22Gx3/4", UTW  
TERSV-23BLK  
23Gx3/4", UTW

TERSV-25BLK  
25Gx3/4", TW  
**Regular, 3-1/2" Tubing,**  
**100/bx, 5bxs/cs**  
TERSV-19BLS  
19Gx3/4", TW  
TERSV-21BLS  
21Gx3/4", UTW  
TERSV-23BLS  
23Gx3/4", UTW  
TERSV-25BLS  
25Gx3/4", TW  
**Intermittent, 3-1/2" Tubing, 100/bx, 5bxs/cs**  
TERSV-23BI  
23Gx3/4", UTW

### **SURFLO Injection Plugs**

Large target area allows for easy needle penetration. Injection Plugs have a tamper resistant seal and easy to open packaging  
TERSR-IP2 Luer Lock,  
100/bx, 10bx/cs

## **Applicators Towelettes & Tongue Blades**

## **Birchwood Laboratories, Inc**



### **Scopettes Comfort Tip" Applicators**

Safe, highly absorbent, uniform and lint free superb quality at an affordable price The 5/8" diameter rayon swab head has a 3/8" cushioned tip

to minimize the possibility of patient trauma

**8" Scopettes, Jr. Ob/Gyn "Comfort Tip"**

**Applicators**

Rayon Swab, Paper Shaft, Non-sterile

BIR34-7021-8

500/cs

BIR34-7021-12

100/bx, 12bx/cs

**16" Scopettes "Comfort Tip" Applicators**

Rayon Procto Swab, Plastic Shaft, Non-sterile

BIR34-7023-12

100/bx, 12bx/cs

## Puritan



**Multi-Purpose Cotton-Tipped Applicators**

Cotton tipped applicators are suited to a wide range of uses in the healthcare field and in industrial settings. High-quality cotton is carded, spun, and applied to northern White Birch, plastic or paper shafts at our integrated facility to ensure consistent quality in every swab. Popular with both our medical and industrial customers, and used everywhere from the exam room to the clean room. Choose the swab right for your application

HW803-WC Wood Shaft, 3" x 1/2", Regular

tip, 100/bg, 5bg/bx, 10bx/cs

HW806-WC Wood Shaft, 6" x 1/12", Regular

Tip, 100/bg, 10bg/bx, 10bx/cs

HW25-806-1WC

Wood Shaft, 6" x 1/12", Regular Tip, Sterile, 1/pk, 100pk/bx, 10bx/cs

**Multi-Purpose Cotton-Tipped Applicators**

HW25-806-2WC

Wood Shaft, 6" x 1/12", Regular Tip, Sterile, 2/pk, 100pk/bx, 10bx/cs

HW806-WCL

Wood Shaft, 6" x 1/12", Large Tip, 100/bg, 10bg/bx, 10bx/cs



**Healthcare Cotton Tipped Applicators**

Different tip sizes and shapes, choice of shaft materials, and sterile packaging options make Puritan the applicator of choice for healthcare professionals worldwide. Many of our swabs can be ordered in non-standard lengths – for use in diagnostic kits, for example

HW826-WC Wood Shaft, 6" x 1/12", Tapered

Mini Tip, 100/bg, 10bg/bx, 10bx/cs



**Dacron Tipped Applicators, Sterile**

Puritan is the world's leading provider of diagnostic grade Dacron polyester tipped applicators. Choose from a range of sizes and configurations with wood, aluminum, or plastic shafts sterile

HW25-806-1WD Wood Shaft, 6" x 1/12", Regular Tip, 1/pk, 100pk/bx, 10bx/cs

HW25-806-2WD Wood Shaft, 6" x 1/12", Regular Tip, 2/pk, 100pk/bx, 10bx/cs

HW25-806-1PD Plastic Shaft, 6" x 1/10", Regular Tip, 1/pk, 100pk/bx, 10bx/cs

HW25-806-2PD Plastic Shaft, 6" x 1/10", Regular Tip, 2/pk, 100pk/bx, 10bx/cs



**Histobrush Cytology Products**



Meets patient comfort and clinical performance needs. 1pk, 10pk/bx, 100bx/cs  
HW2199 Tapered Soft Histobrush



### **Tongue Depressors**

Puritan Tongue depressors come in plain, wrapped flavored, and sterile styles. Customized imprinting and packaging is available for some styles. Industrial grade depressors are suitable for all non-healthcare applications

HW711 Junior, 5-1/2" x 5/8", 500/bx, 10bx/cs



HW710C Junior, 5-1/2" x 5/8", Cherry, 50/bx, 20bx/cs

HW710G junior, 5-1/2" x 5/8", Grape, 50/bx, 20bx/cs



HW704 Standard, 6" x 11/16", 500/bx, 10bx/cs



HW705 Regular, 6" x 3/4", 500/bx, 10bx/cs



### **Applicator Sticks**

Ideal for specimen transfer, for use as mixing rods, for laboratory applications and for general purposes in industrial settings Made of Northern White Birch with Straight-cut ends  
HW807 5-3/4", Thin Diameter, 1000/bx, 20bx/cs



### **Calgiswab Calcium Alginate Applicators**

Tipped with a natural organic fiber derived from seaweed. Calcium alginate contains no fatty acids, and is biodegradable in the human system. Available with wood, plastic, or aluminum shafts, and sterile or non-sterile packaging

HW25-800A Sterile Nasopharyngeal Aluminum Shaft, 5-1/2" x .035", Tightly Wound Tip

HW25-801A Sterile, Urethro-Gential Aluminum Shaft, 5-1/2" x .025", Tightly Wound Tip



### **Proctoscopic Applicators**

Rayon tipped, in a choice of shaft lengths and materials. The paper and polystyrene shaft styles can be autoclaved, and are also suitable as a standard mouth and throat swab. Autoclavable to 250oF

HW808-50 Paper Shaft, 8" x 5/32", Oversized Tip, 50/bx, 10bx/cs

## **CITMED Corporation**

### **Wood Applicators**

Sterile applicators are processed in sealed packets of one's and two's and packaged in both cartons and bulk. Sterile tongue depressors are processed in sealed packets of one and packaged in both cartons and bulk. Non-sterile applicators, packed in sealed bags are autoclavable. Bags have autoclave markings. Our 16" polypropylene swab can also be sterilized by autoclaving. Citmed makes 8" Swabs sterile and non-sterile

### **Sterile Applicators**

#### **Wood Shaft**

**Cotton Tipped, 10 cartons of 100pks/cs**

CTM22-9593 3", 2/pk

**Dacron Tipped, 10 cartons of 100pks/cs**

CTM22-9633 6", 1/pk

**Rayon Tipped, 10 cartons of 100pks/cs**

CTM22-9158 6", 1/pk

### **Nonsterile Applicators, Wood Shaft**

**Cotton Tipped, 10 cartons of 100pks/cs**

CTM22-9610

3", Regular

CTM22-9610

1/pk

**Wood Applicators, 12 cartons of 864/cs**

CTM22-9605 6"

### **Applicators, Paper Shaft**

**8" Jumbo Rayon Swabs, sterile**

CTM22-2006/S

500/cs

**8" Jumbo Rayon Swabs, Nonsterile**

CTM22-9901

100/bx, 10bx/cs

## **Graham-Field**

### **Silver Nitrate Applicators**

Wooden applicator stick tipped with 75% silver nitrate, 25% potassium nitrate. Silver nitrate tips are soluble in all secretions. Control strength by dipping in

water. Convenient dispenser tube vial  
GF1590 6",  
100/vials

## **Topical Sprays**

### **CETYLITE Industries, Inc**



#### **Cetacaine Spray**

Topical anesthetic spray  
Controls pain and gagging Effective only on mucous membrane  
Exclusive Jetco-Spray Cannula provides complete control. 56gm bottle with J-4 Jetco-Spray Cannula

CET00201 56gm  
Cetacaine Spray

## **Soaps & Lotions**

### **Decon Labs, Inc**



#### **BacDown Antimicrobial Handsoap Pump**

BacDown Antimicrobial Handsoap will not dry out your hands, even if you scrub 30,40, or 50 times a day. Triclosan (a broad spectrum Antimicrobial agent), helps to reduce skin-borne microorganisms for a cleaner, safer work environment

DEC21775 16oz.  
Pump, 6/cs

#### **BacDown Antimicrobial Handsoap, 32oz**

DEC28444 32oz,  
12/cs

#### **BacDown Antimicrobial Handsoap, Gallon**

DEC28445 Gallon,  
4/cs



**PROGUARD  
Professional  
Handcream**

PROGUARD Professional Handcream combines soothing skin emollients with a barrier film that seals in moisture and seals out problems. PROGUARD is particularly effective in improving latex glove comfort. The non-greasy formula will not cause the break down of latex gloves. PROGUARD is not perfumed and absorbs rapidly for instant protection

DEC28447 16oz. w/  
dispenser, 6/cs

**Coloplast**



**Isagel Waterless Hand  
Cleaner**

Convenient on-the-spot infection control when water isn't available. In vitro testing demonstrates that this instant hand

sanitizing gel, effective against MRSA and VRE, kills 99.5% of all microorganisms tested in just 15 seconds. Nurtures skin with moisturizers as it disinfects

SWE1644 4oz

**Chester Labs**



**AprilFresh  
Antibacterial Liquid  
Soap**

AprilFresh Antibacterial Liquid Soap contains the active ingredient TRICLOSAN which is an effective germ killer. Strong on bacteria, yet still mild on your hands

CHE000600 128oz.,  
4/cs

**AprilFresh Hand  
Sanitizer**

AprilFresh Hand sanitizer is specially formulated to clean and sanitize hands and skin. This no-rinse formula is ideal for situations where water is not readily available, or extra protection is needed. Special skin conditioners prevent skin from drying out. 70% ethyl alcohol w/conditioners  
Light fragrance

CHE011227 8oz. bottle,  
36/cs

**Huntington**

**BACTI-STAT Liquid  
Hand Soap**

Bacti-Stat provides excellent Antimicrobial activity and gently cleans the skin to avoid affecting its function in providing a natural barrier to outside intrusion of microorganisms. Bacti-Stat contains Aloe Vera and lanolin ethoxylate to help moisten and soothe skin

HNT61017285  
Gallon, 4/cs  
HNT61048722  
540ml, 12/cs  
HNT61011292  
4oz, 72/cs

**Graham-Field**

**Tincture of Green Soap,  
N.F.**

Alkaline cleansing detergent made with pure vegetable oil soap in alcohol and water solution. Distinctive Oil of Lavender scent  
Manufactured in strict accordance with latest edition of National Formulary. Ideal for topical skin cleansing where usage of other detergents containing

U.S.P. Hexachlorophene  
is not desired

GF3375      Gallon,  
4/cs

## HealthLink



### AloeGuard

#### Antimicrobial Soap

“Kills the germs, NOT your hands!” AloeGuard is specially formulated with soothing, moisturizing Aloe Vera to protect your skin’s delicate balance, and Chloroxylenol to provide the high Antimicrobial activity necessary for personal protection and reduction of cross contamination. Feel great, smells great and leaves hands feeling soft and clean

HLK7720      800ml  
bag-in-box, 12/cs

HLK7725      4oz. flip-  
top bottle, 24/cs

HLK7740      Gallon  
refill, 4/cs

HLK7760      18oz.  
pump bottle, 12/cs

HLK7763      18oz.  
empty pump bottle 48/cs

### AloeGuard wall Dispenser

HLK7730      Wall  
dispenser for 7720, 12/cs



### AloeSoothe

#### Moisturizing Lotion

The professional medical grade lotion that soothes and protects! AloeSoothe is the perfect moisturizing lotion for any healthcare professional. Because of its unique blend of ingredients, AloeSoothe moisturizes and helps seal out harsh irritants while allowing the body’s natural healing process to occur. It also helps prevent chapping and irritation caused by frequent hand washing and regular use of latex gloves unlike most commercial lotions. AloeSoothe has no petroleum products added; therefore it’s latex friendly

HLK7750      18oz.  
Pump bottle, 12/cs



### WaterSense Automatic faucet controller

WaterSense eliminates the need for staff or patients to touch potentially contaminated faucets. The unit automatically switches water on as hands move under the spout and switches off when movement stops. The computer-driven sensing system automatically controls water flow enabling individuals to wash without fear of infectious agents. Staff and patients are fully protected from becoming disease-carrying vectors, maximizing infection control in every handwashing environment

HLK7790      12/cs

## Metrex



### **VioNex Antimicrobial liquid Soap**

VioNex Antimicrobial Liquid Soap is intended for hand washing between routine patient contact, procedures, after exposure to potential pathogens and after removal of personal protective equipment. It is formulated to be gentle and non-irritating for hand washing 30-50 times per day. VioNex Antimicrobial Liquid soap is recommended for high-risk areas, critical care areas, nursing stations; washrooms and all patient care areas. 99.945 germ reduction in 30 seconds

MTX10-1500

Gallon, 4/cs

MTX10-1518

18oz. Pump Bottle, 12/cs



### **VioNex Antimicrobial Skin Wipe Towelette**

VioNex Antimicrobial Skin Wipe Towelettes are intended to be used in the field, office, vehicle, clinic, hospital and anywhere else where hand washing with soap and water is not available. The towelettes are formulated with a non-flammable, non-irritating

skin cleansing and moisturizing solution that will leave your skin feeling soft. VioNex Antimicrobial Skin Wipe Towelettes are widely used by firefighters and police officers. Kills germs and cleans skin without running water. Individually pouched for convenience non-flammable and non-irritating

MTX10-1510 5"x7",  
50/pk, 10pk/cs

## **STERIS**



Cal Stat Alcohol Hand Rinse

Cal Stat Alcohol Hand Rinse is a healthcare personnel hand rinse designed for use by personnel involved with routine and high risk patient care when hand washing facilities are not available, or after a routine hand wash when extra protection is desirable. Cal Stat Alcohol Hand Rinse contains 63% isopropyl alcohol in a unique formulation to provide effective topical Antimicrobial activity against MRSA and VRE.

Cal Stat Alcohol Hand Rinse is a waterless hand cleanser, which leaves hands soft and smooth without a sticky after feel

STS1462-K6 4oz, 75/cs



### **CV Medicated Lotion Soap**

CV Medicated Lotion Soap is an Antimicrobial hand wash specifically formulated to be both mild and effective. It contains Triclosan, the same proven Antimicrobial agent used in hospitals, which makes it highly effective against a wide range of bacteria, which may be present on hands. CV Medicated Lotion Soap is superior to harsher formulations, which can cause drying of skin. Its mild formulation is thick, creamy, quick lathering, and easy rinsing. It is the ideal medicated hand wash for controlling bacteria, while leaving hands clean, moist, and healthy

STS6262-24 Quart,  
12/cs



### **Alcure Foamed Alcohol**

Alcure Foamed Alcohol is an alcohol-based surgical scrub and health care personnel handwash with broad spectrum Antimicrobial action against pathogenic microorganisms. The Alcure Foamed Alcohol formulation of 62% ethyl alcohol and emollients provides hand antisepsis and moisturizing at the same time. It spreads easily to facilitate even distribution over all skin surfaces for effective disinfection. Alcure Foamed Alcohol is convenient to use when there are no handwashing facilities available, or as an adjunct to hand washing when extra protection is desirable

STS6395-36 11oz. Can w/ Bracket, 24/cs  
STS6395-57 7oz. Can, 24/cs

### **Alcure Plus Foamed Alcohol**

Alcure Plus Foamed Alcohol is a healthcare personnel handwash used for fast, broad-spectrum reduction of transient

bacteria. It is formulated with 62% ethyl alcohol and a proprietary emollient system to help minimize dryness. Documented efficacy against gram-positive bacteria including MRSA and VRE

STS6399-36 11oz. Can w/ Bracket, 24/cs  
STS6399-57 7oz. Can, 24/cs

### **Septisol NPD With Triclosan**

Septisol NPD with Triclosan is a Healthcare personnel Handwash containing 0.25% Triclosan. Septisol NPD with Triclosan is a mild, natural soap product with broad spectrum Antimicrobial action

STS6281-24 Quart, 12/cs  
STS1357-08 Gallon, 4/cs

## **DEODORIZERS**

### **BARD**



### **Medi-aire Biological Odor Eliminator**

A concentrated unique formula, which chemically eliminates odors, associated with

urine, feces, emesis and necrotic tissue. Available in three scents – Fresh air (A), Lemon (L), and Unscented (U)

BAR7000A Spray Bottle, 1oz. 48/cs  
BAR7000L Spray Bottle, 1oz, 48/cs  
BAR7008A Refiller Bottle, 8oz, 12/cs  
BAR7018A Spray Bottle, 8oz. 12/cs

## **Ultrasound Gels & Accessories**

### **Tyco Kendall- LTP**

#### **Ultrasound Imaging Paper**

Kendall-LTP is the market-leading supplier of diagnostic imaging paper. High quality black and white images with high resolution for ease of review excellent archival properties Designed to fit most popular printers 110mm x 20mm

LTP30751542  
270 prints/roll, 5rolls/box

### **Parker**

#### **Eclipse Probe Cover**

Soft, non-latex, disposable probe cover for medical ultrasound probe/transducers, pre-

gelled inside with Aquasonic 100 Ultrasound Transmission Gel. Fits most endocavity ultrasound probes. 100/bx, 6bx/cs

PAR38-01 Eclipse Probe Cover



### **Transeptic Cleansing Solution**

Multi-purpose spray for professional use in cleaning ultrasound transducer/probe surfaces between procedures. Will not cause transducer membrane swelling. Convenient clear spray dispenser does not contain phenol or glutaraldehyde

PAR09-25 250ml Spray Bottle, 12/cs



### **Thermasonic Gel Warmer**

Warms 3 bottles of ultrasound couplants to body temperature.

Durable wall mountable construction, with warranty Maintains constant temperature Available in 110 and 220

PAR82-04 Thermasonic Gel Warmer



### **Aquasonic 100 Ultrasound Transmission Gel**

World standard for medical ultrasound transmission, AQUASONIC 100 is the most widely used gel for diagnostic and therapeutic medical ultrasound. Recommended for all procedures where a viscous gel is required. Acoustically correct for the broad range of frequencies used in medical ultrasound applications. Gentle wetting of the skin assures excellent transmission of ultrasound waves at all times. Completely aqueous, will not stain clothing, damage transducers or irritate skin

PAR01-02 60g Doppler size tube, 12/bx

PAR01-08 0.25-liter dispenser, 12/bx  
PAR01-34 1 liter SONICPAC with dispenser, 12/cs  
PAR01-50 5 liter SONICPAC with dispenser, 1/bx, 4/cs

### **Sterile Aquasonic 100 Ultrasound Transmission Gel**

World standard for sterile ultrasound gel Sterile Aquasonic 100 is recommended for all medical ultrasound applications where sterility is indicated. Overwrapped, sterilized foil pouches. Acoustically correct for the broad range of frequencies used in medical ultrasound applications Non-irritating, hypoallergenic. Completely aqueous, will not stain clothing, or damage equipment

PAR01-01 20g Overwrapped, sterilized foil pouches, 48/bx



### **Scan Ultrasound Gel**

Our exclusive dual-viscosity formula allows SCAN to be applied as a gel and scan as a fluid for

medical ultrasound procedures where a less viscous gel is desired. Acoustically correct, provides superior transmission of ultrasound waves. Completely aqueous, will not stain clothing, damage equipment or irritate skin

PAR11-28 1 U.S. gallon with dispenser, 4/cs



**Polysonic Ultrasound Lotion**

Multi-purpose ultrasound lotion with superior coupling efficiency in a pleasing cosmetic quality base recommended for all medical ultrasound procedures where a lotion is preferred, also available, Polysonic Ultrasound Lotion with Aloe Vera. Rich, moisturizing formula will not dry or irritate skin. Acoustically correct for the broad range of frequencies used in medical ultrasound. Will not damage equipment or stain clothing

PAR21-28 1 U.S. gallon with dispenser, 4/cs



**Aquasonic Clear Ultrasound Gel**

CLEAR ultrasound gel, Based on the same “can’t be copied” formula of the world standard Aquasonic 100 Ultrasound Transmission Gel without perfume or color For diagnostic and therapeutic medical ultrasound procedures Acoustically correct for the broad range of frequencies used in medical ultrasound. Water-soluble

PAR03-50 5 liter SONICPAC w/ disp, 1/bx, 4/cs

PAR03-54 5 liter ECONOPAC contains: FOUR 5-liter, 2/cs

**Spacelabs Medical**  
Burdick Products.

**Lectro-Sonic Ultrasound Transmission Gel**

BUR027851 Gallon, 4/cs

**Chester Labs**



**LiquaSonic Ultrasound Gel**

LiquaSonic Ultrasound Transmission Gel has a thicker, no-slip formula that remains viscous even when applied to salty skin. Our salt free, formaldehyde-free formula eliminates the air between the transducer and skin, offering exceptional coupling efficiency for more effective treatment/diagnosis. LiquaSonic gel has been formulated to be acoustically correct for the broad range of frequencies used in medical ultrasound applications. It is water soluble and non-greasy preventing damage to expensive transducers and clothing

CHE001222 5-liter collaps, container w/ dispens. Cap, 4/cs

CHE001205 8 fl. Oz. dispensing, 12/cs

**Graham-Field**

**Ultrasound Transmission Gel**



Designed for diagnostic and therapeutic medical ultrasound, the transmission gel is ideal for procedures requiring a viscous or nonspicidal gel. Skin problems derived from prolonged skin contact have been eliminated by the transmission gel's hypoallergenic and completely aqueous formulation. The gel will not stain clothing and is applied and removed with ease. Salt and formaldehyde-free

GF4001GF 8oz., 12/cs  
 GF4002GF 5 liters,  
 1ea

## CHEMO SUPPLIES

**Tyco**  
 Healthcare

### Kendall-LTP

#### Chemo + Plus Gloves

Glove length from tip to cuff is approximately 12" to protect wrists. Textured surface enhances grip. Gloves are hand specific to reduce fatigue. Powder-free with beaded cuff to resist tearing. Tested for resistance to permeation of commonly used

chemotherapy drugs.  
 50/bx, 10bx/cs

LTPCT0192-1 Medium

## Lubricating Jelly

**Johnson Johnson Medical**



#### K-Y Brand Lubricating Jelly

K-Y Brand Lubricating Jelly is a sterile, greaseless, water-soluble jelly used to lubricate gloves, devices or instruments. K-Y Brand Lubricating Jelly spreads easily over rubber, metal or a synthetic surface, forming a transparent lubricating film that adheres well, yet washes away easily. Available in multi-use and single-use packaging

JJ8919  
 4oz. Tube, 12/bx, 6bx/cs  
 JJ8946  
 5g single-use tube, 48/bx,  
 6bx/cs  
 JJ8917  
 2oz. single-use tube,  
 12/bx, 6bx/cs  
 JJ8942  
 2.7g Foil Pack, single-  
 use, 150/bx, and 12bx/cs

JJ8944  
 5g Foil Pack, single-use,  
 150/bx and 4bx/cs

## Triad Medical Products

#### Sterile Lubricating Jelly

Never tacky, our clear, sterile jelly instantly lubricates orifices, preventing friction that can tear delicate tissues. During routine exams, healthcare professionals can easily insert instruments and gloved fingers into the body without discomforting or irritating the patient. Water soluble, non-staining and non-greasy, our lubricating jelly is easily removed with tissues or a towel

TRI10-8344 3gram,  
 144 pk/bx, 12bx/cs

**Tyco**  
 Healthcare

### Kendall



#### VASELINE Pure Ultra White Petroleum Jelly

Bland emollient dressing  
 Soothing Ointment  
 Stable, non-water soluble,  
 topical lubricant

Ointment base for many topical therapeutic agents

KEN8884430200

1oz. Tube, 144/cs

KEN8884430300

3.25oz. Tube, 72/cs

KEN8884433200

5gm Packet, Peel pack,

576/cs

## Thermometers

### 3M Health Care



#### 3M Tempa DOT Single-Use Clinical

##### Thermometers

Accurate, consistent, reliable readings proven in an independent testing laboratory Sterile to reduce the risk of cross contamination Oral readings in just 60 seconds Soft, flexible and Non-toxic, perfect for oral and axillary use, even with newborns and children

#### 3M Tempa DOT Thermometers


MMM5122 Oral and axillary, 100/bx, 20bxs/cs

#### 3M Tempa DOT B Thermometers

MMM5147 Oral and axillary, Non-Sterile 250/bx, 24bxs/cs

## Alaris

Medical Systems




#### P850 Probe Covers

Disposable probe covers for use with models 2080, 2085, 2180, 2185, 4200, 4410, and 4415, 20/bx, and 10bx/dispenser, 25 dispensers/case, for a total of 5,000 probe covers/cs

IVAP850A P850

Probe Covers



#### P90 Probe Covers

Disposable probe covers for use with the CORE CHECK Thermometer System Model 2090. 100/stack, 5 stacks/bx, 10bx/cs, for a total of 5,000 probe covers/case

IVAP90 P90 Probe Covers

## OMRON

### Lens Covers

Lens covers for the MC505 Ear Thermometer for hygienic purposes it is recommended to replace lens covers after each use

OMRMC-505LC

40/bx

## TERUMO

### Probe Covers for Electronic

#### Thermometers

Latex-Free. 1000/box, 5 boxes/case

TER8ETM73 Probe Cover

## Becton Dickinson

### Digital Thermometers and Probe Covers

Delivers fast, easy readings in about 60 seconds. Beeper signals when ready. Auto-shut off. Water-resistant

BD522860 Oral Thermometer, 12/bx

BD522870 Oral Probe cover, 50/bx, 1000/cs

BD522861 Rectal Thermometer, 12/bx

BD522871

Prelubricated Rectal probe Covers, 50/bx, 1000/cs

### Temp-Away Probe Covers

**Oral**  
 BD003700 50/bx,  
 5000/cs  
 BD003715 1000/cs  
**Rectal**  
 BD003700 50/bx,  
 5000/cs

## Welch Allyn



### Welch Allyn Braun ThermoScan Pro 3000 Ear thermometer


Welch Allyn offers the simple, accurate solution for your thermometry testing needs. The Braun ThermoScan Pro 3000 incorporates a new and improved measurement technology. The new ergonomic design provides a small profile, which is comfortable to hold. Built-in-safeguards reduce the risk of cross-types appreciate the quick, non-invasive test procedure

WA03000-200  
 Pro 3000 Ear  
 Thermometers  
 WA03679-200  
 SureTemp 679  
 Thermometer and Braun

ThermoScan  
 Thermometer (Comb.  
 Package)  
 WA05074-800  
 Disposable Probe Covers  
 WA023001-0000  
 Wall Mount for Pro 3000

## EYE Examination & Protection

### Graham-Field



### Snellen Eye Chart

Plastic matte finish, non-reflective with green and red color bars  
 Grommeted hole  
 Washable

GF13-1240 22"x11"

### Eye Test Charts

Non-reflecting, mat finish, eye charts features brass eyelet reinforced holes. Designed for 20 feet test distance. Individually polybagged. Plastic with green and red color bars 22" x 11"

GF13-1241 literate

GF13-1243  
 Kindergarten

### Pocket size Eye Chart, Plastic

Precise, clear, legible, provides 20/800-distance equivalent testing at 14"

GF13-1243-1 6-3/8"x3-1/2", 1ea

### Protective Goggles

Flash Resistant  
 GF9675 24/bx

### Ishihara Color Deficiency Testing

**Books**  
 GF1257 10 Plate  
 Books  
 GF1254 14 Plate  
 Books



### Standard Eye Occluder

Shaped to occlude either eye without pressure. Constructed of black plastic  
 GF1266 6", 1ea

## Stethoscopes

### 3M Health Care



### 3M Littmann Master Cardiology Stethoscopes

The ultimate acoustical stethoscope Patented single-sided “Master Chest-piece with tunable diaphragm lets you monitor both high and low frequencies without having to turn over the chest-piece. Features handcrafted, solid stainless steel chest piece, two-tubes in one design, lightweight comfortable headset, special procedures adaptor for difficult site and pediatric auscultation, nonchill rim, double-leaf binaural spring and patented soft-sealing ear tips. Seven-year warranty

MMM2159 22”, Black  
MMM2160 27”, Black  
MMM2163 27”, Burgundy  
MMM2164 27”, Navy Blue  
MMM2165 27”, Hunter Gm  
MMM2161 27” Black Edition, (Black tubing, black-plated Chestpiece and earpieces)



### 3M Littmann Master Classic II Stethoscope

Outstanding high frequency response Single-sided “Master technology” Chestpiece for excellent acoustics. Dual-frequency, tunable diaphragm Top of the line single lumen for portability. Three-year warranty, plus lifetime service and repair policy

MMM2144L 27”, Black  
MMM2146 27”, Burgundy  
MMM2147 27”, Navy Blue

### 3M Littmann Classic II Infant Stethoscopes

Patented floating diaphragm Solid stainless Steel Chestpiece is optimally sized for little patients. 3/4”inch bell and diaphragm conform to the contours of infant or pediatric patients. Soft-sealing eartips for maximum comfort and acoustical seal Three-year warranty plus lifetime service and repair policy

MMM2114 28”, Black



### 3M Littmann Classic II SE Stethoscopes

The dual-frequency, “tunable” diaphragm allows you to listen to both low and high frequency sounds by alternating pressure on the Chestpiece. Firm pressure for high frequency sounds (Previously heard using the traditional bell). Other features include tough Elastomer polyurethane rim and crack-resistant, reinforced PVC tubing, internal, double-leaf binaural spring helps assure a long life and allows easy adjustment of eartube tension

MMM2201 28”, Black  
MMM2203 28”, Gray  
MMM2205 28”, Navy Blue  
MMM2211 28”, Burgundy  
MMM2214 28”, Teal



### 3M Littmann Classic II Pediatric Stethoscopes

Patented floating diaphragm Solid stainless steel Chestpiece is optimally sized for little patients. 1-inch bell and diaphragm conform to the contours of infant or pediatric patients. Soft-sealing eartips for maximum comfort and

acoustical seal Three-year warranty, plus lifetime service and repair policy

MMM2113 28", Black  
MMM2113R 28", Red



### **3M Littmann Lightweight Stethoscopes**

Patented floating diaphragm assures high acoustic sensitivity and consistent, high frequency responses. Large open bell provides greater sensitivity of low-frequency sounds, Non-chill rim and diaphragm and bell for patient comfort. Lightweight comfort, yet extremely durable with aluminum Chestpiece 1/box, 3bxs/cs

MMM2190 28", Black  
MMM2193 28", Burgundy  
MMM2197 28", Light Blue  
MMM2198 28", Royal Blue



### **3M Littmann Select Stethoscopes**

Single-sided listening for both low and high frequencies (bell and diaphragm combined). Reliable and clear acoustics Lightweight, durable design with scratch resistant finishes on Chestpiece Non-Chill Chestpiece for patient comfort. Backed by a two-year warranty, plus a lifetime service and repair policy

MMM2290 28", Black  
MMM2291 28", Caribbean Blue  
MMM2293 28", Burgundy  
MMM2294 28", Purple



### **3M Littmann Cardiology III Stethoscopes**

The exceptionally versatile adult and pediatric-sided stethoscope also features outstanding acoustics, two patented tunable diaphragms, solid stainless steel Chestpiece, and two tubes in one design, nonchill rims,

adjustable double-leaf binaural spring, and patented Littmann soft-sealing eartips. Five-year warranty

MMM3127 22", Black  
MMM3128 27", Black  
MMM3129 27", Burgundy  
MMM3130 27", Navy  
MMM3134 27", Hunter Green

## **OMRON**



### **Sprague Rappaport- Type LATEX FREE Stethoscope**

5-in-1 multipurpose Quintscope has approximately 22" long tubing. Includes silver Chestpiece, three size open bells, two size replacement diaphragms, and two sets of eartips in a vinyl storage case

OMR416-22-BLK  
Black Tubing

### **PureTone Classic Stethoscope**

Titanium version offers 40% less weight. ENR tube technology reduces external noise and optimizes internal acoustics.

AUTOCLAVABLE.  
Cardiology, Traditional,

Pediatric, and Infant models for a variety of clinical needs 27" and 22" lengths for optimal sound transmission

OMRPT2B27  
27" Black Tubing

### **Echo + Digital Recording Stethoscope**

Digital and conventional acoustics in one, Up to 400% increase in amplification in proprietary Digital Mode. Record, replay, and slow replay of recorded sounds (select models) Data-Link Storage and internet-ready (select models) for telemedicine applications 27" and 22" lengths for optimal sound transmission

OMREC2B27  
27" Black Tubing

## **Welch Allyn**



### **Medallion Gift Stethoscope**

Flawless function in a design of distinction comes in beautiful gift box comes with gold-colored diaphragm and bell. Made of single-lumen tubing with rotating binaurals and

color-coded eartips. Single head. Replacement chestpieces

WA5079-340  
Medallion Gift  
Stethoscope with Black Tubing



### **Tycos Elite Double-Head Stethoscope**

Provides superior acoustical performance in the most important range of 200 to 500Hz. Dual channel PVC tubing provides separate sound pathways to each ear for unrivaled accuracy. Quality stainless steel chestpiece and binaurals designed for optimal weight and durability. Lightweight and comfortable to use and carry. Triple-leaf binaural spring encased in polyurethane permits size adjustment while preventing breakage. Offers flexibility of pediatric capability in the same instrument used on older patients – at no additional cost. Free interactive multimedia CD for learning heart sounds

WA5079-122 22"  
WA5079-125 25"  
WA5079-120  
Tubing for 25" Model

WA5079-115  
Pediatric Diaphragm Assembly

### **Stethoscope Accessories**

WA5079-233 Comfort Sealing Eartips, Large