SAFETY

Blades & Scalpels





Safety Scalpel with Retractable Blade

Designed by a surgeon and a safety engineer the blade can be fully retracted using one hand Safe and easy to open and close ribbed, slip-resistant finger grips ensure a secure hold in gloved hands comfort styling allows for a variety of hand positions during surgical cutting latex safe and sterile

STAD4510 #10 Safety Scalpel STAD4511 #11 Safety Scalpel STAD4515 #15 Safety Scalpel

Futura

The New Thinking Safety



Futura Safety Scalpel

The Futura Safety Scalpel provides superior safety features and retraction technology the user controls in-use or hidden position of blade via locking inset retraction feature hides blade in safety position thereby reducing accidental cuts the device is well balanced, compact and completely disposable

FTAL-342-20 #10, Green, 20/bx FTAL-343-20 #11, blue, 20/bx FTAL-344-20 #15, yellow, 20/bx

Becton Dickinson



BD Bard-Parker SafetyLock Surgical Blade System BD Bard-parker SafetyLock Surgical Blade System offers ease of use, convenience, and

a protective design to help minimize accidental blade injury, especially when disarming simple downward pressure on the SafetyLock handle lever releases the used blade safety and efficiently without the need for a hemostat **BD Bard-Parker** SafetyLock Surgical Blade System, 50/bx, 150/cs BDA371150 #10 Blade BDA371151 #11 Blade BDA371152 #12 Blade BDA371154 #20 Blade BDA371155 #21 Blade BDA371156 #22 Blade BDA371157 #23 Blade **BD Bard-Parker** SafetyLock Handles, 5/cs BDA371036 #3 handle BDA371046 #4 Handle **BD Bard-Parker** SafetyLock Sterile **Disposable Scalpels**, 25/bx, 100/csBDA371670 #10 Handle BDA371672 #12 Handle BDA371675 #15 Handle

IV PRODUCTS

Becton Dickinson



BD Posiflow IV Access **Products**

BD Posiflow eliminates blood drawback into the catheter upon luer disconnect. Compatible for use with lipids and chemotherapy agents Intuitive, easy-to-use and patient-friendly design Low-priming volume, dead space. Straight, twoway fluid path Flow rate capability to meet clinical needs. Flat swabable surface. Latex free. Contains no metallic materials

6" Luer BD385302 Slip Extension Set, 50/bx, 200/cs BD385303 6" Bi Luer-Lok Extension Set. 50/bx, 200/cs BD385304 6" Bi Luer Slip Extension Set, 50/bx, 200/cs 8" "Y" BD385305 Extension Set, 25/bx, 200/csBD385306 34" Dual "Y" Extension Set. 100/cs 100" Dual BD385307 "Y" Drip Set, 25/bx, 100/csBD385308 Multi Dose Vial Adapter, 25/bx, 100/cs BD385309 Port Converter, 50/bx, 200/cs 100" "Y" BD385323 Drip Set, 25/bx, 100/cs BD385324 6" "Y" Extension Set, 25/bx, 200/cs



BD InterLink Blunt Plastic Cannula

Designed for improved performance, the BD **Blunt Plastic Cannula** enhances the Interlink IV Access System. The BD Blunt Plastic Cannula is a truly universal cannula, and can be used with all types of split septum IV injection sites and vials designed for needleless access. It provides simplified access to the interlink injection site with smoother insertion. reduced penetration force, and greater fluid dispersion

BD303367 Vial access Cannula 100/bx, 1000/cs BD303381 InterLink **Blood Collection** Assembly 25/bx, 200/cs BD303400 3cc Syringe with Syringe Cannula 100/bx. 800/cs BD303401 3cc Syringe with Vial Access cannula 100/bx. 800/cs BD303402 5cc Syringe with Syringe Cannula 100/bx, 400/cs BD303403 5cc Syringe with Vial Access Cannula 100/bx, 400/cs

BD303404 10cc Syringe with Syringe Cannula 100/bx, 400/cs BD303405 10cc Syringe with Vial Access Cannula 100/bx, 400/cs



BD Saf-T-Intima Integrated Safety IV Catheters

It's the Complete Safety Catheter. BD Saf-T-Intima integrated safety IV catheters combine advanced needle-stick safety with an innovative closed catheter system design. The result is a virtually bloodless insertion and greater safety for you. The BD Saf-T-Intima integrated safety IV catheter with "Y" adapter allows for pre-priming and flushing during venipuncture. Proven to optimize patient outcomes **BD Saf-T-Intima Integrated safety IV** Catheters, 25/bx, 200/cs DES383312 24ga x 0.75in (0.7 x 19mm) Yellow DES383322 22ga x 0.75in (0.9 x 19mm) Blue DES383335 20ga x 1.00in (1.1 x 25mm) Pink **BD** Saf-T-Intima **Integrated Safety IV** Catheters with Y Adapter, 25/bx, 200/cs

DES383313 24ga x 0.75in (0.7 x 19mm) Yellow DES383323 22ga x 0.75in (0.9 x 19mm) Blue DES383336 20ga x 1.00in (1.1 x 25mm) Pink DES383346 18ga x 1.00in (1.3 x 25mm) Green




Shielded IV Catheters BD Insyte Autoguard Shielded IV Catheters mean added safety for you. They're the safetyengineered version of the popular BD Insyte IV Catheter. The Patented BD Autoguard pushbutton shielding technology is proven to effectively safeguard healthcare workers from accidental needle stick iniuries **BD** Insyte Autoguard Shielded IV catheters, 50/bx. 200/cs DES381412 24ga x 0.75in (0.7 x 19mm) Yellow DES381423 22ga x 1.00in (0.9 x 25mm) Blue DES381433 20ga x 1.00in (1.1 x 25mm) Pink DES381434 20ga x 1.16in (1.1 x 30mm) Pink DES381437 20ga x 1.88in (1.1 x 48mm) Pink DES381444 18ga x 1.16in (1.3 x 30mm) Green

DES381447 18ga x 1.88in (1.3 x 48mm) Green DES381454 16ga x 1.16in (1.7 x 30mm) Gray DES381457 16ga x 1.77in (1.7 x 45mm) Gray DES381467 14ga x 1.75in (2.1 x 45mm) Orange **BD** Insyte Autoguard Winged Shielded IV Catheters, 50/bx, 200/cs DES381512 24ga x 0.75in (0.7 x19mm) Yellow DES381523 22ga x 1.00in (0.9 x15mm) Blue DES381533 20ga x 1.00in (1.1 x 25mm) Pink DES381534 20ga x 1.16in (1.1 x 30mm) Pink DES381537 20ga x 1.88in (1.1 x 48mm) Pink DES381544 18ga x 1.16in (1.3 x 30mm) Green DES381547 18ga x 1.88in (1.3 x 48mm) Green DES381554 16ga x 1.16in (1.7 x 30mm) Gray DES381557 16ga x 1.77in (1.7 x 45mm) Gray

LANCETS

Becton Dickson



Microtainer safety Flow Lancets

No assembly or disassembly prevents injuries. Completely disposable eliminates. Cross contamination. Easy to use Pink Safety Flow Lancet is specifically designed for low blood volume when repetitive punctures are required

BDV366359 Pink, Blade 0.5mm W x 1.4mm D 50/bx, 500/cs BDV366356 Green, Blade 1.0mm W x 1.4mm D 50/bx, 500/cs BDV366357 Blue, Blade 1.0mm W x 1.9mm D 50/bx, 500/cs BDV365759 Yellow, Blade 1.0mm W x 2.2mm D 50/bx, 500/cs

Futura

The New Thinking Safety



Safe-T-Lance Plus Lancets The Futura single use lancets feature automa

lancets feature automatic needle retraction, which inhibits needle stick injures and cross contamination the needle is shielded prior to and after use needle depth control provides for consistent and accurate sampling results

FTA1014 Green 21g, 150/bx, 24bx/cs FTA1020 Blue 25g, 150/bx, 24bx/cs

CHRONIMED

Haemolance Safety Lancets

Single use, disposable devices used to obtain a capillary drop of blood, lancet automatically retracts after use

CHM990300 21 gauge, 150/bx

Bayer

Singlelet Safety Lancets AME6568 Needle 23G, 200/bx



Genie Lancet Ease of use Quick one step blade activation Self-contained, singleuse, and supplies ready to use Full range of blade depths for sufficient blood volume in multiple sample requirements

Advanced protection and safely

Automatic blade retraction minimizes the possibility of injury and reuse No loading or special blade preparation Rapid, precise, consistent incisions providing quality blood samples **Optimum patient comfort** One millisecond lancet penetration and retraction Ultra-sharp surgical stainless steel blades and needles 200/bx, 200/cs

BDV366583 Needle, Orange (Single Drop), 2.25mm, 23 Gauge BDV366582 Lancet, Blue (High Blood Flow), 2.0mm deep, 1.5mm Long BDV366581 Lancet, green (medium to High Blood Flow), 1.5mm deep, 1.5mm long BDV366580 Lancet, pink (Low Blood Flow), 1.0mm deep, 1.5mm long



Quikheel Safety Lancet A gentle, safe bloodsampling device for premature infants, newborns and babies, The BD Quikheel safety Lancet's advanced design technology reduces the risk of injury (cuts or punctures) to the health care worker and provides a virtually pain-free incision BDV368100 Quikheel Preemie Lancet, 0.85mm, 1.75mm, 50/bx, 200/cs BDV368101 Quikheel Infant Lancet, 1.0mm, 2.50mm, 50/bx, 200/cs

> Needles & Syringes

Retractable

Technologies, Inc

VanishPoint Automated Retraction Syringes

- ?? Automated Retraction
- ?? Safe Disposal

?? Non-reusable VanishPoint syringes offer engineered sharps injury protection and feature patented technology which automatically retracts the needle directly from the patient upon activation automated retraction virtually eliminates exposure to contaminated sharps in addition VanishPoint syringes are non-reusable and contain no natural rubber latex Disposal is safe, cost effective 100/bx



| 3CC Vanis | shPoint Syringe |
|------------------|-----------------|
| RTI10301 | 25Gx5/8" |
| RTI10391 | 25Gx1" |
| RTI10311 | 23Gx1" |
| RTI10331 | 22Gx1" |
| RTI10341 | 22G1x-1/2 |
| RTI10351 | 21Gx1" |
| RTI10361 | 21Gx1-1/2 |
| RTI10371 | 20Gx1" |
| RTI10381 | 20Gx1-1/2 |

VanishPoint Automated Retraction Syringes

| iteraction of | |
|--|--|
| 5cc VanishPoin | nt Syringes |
| RTI10531 | 22Gx1" |
| RTI10541 | 22Gx1-1/2 |
| RTI10551 | 21Gx1" |
| RTI10561 | 21Gx1-1/2 |
| RTI10571 | 20Gx1" |
| RTI10581 | 20Gx1-1/2 |
| | |
| 10cc VanishPo | int Syringe |
| 10cc VanishPo RTI11031 | int Syringe 22Gx1" |
| | • • |
| RTI11031 | 22Gx1" |
| RTI11031 RTI11041 | 22Gx1" 22Gx1-1/2 21Gx1" 21Gx1-1/2 |
| RTI11031 RTI11041 RTI11051 | 22Gx1" 22Gx1-1/2 21Gx1" |
| RTI11031 RTI11041 RTI11051 RTI11061 | 22Gx1" 22Gx1-1/2 21Gx1" 21Gx1-1/2 |



VanishPoint Blood Collection Tube Holder

?? Automated Retraction ?? Safe disposal ?? Non-reusable The VanishPoint blood collection tube holder offers engineered sharps injury protection and features patented technology which automatically retracts the needle directly from the patient upon activation automated retraction virtually eliminates exposure to contaminated sharps by protecting the user from both ends of the blood collection needle the tube holder accommodates conventional multisample blood collection needles up to 1-1/2" (38mm) in length and is non-reusable not for use with winged needle sets or luer adapters a small diameter tube adapter is available for use with 10.25mm tubes in addition, the VanishPoint blood collection tube holder and small diameter

tube adapter contain no natural rubber latex

RTI22701BloodCollection Tube HolderRTI22711SmallDiameter Tube Adapter

Futura

The New Thinking Safety



Pro-Ject Blood Needle Holder

Pro-ject needle holders provide a safe and economical mechanism for the disposal of contaminated blood collection needles, Easy one-handed procedure. No need to recap or unscrew from holder

FTAPJ200-10 10/box FTAOJ200-25 25/box FTAPJ200-250 250/bx FTAPJ200-10H For use w/ Hemogard 10/bx FTAPJ200-25H For use w/ Hemogard, 25/box FTAPJ200-250H For use w/ Hemogard, 250/bx

Futura Safety Syringe

Automatic needle retraction requires minimal change to injection technique and virtually eliminates exposure time. The barrel harbors needle postretraction shielding clinician and patient Complies with needle safety regulations and guidelines

FTA103213 3cc, 21g x 1", 50/bx, 10bx/cs FTA103214 3cc, 21g x 1.5", 50/bx, 10bx/cs FTA103223 3cc, 22g x 1", 50/bx, 10bx/cs FTA103224 3cc, 22g x 1.5", 50/bx, 10bx/cs FTA103233 3cc, 23g x 1", 50/bx, 10bx/cs FTA103252 3cc, 25g x 5/8", 50/bx, 10bx/cs FTA103253 3cc, 25g x 1", 50/bx, 10bx/cs

TERUMO

SurGuard Needle Protection System

Unique needle locking mechanism to help reduce the risk of needlestick injuries one-handed activation meets key criteria for safety and ease of use Available with or without attached needle compatible with luer lock and luer slip syringes that conform to ISO 594, ISO/ANSI fittings cost-effective, easy to use positive lock design latex-free 100/bx. 8boxes/case

TERSG1825 SurGuard sheath with 18G x 1"needle TERSG1838 SurGuard sheath with 18G x 1-1/2"needles TERSG1925 SurGuard sheath with 19G x 1"needle TERSG1938 SurGuard sheath with 19G x 1-1/2"needle TERSG2025 SurGuard sheath with 20G x 1"needle TERSG2038 SurGuard sheath with 20G x 1-1/2"needle TERSG2125 SurGuard sheath with 21G x 1"needle TERSG2138 SurGuard sheath with 21G x 1-1/2"needle TERSG2225 SurGuard sheath with 22G x 1"needle TERSG2238 SurGuard sheath with 22G x 1-1/2"needle TERSG2325 SurGuard sheath with 23G x 1"needle TERSG2516 SurGuard sheath with 25G x 5/8"needle TERSG2525 SurGuard sheath with 25G x 1"needle TER129331 SurGuard Short sheath used with 5/8", ³/₄", 7/8" and 1"needles TER129349 SurGuard Long sheath used with $1^{-1}\frac{1}{2}$ needles

Becton Dickinson

BD Saf-T E-Z Set

Winged Infusion Set 50/bx, 200/cs Standard 12in (305mm) Tubing DES387312 27g x 0.38in needle (0.4mm x 9.5mm), Tan DES387323 25g x 0.50in needle (0.5mm x 13mm). Blue DES387326 25g x 0.75in needle (0.5mm x) 19mm), Blue DES387336 23g x 0.75in needle (0.6mm x 19mm), Green DES387346 21g x 0.75in needle (0.8mm x) 19mm) yellow Standard 12in (305mm) **Tubing**, Cont DES387359 19g x 0.88i Needle (1.1mm x 22mm). Pink Shorty 3in (76mm) tubing DES387321 25g x 0.38in needle (0.5mm x 9.5mm) Blue DES387334 23g x 0.75in Needle (0.6mm x 19mm), Green DES387344 21g x 0.75in Needle (0.8mm x 19mm), yellow DES387357 19g x 0.88in Needle (1.1mm x 22mm), Pink PRN 3in (76mm) Tubing DES387358 19g x 0.88in needle (1.1mm x 22mm), Pink



BD Safety-Lok Syringe

BD329464 1ml Insulin Syringe with Permanently Attached 29gx 1/2in. ULTRA-FINE Needle 9U-100) 100/bx, 500/cs BD305553 1ml Tuberculin Syringe with Permanently Attached 27gx1/2in. Needle, 100/bx, 500/cs BD305554 1 mlTuberculin Syringe with permanently Attached 25gx5/8in. Needle, 100/bx, 500/cs BD309606 3ml Syringe only 100/bx, 800/cs BD309592 3ml syringe with 25gx 5/8in. Needle, 100/bx, 800/cs BD309564 3ml Syringe with 23g x 1in. Needle, 100/bx, 800/cs BD309593 3ml Syringe with 22g x 1-1/2in. Needle, 100/bx, 800/cs BD309596 3ml Syringe with 22g x 1in. Needle, 100/bx, 800/cs BD309595 3ml Syringe with 21g x 1-1/2in. Needle, 100/bx, 800cs BD305558 5ml Syringe only 50/bx, 400/cs

BD305559 10ml Syringe only 50/bx, 400/cs BD305561 5ml Syringe with 21g x 1-1/2in. Needle, 50/bx, 400/cs BD305564 10ml Syringe with 21g x 1-1/2in. Needle, 50/bx, 400/cs



Shielding Hypodermic Needles & Syringes Activation-assist technology allows for fast and easy, single-handed needle shielding. It can help prevent needle shielding. It can help prevent needle stick injuries by allowing easy and secure shielding of the needle tip immediately after injection. BD SafetyGlide Syringe a line extension of the BD SafetyGlide Shielding Hypodermic Needles, is similar to the larger BD SafetyGlide sizes, BD SafetyGlide Syringe offers singlehanded immediate activation with minimal change in technique

BD305930 1ml Insulin Syringe with 28g 1/2in Needle, 100/bx, 400/cs BD305932 1/2mlInsulin Syringe with 28g 1/2in Needle, 100/bx, 400/cs BD305935 3/10ml Insulin Syringe with 29g 1/2in Needle, 100/bx, 400/cs BD305945 1 mlTuberculin Syringe with 27g 1/2in Needle, 100/bx, 400/cs BD305946 1ml Tuberculin Syringe with 26g 3/8in Needle, (IB) 100/bx, 400/cs BD305901 25g 5/8in Subcutaneous Needle, 50/bx, 500/cs BD305903 1ml Syringe with 25g 5/8in Sub cut. Needle, 50/bx. 400/csBD305904 3ml Syringe with 25g 5/8in Sub cut Needle, 50/bx, 400/cs BD305916 25g 1in IM Needle, 50/bx, 500/cs BD305902 23g 1in IM Needle, 50/bx, 500/cs BD305900 22g 1-1/2in. IM Needle, 50/bx, 500/cs BD305915 21g 1in IM Needle, 50/bx, 500/cs BD305917 21g 1-1/2in. IM Needle, 50/bx, 500/cs BD305905 3ml Syringe with 23g 1in IM Needle, 50/bx, 400/cs BD305906 3ml Syringe with 22g 1-1/2in IM Needle, 50/bx, 400/cs

BD305907 5ml Syringe with 22g 1-1/2in IM Needle, 50/bx, 400/cs BD305908 10ml Syringe with 22g 1-1/2in IM Needle, 50/bx, 400/cs BD63190 1ml Allergist Tray with 27g 1/2in Needle, (RB) 25/Tray 1000/cs BD66235 1 mlAllergist Tray with 26g 3/8in Needle, (IB) 25/Tray 1000/cs

Shielding Blood Collection Needle

Integrated, single-use blood collection needle for fast and easy needle protection **Advanced Protection** Technology -Single handed activation -Immediate protection -No change in technique -Unique double locking system -Audible click indicator Ease of use -Integrated safety needle shield design -Intuitive in use -Single-use -Bevel up orientation

BDV368607 21G, 1-1/4", 48/bx, 480/cs BDV368608 22G, 1-1/4", 48/bx, 480/cs

Sharps Containers

Becton Dickinson



Multi-Use One-Piece Sharps Containers

Disposal containers of puncture resistant plastic ideal for needles, syringes, and other contaminated sharps features self-sealing locktight cap and safety neck to prevent contents from being extracted contains biohazard caution label on front mounting bracket, lock, and quick release nylon straps available

BD300464 Small. 3.3 gts non-vented, 24/cs BD300467 Medium, 6.9 gts non-vented, 12/csBD300470 Large, 8.2 qts non-vented, 12/cs BD300473 X-large, 5 gal w/open top, nonvented BD305486 Bracket for Phlebotomy & tray size, C-clip, 12/cs BD305497 Strap for 3.3 qts thru 5 gallon, 12/csBD305494 Security lock. 12/csBD305496 Bracket for 3.3 qts thru 9.2 qts, 12/cs

BD305487 Tray size, 1.4 qts, 36/cs BD300458 Tray size, 1.4 qts non-vented, 36/cs BD305488 Small. 3.3qts, 24/cs BD305489 Medium, 6.9 qts, 12/cs Large, 8.2 BD305490 qts, 12/csBD305478 9.2 qts, w/offset neck, 12/cs BD305491 X-Large, 5 gallon w/open top 8/cs BD305477 X-Large, 5-gallon w/large funnel, 8/cs



Point-of-Use Sharps Collectors

Puncture resistant disposal containers for needles, syringes, and other contaminated sharps they feature a counterbalanced door/horizontaldrop design, providing maximum capacity to the fill line with side or horizontal entries 5.4

quart **Horizontal**

BD305445 Red, 20/cs BD305446 Clear, 20/cs BD305452 pearl, 20/cs Horizontal, Next Generation BD305426 Red, 12/cs BD305428 peal, 12/cs Side Entry BD305443 Red BD305444 Ivory BD305425 Ivory Bracketry BD305447 5.4 qt locking bracket, 10cs BD305448 Glove Box Cabinet. 6/cs BD305550 5.4qt Wall mate 1/ea

Multi-Use Nestable Sharps Containers

With an exclusive funnelshaped opening and safety neck, these sharps collectors are Nestable for easy storage they allow point-first drop and reduce the risk of tampering and spilling **3.2 Quart, Vertical Entry** BD305469 Clear, 36/cs BD305471 Red. 36/cs BD305459 Clear, 6/cs Red. 6/cs BD305468 8 Quart, Regular Funnel BD305460 Clear/Red. 24/csBD305479 Red/Red, 24/csBD305463 Clear/Pearl 24/cs**14 Ouart, Regular Funnel** BD305464 Clear/Red. 20/cs14 Quart, Large Funnel BD305480 Clear/Red. 20/cs6 Gallon, Large Funnel BD305465 Clear Red. 12/cs6 gallon, Open Top BD305457 Clear/Red. 12/cs

BD305481 Red/Red, 12/cs

Chemotherapy Sharps Collectors

Chemotherapy Sharps Collectors feature a wide entry port design for disposal of bulky chemotherapy waste in addition; these collectors are designed to offer increased leak resistance

BD305492 9.2qt, Plug Cap, 12/cs BD305493 5 Gallon, Plug Cap, 8/cs

Bemis



Wall Safe Sharps Container

Ideal for patient rooms, Bemis wall safe mounts unobtrusively at any point-of-use location. The rotating cylinder top permits one-handed disposal of sharps, dropping needles horizontally to maximize capacity.

BEM150020 5 qt

BioMed



Sharpx Needle Destruction Unit

Sharpx needle destruction unit saves \$25.000 per 100.000 needles, assisting in compliance with the new, revised OSHA regulations. Incinerate existing needles eliminating blood borne pathogens instead of purchasing expensive safety needles

8001 Sharpx NeedleDestruction Unit0024 Sharpx One Yearservice contract0012 Sharpx Two Yearservice Contract

Sharpx Needle Destruction Accessory Pack

Includes 10 wipes and 10 cups for use with the Sharpx Needle Destruction Unit

5009 Sharpx Accessory Pack

Futura

The New Thinking Safety

Pro-Tec Sharps

Containers For safe disposal of sharps, Pro-Tec containers are constructed of rigid, punctureresistant polypropylene, supplied pre-assembled and ready to use with grip handles and locking lids for safe transport. Containers are available for both biohazard us (red) and chemotherapy (yellow) materials

FTAC1GRED1gallon(150-180 syringecapacity), Red, 18/csFTAC2GREDFTAC2GRED2gallon(260-280 syringecapacity), Red, 12/csFTAC3GREDFTAC3GRED3gallon(460-500 syringecapacity), Red, 8/cs

Tyco Healthcare

Kendall



DEVON Large Volume Sharps Containers Designed for the convenient disposal of large instruments and high volumes the versatile lid on the 7 and 10 gallon containers features both a wide opening for large instruments and an integrated sharps port for smaller items gasketed containers for spill resistance clear lid keeps contents and fill level visible for increased safety convenient molded carrying handles to aid in lifting heavy, filled containers for final disposal

KEN31143533 7.5 Gallon, Red, Round, 6.25"H x 4.75"D x 10.75"W, 6/cs KEN31143665 10 gal, Red, Split, Lid, Sharps Port, 6.25"H x 4.75"D x 10.75"W, 6/cs



Chemosafety Containers with Sliding Lids

Sealing gasketed lids provide leak resistant containment of fluids absorbent pad inside of container helps absorb excess liquid inside of container adjustable sliding lid reduces the possibility of aerosolization and accommodates chemotherapy waste of various types and sizes available carts with wheels and built-in carrying handles for easy transport and lifting of heavy, full containers. Containers provide

temporary and final locking closure for disposal KEN8939 18 Gallon Yellow w/ sliding Lid, 25.25"H x 18.5"D x 13"W, 5/cs KEN8934 12 Gallon Yellow w/ Sliding Lid, 18.5"H x 13"D x 18.5"W, 10/cs



Large Volume Containers with Sliding, Sealing Gasket, and Hinged Lids

Adjustable rotor opening and hinged lids accommodate a variety of sizes of sharps and provide temporary and permanent closure carrying handles and available carts with wheels for easy transport and lifting of heavy, full containers optional gasketed lids provide leak resistant containment of contents. Containers with gaskets have DOT approved UN3291 for transport of regulated medical waste except for cultures and stocks easy simple to use assembly containers lock for final disposal Large Volume containers with sliding Sealing **Gasket lids**

KEN8935 12 Gallon, Red, 18.75"H x 12.75"D x 18.25"W, 10/cs **Large Volume Containers** with Hinged, Sealing **Gasket Lids** KEN8997 8 Gallon, Red Biomax, 17.5"H x 11"D x 15.5"W, 10/cs 8 Gallon. Ken8980 Red, 17.5"H x 11"D x 15.5"W, 10/cs KEN8991 18 Gallon, Red, 26"H x 13"D x 18"W, 5/cs KEN8998 18 Gallon. Red Biomax, 26"H x 13"D x 18"W, 5/cs



DEVON Sharps Containers

Clear, locking tops with horizontal load openings large horizontal openings accommodate larger syringes easy, simple assembly containers lock for final disposal Nestable containers save valuable storage space full range of brackets Plastic and pigments are safe to autoclave or incinerate **Red with Clear Lid** KEN31143699 1 Gallon, 32/cs

KEN31143731 2 Gallon, 30/cs KEN31143681 4 Gallon, 15/cs



Phlebotomy Sharps Containers

Adapts to present blood drawing trays and techniques, including the **AUTO-DROP Blood** Needle Holder System (1510SA) facilitates the disconnection of multisample needles from needle holders a variety of container sizes and openings allow for the disposal of other sharps containers lock for final disposal. For disposal of blood, hypodermic and IV needles, phlebotomy containers are portable and fit into present blood drawing trays. These containers facilitate the disconnection of multisample blood needles from holders

KEN8900SA 1 Qt, Red Container, 100/cs KEN1522SA 2.2 Qt, Red Container, 60/cs KEN1525SA 8 Qt, Red Container, 20/cs

DEVON Phlebotomy Containers

KEN31143657 1 Quart Red Container, 80/cs



The IN-ROOM System Wall Enclosures and Glove Boxes

Sharps containers are secured in a sturdy, locked enclosure. All cabinets are conveniently keyed alike. Made with sturdy ABS plastic for extra protection and long lasting use. Easily monitor fills level through viewing windows on front and sides. Keep gloves readily available with optional Glove Dispenser. Mounted Wall enclosures secure sharps containers, free up counter space and help prevent access to contents in high traffic areas. Hospital personnel can view contents, but unsightly items are not fully apparent to patients and visitors. These sturdy enclosures are made with **ABS** Plastic-the same material used in football helmets

KEN85161H 2 & 5 Qt Wall Enclosure, 1/cs KEN8550LG Glove Dispenser Attaches to 12 and 5 Quart Wall Enclosures, 10/cs KEN8556H 5 Quart Wall Enclosure w/ Glove Dispenser and Sharpstar Container, 1/cs



The IN-ROOM System w/ Mailbox-Style Lids Mailbox-style lid can be kept closed when container is not in use. Economical horizontal drop maximizes use of container volume. Containers lock for final disposal. Wall enclosures. brackets and adhesive backed holders are available. Secure containers wherever needed. IN-ROOM containers for disposal of needles and syringes can be mounted at the bedside or in examination rooms, etc, for point-of-use sharps disposal and to free up counter space

KEN85021 2 Qt, Clear, 20/cs KEN85031 2 Qt, Transparent Red, 20/cs KEN85121 5 Qt, Clear, 20/cs KEN85131 5 Qt, Transparent Red, 20/cs

The IN-ROOM System with Sharpstar Lids Minimizes the potential for overfilling by stopping in the "full" position when it reaches maximum capacity, also allows visual monitoring. Economical horizontal drop maximizes use of container volume. Unique lid design limits accidental or intentional access to containers contents for pediatric, geriatric, or high risk areas. Containers lock for final disposal for use in wall enclosures. Secure containers wherever needed. Allows for visual monitoring. Sharpstar container is the only sharps disposal container awarded the Seal of Accept-ance from the Alliance of Children's Hospitals Inc. Sharpstar System is ideal for pediatric, geriatric or high-risk areas **In-Room Sharps Cont.** w/Sharpstar Lid KEN8508SA 5 Qt, Clear, 20/cs KEN8509SA 5 Qt, Transparent Red, 20/cs **IN-ROOM Sharps Container with Sharpstar** Lid & Counter-Balanced Door KEN8536SA 12 Qt, Clear. 10/cs



The IN-ROOM System w/ Always-Open Lids

Always-Open lids for provide quick and easy disposal. Accommodates larger syringes. Economical horizontal drop maximizes use of container volume containers lock for final disposal. Wall enclosures, brackets and adhesive backed holders are available. Secure containers wherever needed. Accommodates larger syringes. Always open IN-ROOM containers for disposal of needles and syringes can be mounted at the bedside or in examination rooms etc., for quick point-ofuse sharps disposal

KEN851301 5 Qt, transparent Red, 20/cs KEN85221R 12 Qt, transparent Red, 10/cs



DEVON Tortuous Path Sharp Containers One-hand, one-step, notouch sharps depositing Horizontal loading, nonmechanical sharps opening Tortuous-path design impedes hand entry. Clear top standard, translucent base offerings. Plastic and pigments are safe to incinerate. Various mounting brackets and glove box holders are available

KEN31143897 5 Qt, Clear, 30/cs KEN31144010 5 Qt, Translucent Red, 30/cs



Chemosafety Containers with Hinged Lids

Chemosafety containers with hinged lids in appropriate sizes can accommodate larger sizes of chemotherapy waste and sharps including IV bags, tubing, gowns, etc. The sealing gasket and absorbent pad inside help contain excess liquid and guard against spills

KEN8982 2 gallon, yellow w/ hinged lid, 20/cs KEN8985 8 Gallon, yellow w/ Hinged Lid, 10/cs KEN8989 18 Gallon, yellow w/ Hinged Lid, 5/cs



AUTODROP System

Drops needles easily into disposal container in a true one-handed, one-step needle removal and disposal technique. No need for manipulation of blood needles by disposing of the holder with attached needle into proper container AUTODROP System combined with Point-of-Use Sharps Disposal **Containers** Provides healthcare facilities the flexibility to comply with new OSHA regulations

KEN1510SA

AUTODROP System 1 Needle Holder, 60 threads, 50pkg/cs



Multi-Purpose Sharps Containers

Adjustable rotor opening or hinged lids accommodate a variety of sizes of sharps and provide temporary and permanent closure. Easy, simple to use assembly Containers lock for final disposal Nestable containers save valuable storage space available holders and brackets secure containers wherever needed Multi-Purpose, w/Rotor **Opening Lid** KEN8920SA ¹/₂ Gal Needle Syringe collection, 60/cs KEN8950SA 5 Ot, Red, round. 40/cs KEN8970 2 Gal Red. 20/cs**Multi-Purpose Container** with Hinged Lid KEN8526 3 Gal. Red, 10/cs

Transportable Sharps Containers

Simple, one-handed operation opens and closes the 1 Quart disposal openings convenient size fits into carry bags for nursing, EMS and Fire/ Rescue. Leak-resistant screw cap offers secure closure between uses when properly closed. Large disposal opening easily accommodates small syringes and other sharps waste 2 Gallon Rotor Lid adjusts to accommodate small and large sharps Built-in handles for transport to usage site and final disposal

KEN8961 8 Quart, Transparent Red, 20/cs

DEVON Brackets and Accessories KEN31139747 Extra Keys for Hidden Bracket Set 4841-K, 4841-HK. 10/cs

Horizontal-Drop Opening Lids

For use on countertops or mounted on wall Economical horizontal drop maximizes use of container volume large horizontal openings accommodate larger syringes easy, simple to use assembly containers lock for final disposal Nestable containers save valuable storage space adhesive-backed holders and wall mounting brackets are available to secure containers wherever needed plastic and pigments are safe to autoclave or incinerate

KEN89651 2 Gal Container, Red with Clear Lid. 20/csKen89671 2 Gal Container, Transparent Red with Clear Lid, 20/cs





Maxxim Sharps Collectors Red sharps collector with clear locking lid

MXX8702 1 quart, 72/cs

MXX8705 16 quart, 10/csMXX8703 24/csMXX8704 24/cs

4 quart, 8 quart,

Post Medical



Automatic Needle Disconnect Sharps Collector

The Automatic Needle Disconnect (A.N.D.) phlebotomy container from Post medical is the ideal system for safe disposal of blood collection needles, lancets and butterflies. Molded of durable, punctureresistant plastic, this unit's compact size fits easily into phlebotomy trays. Separate A.N.D. units are designed to work with BD Vacutainer, Sherwood Monoject, Terumo, Nissho, and Greiner blood collection needles. Simply push holder downward with ONE HAND – the blood collection needle is automatically unwound and deposited safely into the container. 24 per case

PMIBE-100 200 **Needles** Capacity

BioMed



NICSx Compliance **Assistance Program** NICSx web-based software publishes your Sharps Exposure Control Plan. Trains and certifies caregivers on plan, records safety product evaluations, and logs any needle stick incidents. Review data preventing future needle stick incidents

N005 NICSx for 1 to 5 Employees N010 NICSx for 6 to 10 Employees N10P NICSx for Additional 10 Employees N100 NICSx for up to 100 Employees N025 NICSx for up to **25** Employees N050 NICSx for up to 50 Employees NMAN NICSx paper Manual

M & I Medical Sales, Inc., 4711 SW 72 Ave, Miami, FL 33155, U.S.A. Ph 305-663-6444, Fax 305-663-6201, contact@mimedical.com, www.mimedical.com