Urology

BAGS

Allegiance

Urinary Drainage Bags with Antireflux Chamber with SPLASHGUARD II Drain Spout

(Kendall) These drainage bags offer a SPLASHGUARD II spout that makes drainage easy and helps protect you from contact with potentially hazardous body fluids the Antireflux chamber minimizes urine backflow each bag includes nylon hook handing device and string hanger that fit securely on bedrails and wheelchairs, urine sampling port and 11/32"tubing sterile bags packaged individually 20 bags/cs

BGHKB5900 poly bag PK only



Drainage Bag w/ Mono-Flo Antireflux II Drain Spout

(Kendall) This 2L drainage bag features a Mono-Flo Antireflux device that automatically closes when you tilt, lift or lay the bag on its side, helping prevent urine backflow SPLASHGUARD II drain spout makes drainage easy and helps protect you front contact with potentially hazardous body fluids universal hanging device fits easily over the bedrail and doubles as a carrying handle for ambulatory patients each sterile, individually packaged bag is supplied with 11/32" tubing, nylon hook handing device and string hanger 20 bags/cs

BGHKB5912 Drainage bag



Leg Bags

Our leg bags feature an Antireflux valve, universal connections, a twist open drainage valve and adjustable leg straps with buttons bags are individually packaged with a sterile fluid path

BGHP4015M Medium, 500cc, 40/cs BGHP4017L Large, 1000cc, 40/cs

BARD

Dispoz-a-Bag Leg Bags

Made from heavy duty plastic film for long wear time, odor containment and leak resistance each leg bag contains an antireflux valve reusable, sterile

BAR150102 19oz, 48/cs BAR150103 32oz, 48/cs



Center-Entry Urine Drainage Bag

- -Bard Safe Sampling Port allows for needle-free aspiration
- -Anti-Reflux device
- -Double hook hanger
- -Single-use, sterile

BAR153504 2000ml



Urine Drainage bag with Anti-Reflux Chamber

- -Bard Safe sampling Port allows for needle-free aspiration
- -Single hook design with string hanger single-use, sterile

BAR154002 2000ml, 20/cs

Tyco Healthcare

Kendall



CURITY Urine Leg Bags

Anti-Reflux valve longer catheter adapter for greater security with all types of catheters twist valve for easier emptying and less spill, no bottom cap to lose or misplace soft, cloth straps latexfree leg straps

KEN145516 Large 32oz/ 1000ml, 50/cs



KENGUARD Anti-Reflux Chamber Urinary Drain Bags

Contoured design for more even filling and complete emptying 11/32" drainage tube provides quick drainage, eliminates standing column of urine universal hanger fits securely on bedrails and wheelchairs convenient carrying handle for ambulatory patients

KEN3512V 2000ml bag with Anti-reflux chamber and VELCRO hanger KEN3502 2000ml bag without anti-reflux chamber, 20/cs



CURITY mono-Flo Anti-reflux Device Drain Bag with SAFEGUARD II Drainage Device

Mono-Flo Anti-Reflux valve automatically closes when bag is tilted, lifted, or laid on side large diameter drip chamber provides dry barrier against bacteria, facilitates passage of clots Splashguard II Drainage Device prevents backsplash of urine SAFEGUARD Needles less sampling port prevents needle sticks

KEN6300 2000ml MONO-FLO, Polybag package, 20/cs KEN6206 2000ml Anti-reflux chamber with SPLASHGUARD II, CSR Wrap, 20/cs

Catheters

Allegiance



Siliconized Latex Foley Catheters

Our catheters feature high-quality natural latex that provides patient comfort while remaining sufficiently firm for insertion 12 catheters per case

2-way catheters with 5cc balloon

BGHPE512 12 14 BGHPE514 BGHPE516 16 BGHPE518 18 BGHPE520 20 BGHPE522 22 2-way catheters with 30cc balloon BGHPE616 16



18

Red Rubber urethral Catheters

BGHPE618

Our red rubber urethral catheters are made of X-ray-opaque latex, which is soft and flexible to minimize trauma to delicate mucosal tissue single-use catheters are sterile 50 catheters per case

BGHPRU110 10 BGHPRU112 12 BGHPRU114 14 BGHPRU116 16 BGHPRU118 18

Silicone Elastomer-Coated Two-way Foley Catheters

Our catheters feature permanent silicone-Elastomer coating bonded to latex for a smooth surface, helping minimize encrustation the inert formulation protects against irritation, which allows the catheters to remain in place for a longer time period with comfort two-way catheters are sterile and individually packed caution: Federal (USA) law restricts this device to sale by or on the order of a physician 12 catheters per case

Catheters with 5cc balloon

BGHPS516 16 BGHPS518 18 BGHPS520 20 Catheters with 30cc balloon PS620 20

BARD

5cc Coude Bardex Lubricath Latex Foley Catheter (Carson Model)

-Coude tip design facilitates negotiation of the prostatic urethra -Medium length – olive Coude tip with a single drainage eye -30cc balloon for postoperative hemostasis (0100LXX)

BAR0168L14 14Fr, 5cc, 12/cs BAR0168L16 16fr, 5cc, 12/cs BAR0168L18 18Fr, 5cc, 12/cs BAR0168L20 20Fr, 5cc, 12/cs

5cc Coude Bardex Lubricath latex Foley Catheter (Tiemann Model)

-Coude tip design facilitates negotiation of the prostatic urethra -Medium length – olive Coude tip with a single drainage eye -Red latex construction for added catheter stiffness and Radiopacity

BAR0102L16 16Fr, 5cc, 12/cs BAR0168L18 18Fr, 5cc, 12/cs


Util-Cath Vinyl Catheters

Sterile, whistle tip BAR015914 14Fr, 6in, 200/cs



Bardex Lubricath Latex Foley Catheters

-Lubricath hydrophilic polymer coating which helps reduce irritation of the urethral mucosa and which helps minimize encrustation -High-quality Bardex construction -Single-use, sterile Short-Round tip w/ two opposing eyes BAR0165L12 12fr, 5cc, 12/cs BAR0165L14 14Fr, 5cc, 12/cs BAR0165L16 16Fr, 5cc, 12/cs BAR0165L18 18Fr, 5cc, 12/cs BAR0165L20 20Fr, 5cc, 12/cs BAR0165L22 22Fr. 5cc. 12/cs BAR0165L24 24fr, 5cc, 12/cs



Medium length – Round tip with two staggered eyes BAR0166L16 16fr, 30cc, 12/cs BAR0166L18

Bardex All-Silicone 2-Way Foley Catheters

18Fr, 30cc, 12/cs

-100% Silicone
constructions
-Single-use, sterile 5cc
Bardex all-silicone Foley
Catheter
-Short length – round tip
with two opposing eyes

BAR165816 16Fr, 5cc, 12/cs

Bardex Silicone-Elastomer Coated latex Foley Catheter

-Silicone Elastomer coating provides protection against urethral irritation -5cc balloon Foley catheter features a short length – round tip and two opposing eyes (0165VxxS) -Single-use, sterile

BAR0165V16S 16Fr, 5cc, 12/cs

Urethral Catheters

-Semi-rigid construction to facilitate negotiation of the prostatic curve -Coude olive tip with a single drainage eye located on the top of the prostatic curve -Raised indicator to identify direction of Coude tip -Semi-rigid Red Rubber

Bard Coude Tip Latex Urethral Catheters

One Eye (Tiemann Model)

-Single-use, sterile BAR010114 14Fr, 12/cs BAR010116 16Fr, 12/cs BAR010118 18Fr, 12/cs

Two Eyes (Tiemann Model)

-Coude olive tip with two staggered drainage eyes to provide for higher drainage flow -Single-use, sterile BAR120614 14Fr, 12/cs BAR120616 16Fr, 12/cs

Bard Red Rubber All-Purpose Urethral Catheter

- -Radiopaque red rubber catheter with a round, hollow tip can be used as a Robinson or a Nelaton catheter
- -Two opposing drainage eves
- -16"long
- -Single-use, sterile

BAR177714 14Fr, 100/cs BAR177716 16Fr, 100/cs

Tyco
Healthcare

Kendall



KENGUARD Silicone- Coated Foley Catheter

Silicone coating provides a smooth exterior shaft, facilitates insertion

5cc Balloon, 2-way

KEN3560 16Fr 17", 10/cs KEN3563 18Fr 17", 10/cs

30cc Balloon, 2-way KEN3607 18Fr 17", 10/cs



ULTRAMER Foley Catheter

100% silicone Foley catheter Add-A-Foley Tray

KEN8256 10cc Prefilled syringe, 10/cs

Mentor



Self-Cath Female

Self-Cath 240 is designed for women who practice intermittent selfcatheterization self-cath is made of medical grade PVC, then siliconized for friction-free insertion Mentor's unique fire polished eyelets enhance patient comfort

MNT240 14Fr, 6" Long, 30/bx



Self-Cath Male

The self-cath family of latex-free catheters offers patients a variety of options for intermittent self-catheterization each catheter is made with the highest quality standards dedicated to patient comfort and ease of use the Self-Cath Straight Tip is 100% latex-free with fire-polished eyelets and a siliconized surface for smooth insertion

MNT450 14Fr, 16" Long, Funnel End, Curved Packaging, 50/bx



Freedom Cath Male External Catheter

Mentor offers a complete line of latex and non latex male external catheters a variety of adhesive styles, sheath lengths and one two-piece configurations are designed to suit a variety of patient needs a complete incontinence system can be achieved by combining a Mentor's male external catheter and Leg Bag Freedom Cath is a self-adhering, one-piece catheter for secure everyday use, combining a comfortable sheath with a wide, watertight adhesive seal

MNT8400 Large, 100/bx

RUSCH



All purpose Robinson/ Nelaton Catheters

Red rubber latex, 2 opposed eyes, X-ray opaque length 40cm (16") Sterile, single use 100/bx

RUS3510-08 8Fr RUS3510-10 10Fr RUS3510-12 12Fr RUS3510-14 14Fr RUS3510-16 16Fr RUS3510-18 18Fr RUS3510-20 20Fr RUS3510-22 22Fr



Catheters

RUS564000 Cath Word Bartholin Gland Catheter, 5cc Balloon, 5cm Length, Latex 6/bx RUS220005140 Female ERU Siliconized PVC, 14Fr, 6.5" Single Eye, 50/cs

RUSCH



Foley Catheters RUS318116 PTFE Teflon Coated, 2 ways, 5cc, 16Fr, 16" Coude

Single Eye Tiemann Model, 10/cs RUS318118 PTFE Teflon Coated, 2 Way, 5cc, 18Fr, 16" Coude Single Eye Tiemann Model, 10/cs



RUS180705160 Silicone Coated, 2 Way, 5cc, 16Fr, 16" 2 Opposed Eyes, 10/cs RUS180705180 Silicone Coated, 2 Way, 5cc, 18Fr, 16" 2 Opposed Eyes, 10/cs

Catheter Kits

BARD

Foley Catheter kits

Bard Foley Kits feature high quality Bardex Foley catheters in a selection of coatings: Lubricath hydrophilic polymer or Silicone-Elastomer these coatings provide protection against irritation of the urethral mucosa, help to minimize encrustation and enhance patient comfort -Each kit contains: -5cc Latex Foley catheter – Prefilled water inflation syringe –5 gm Lubricant packet -The Bardex Lubricath Foley Catheter provides for smoother insertion

and enhanced patient comfort through the hydrogel coating's unique ability to absorb 800% of its own volume in fluid, thus creating a hydrophilic "cushion" between the catheter surface and the urethra -Single-use, sterile

BAR710012S-26S 12-26 Fr, Silicone-Elastomer Coated Foley Catheter, 20/cs



Female, Infant and Pediatric Catheter kits

-Closed system reduces the risk of specimen contamination and associated false readings -Collection tube is centrifuge-ready -Kit includes: pair of vinyl gloves, towel, 3 iodophor swabsticks, 3g lubricant packet and label -Single-use, sterile

BAR0035380 8Fr, Basic Female Cath kit/rigid (gloves, swabs, and label), 100/cs DAV0035630 5Fr, Infant, Soft, 50/cs DAV0035640 8Fr, Pediatric, Soft, 50/cs



FAST-CATH Female Urethral Catheterization Kit

Kit Contents: Underpad, Vinyl Gloves, PVI swab sticks (3), 15ml Specimen Tube and Flipcap, Sterile Lubricating Jelly, 10gm, Patient Label KEN3411 female urethral Cath kit, 25/cs

TRAYS

Allegiance



Basic Foley Catheterization trays With Syringe and PVP

Sterile Foley catheterization trays are available with either a 10 or 30cc printed syringe with sterile water and Luer tip no catheterincluded trays include:

- -Outer basin, 1000cc
- -Water-resistant underpad
- -One pair latex gloves
- -Fenestrated drape
- -Three povidone-iodine swab sticks (1% available iodine)
- -Lubricating jelly packet

BGH3T5786A W/ 10cc Syringe BGH3T5785A W/30cc Syringe

Basic Urethral Catheterization Trays

With BZK and Specimen Container

Our sterile catheterization trays are available with a red rubber latex or plastic 14Fr catheter or without catheter trays include:

- -Outer basin, 1200cc
- -Water-resistant underpad
- -One pair latex gloves
- -Fenestrated drape
- -Three Benzalkonium chloride (1:750) swab sticks
- -Lubrificant jelly packet
- -Specimen container 20 trays per case

BGH3T6033 With red rubber latex catheter BGH3T6030 With plastic catheter BGH3T6031 Without catheter



Universal Foley Catheterization Trays With Hospital Wrap

(Kendall) These convenient sterile trays come with or without catheter 20 tray per case trays include:

- -CSR wrap
- -Anatomically shaped basin, 1200cc
- -Water-resistant absorbent underpad
- -One pair latex exam gloves
- -Fenestrated drape
- -Prepping balls

- -Plastic forceps
- -10% povidone-iodine solution, 1oz
- -Lubricant jelly, 2.7g
- -Sterile specimen container

Trays without catheter BGHKB4708B

Prefilled sterile water syringe, 10cc

Trays with Elastomercoated latex Foley catheter BGHKB4709S

14Fr catheter Prefilled sterile water syringe, 5cc BGHKB4710S

16Fr catheter Prefilled sterile water syringe, 5cc



Urethral Catheterization trays

20 trays per case trays include:

- -Anatomically shaped basin, 1200cc
- -Rayon balls
- -Removable prep cup
- -Water-resistant underpad
- -Fenestrated drape
- -Specimen container
- -One pair of latex gloves
- -Lubricating jelly, 2.7g
- -10% povidone-iodine prep solution, 1oz
- -Forceps
- -Hospital wrap

BGHKB4731B

Tray w/out cath

BGHKB4730B

Tray with 14Fr vinyl

catheter

BGHKB4733B

Tray w/ 14Fr red rubber latex catheter



Foley catheterization trays with Antireflux Device drainage Bag and Drain Spout

(Kendall) These sterile trays offer a bag with Splashguard II drain spout that makes drainage easy and helps protect you from contact with potentially hazardous body fluids trays include:

- -Patient information pamphlet
- -CSR wrap
- -Water-resistant absorbent underpad
- -One pair latex exam gloves
- -Fenestrated drape
- -10% povidone-iodine solution, 1oz
- -Prepping balls
- -Plastic forceps
- -Sterile lubr. jelly, in bellows syringe, 10g
- -Prefilled sterile water syringe, 10cc
- -Tape strips
- ape surps
- -Sterile specimen bottle
- -Drainage bag with

Mono-Flo Antireflux valve, 2000ml

- -Tubing, 11/32" diameter
- -Splashguard II drain spout
- -Bedsheet clamp

Trays with 100% silicone cath 10 trays/cs

BGHKB5936A 16Fr Cath W/



Foley Catheterization Trays With Antireflux Device Drainage Bag Drain Spout & Lube Jelly

(Kendall) these sterile catheterization trays offer a bag with Splashguard II drain spout that makes drainage easy and helps protect you from contact with potentially hazardous body fluids a "Wash your Hands" reminder encourages procedural asepsis trays include:

- -Patient information pamphlet
- -Water-resistant absorbent underpad
- -One pair latex exam gloves
- -Fenestrated drape
- -10% povidone-iodine solution, 1oz
- -Prepping balls
- -Plastic forceps
- -Sterile lubricating jelly, 5g
- -Prefilled sterile water syringe, 10cc
- -Tape strips
- -Sterile specimen bottle

-Drainage bag with Mono-Flo Antireflux device, 2000ml -Tubing, 11/32" diameter -Splashguard II drain spout -Bedsheet clamp 10 trays per case Tray without catheter and **CSR** wrap BGHKB5970 Trays with silicone **Elastomer-coated latex** catheter and CSR wrap BGHKB5974S With 14Fr catheter BGHKB5976S With 16Fr catheter **BGHKB5978S** With 18Fr catheter

Tyco
Healthcare

Kendall



KENGUARD Foley Catheterization Tray

Includes #3512 Drain
Bag-2000ml capacity,
absorbent underpad, latex
gloves, fenestrated drape,
PVI swab sticks (3),
lubricating jelly, 5 gram
pouch, CSR wrapped,
specimen container with
cap, 10cc pre-filled
syringe

ADD-A-Foley Tray KEN3531 30cc Prefilled Syringe with 3502 Drain Bag



CURITY MONO-FLO Antireflux Device Foley Trav

Ultramer Foley Catheter
KEN6014 16Fr, 5cc,
10/cs
100% Silicone Foley
Catheter
KEN6140 16Fr, 5cc,



KENGUARD Universal Catheter Tray (Peel Back Lid)

All trays contain (unless otherwise indicated):
Lidded tray, underpad otherwise indicated):
Lidded tray, underpad, latex gloves, fenestrated drape, PVI solution, prepping balls, forceps, lubricating jelly, graduated basin

KEN3505 30cc Prefilled Syringe, 20/cs

KENGUARD Universal Catheter Insertion Tray (Peel Back Lid)

All trays contain (unless otherwise indicated): Lidded tray, underpad, latex gloves, fenestrated drape, lubricating jelly, and graduated basin KEN76020 30cc Prefilled Syringe, BZK Swab Sticks (3), 20/cs



CURITY Urethral Catheter Tray (Open)

All trays feature (unless otherwise indicated): CSR wrap, underpad, latex gloves, PVI solution 10%, prepping balls, forceps, lubricating jelly, specimen container

KEN3217 14Fr, ULTRAMER Red Rubber, 20/cs

RUSCH



Catheter Tray

RUS75735 urethral Cath Tray w/ Red Rubber Cath, 14Fr (1100ML Basin PVI Swabs), 20/cs

BARD



Bardex Lubricath Foley Tray

-5cc Bardex Lubricathlatex foley catheter-Bard safe sampling port

- -Patented Tamper-Evident seal
- -Uro-Prep tray
- -Single-use, sterile

BAR899614 14Fr, 10/cs BAR899616 16Fr, 10/cs BAR899618 18Fr, 10/cs

Bard Lubricath Foley Trav

- -5cc Standard Bard Lubricath latex foley catheter
- -Bard safe sampling port
- -Patented Tamper-Evident seal
- -Uro-prep tray
- -Single-use, sterile

BAR898316 16Fr, 10/cs



Bardex Lubricath Foley Tray

- -5cc Bardex Lubricath latex foley catheter
- -Bard safe sampling port
- -Patented Tamper-Evident seal
- -Uro-prep tray
- -Single-use, sterile

BAR899716 16Fr, 10/cs



Infection Control Add-A-Foley Tray

-Patented Microbicidal outlet tube

-Complements Bardex I.C. system

BAR899100 10/cs



Bard Bilevel Universal Travs

-Uro-prep Tray with 5g lubricant pack -Single-use, sterile BAR792116 With 16Fr, 5cc Bardex Lubricath Latex Foley Catheter, 20/cs BAR781600S With 16Fr, 5cc Bard Silicone-Elastomer Coated Foley Catheter, 20/cs BAR782102 With 30cc pre-filled water syringe, 20/cs

Bard Bilevel Urethral Catheterization Procedure Trays

Basic tray contents:
-Sequence packed in a
bilevel configuration,
which separates prepping
components from
catheterization
components
-Prepping components
include; drape, underpad,
latex gloves, forceps,
rayon balls, 5g lubricant

rayon balls, 5g lubricant packet, specimen container, label, povidone-iodine solution -Single-use, sterile

BAR772100 Basic Tray without either a catheter or a collection bag, 20/cs

Urine Collectors

BARD

Bard Midstreams with Funnel Collector

-Pre-connected large target collection funnel -Instructions for male and female patients -Lid protector minimizes contamination potential -Single-use, sterile fluid path

BAR842801 3oz. container with BZK Wipes, 50/cs

TycoHealthcare

Kendall



Pediatric Urine Collectors

Non-Sterile
Ken145500 100ml,
50/bx, 10bx/cs
Sterile
KEN14501 100ml,
30/bx, 10bx/cs



SAGE Midstream Catch Kits

KEN2110SA Calculi Strainer 100



Commode Specimen Collection System

Ideal collection unit for urine and stool samples or kidney stones Wipe lip and firm base helps prevent splashing, spilling, tipping clearly marked gradients are easy to read KEN2500SA Commode Specimen Measuring System, 100/cs KEN4014 CURITY SPECIPAN with clearly marked gradients, 100/cs



PRECISION 24 hour Specimen Container 3000ml capacity lined

plastic screw cap minimizes leakage easyto-ready graduations

KEN5000SA 3000ml, 40/cs



PRECISION Specimen Containers

Complete line of specimen containers allow for sterile specimen collection protects against leakage clearly marked gradients

17000 Series: Polyethylene Screw Caps, Sterile, individually wrapped with screw-on caps KEN17000 4.5oz Graduated Containers, 200/cs

KEN17080 4.5oz Graduated Containers, Attached ID Label, 100/cs

17000 Series: Polyethylene Screw Caps, Non-Sterile, individually wrapped with screw-on caps KEN17500 4.5oz Graduated Containers, 400/cs

13000 Series: Snap Caps, Sterile, individually wrapped containers KEN13000 5oz Graduated Containers, 200/cs

Precision Dynamics



Lil-Katch Pediatric Urine Collector

Versatile design provides a comfortable fit for infant or pediatric use, male or female medically approved adhesive protects sensitive skin adhesive completely surrounds opening, assuring a leak-proof fit and a clean, usable urine specimen enlarged foam area provides ease of application and comfortable fit, while keeping walls of collector open available sterile or clean, 50/bx capacity 7fl. Oz. 200ml

PRE05001-00MCF Non-Sterile PRE05002-00MCF Sterile

Premium Plastics



Specimen Collection Commode and Toilet Collection Pan

Securely fitting standard toilets and commodes, this specimen collector's winged edges facilitate application, handling and removal each Gent-L-Pan has easy-to-read graduations in 50cc and 1oz increments up to 800cc's a handy pour permits specimens to be discarded quickly and

safety Bulk (2 bags of 50ea) 100/cs

PP4707 800cc, 100/cs



Calculi Strainer

Fine mesh filters stones or other foreign matter quickly and easily strainers fit securely on to 24 hour specimen containers, remaining stable throughout the entire filtering process built-in ridge serves as a handle for simple application and removal

PP4712 White, 25/bg, 4bg/cs



Specimen Containers

6.5oz containers are available in polystyrene and polypropylene. Polystyrene enhances clarity and polypropylene improves durability both come with a convenient pour spout for simplified specimen retrieval each container is graduated in 10cc (ml) and 5oz increments up to 6.5oz printed lids can be

purchased separately or with specimen containers

Polypropylene
PP4629 Specimen
Container w/o lid, 20/bg,
25bg/cs
PP4646 Specimen
Container w/ Printed
Paper Lid, 20/bg, 25bg/cs
Polystyrene
PP4630 Specimen

Container w/o Lid, 20/bg, 25bg/cs
PP4647 Specimen
Container w/ Printed

Container w/ Printed Paper Lid, 20/bg, 25bg/cs

Specimen Containers w/ Liqua Guard Lid

Available with several packaging options, each specimen container is marked in 10cc (ml) and .25oz increments up to 120cc's containers are available with Liqua Guard lids, ensuring leak resistance and making them perfect for transport each Liqua Guard lid is engraved and has space for writing date, name, doctor and room number depending on package on package lids can be purchased separately or with containers

PP4651 4oz Specimen Container w/o Lid, 20/bg, 25bg/cs PP4653 4oz. Specimen Container w/ Liqua Guard (Turquoise), 20/bg, 25bg/cs



Specimen Containers W/ Liqua Guard Lid and Label

4 oz sterile specimen container w/ gray Liqua Guard lid & label each specimen container is marked in 10cc (ml) and .25oz increments up to 120cc's Liqua Guard lids ensure leak resistance and make them perfect for transport each Liqua Guard lid is engraved and provides ample space for writing the date, name, doctor and room number polybagged

PP4928 4oz. Specimen Container w/ Gray Liqua Guard Lid & Label, 100/cs

Vollrath



Sterile specimen Containers

Vollrath Group's 4oz. specimen container features full-turn screw thread, provides excellent leak-resistant seal tough, durable inert material. Integrity seal includes

patient ID label on body Leak-resistent seal meets NCCLS specimen handling guidelines for primary containers. (M29-P, Vol7, No.9). Integrity seal indicates internal sterility when intact. Cup is graduated in both ounces and ml on two sides

VOL01055 100/cs



Specimen Containers without Cap

Vollrath Group's nonsterile specimen container without cap is graduated for easy measuring in cc and oz

VOL02053 4.5oz. Polypropylene 25/sleeve VOL02263 6oz. Styrene, 25/sleeve, 500/cs VOL02265 6oz. Polypropylene, 25/sleeve



Specimen Containers with Cap

Vollrath Group's screw Cap Specimen Containers are graduated for easy measuring in cc and oz. Area on white cap allows for patient identification. Caps are sleeved and cups are stacked VOL02051 4.5oz Polypropylene, 25/Sleeve VOL02052 4.5oz. Styrene, 25/Sleeve, 500/cs



Measurite Specimen Pan With Black Graduations

Vollrath Group's improved injection molded MeasuRite Specimen Pan offers accurate urine output measurements, even at volume as low as 50cc's no need for additional measuring black printed graduations in oz and cc. Ounce graduations are 2,4,6,8,10,12,16,20,24,28, 32, cc's graduated at 50, then 100 through 1000 in 100cc increments increased wall thickness to provide improved stability during transport large capacity to reduce the risk of accidental spills and permit urine output measurements up to 1000cc's a molded area on the top front with space to record the doctor, patient, and room number a grooved pour spout for easy and accurate pouring of contents from the front or back. The specimen pan also has a hole for

convenient hanging a wide flange with reinforcing ridges for strength and stability when emptying pan or carrying pan with contents a sturdy, stackable one-piece design that can be nested for compact storage graduations in cc's and oz's have large, angled, easy-to-read numbers

VOL02071 100/cs

Miscellaneous

Allegiance

Catheter plug/Protector

(Kendall) This catheter plug/protector is sterile and individually packaged 200 plug/protectors per case

BGHKB9660 catheter Plug/Protector

BARD

Urinary Cunningham Clamp

Used to control stress incontinence or dribbling of urine constructed of soft sponge rubber and can be adjusted to fit comfortably reusable, non-sterile

BAR004053 Regular, 2in

Bard Penile Clamps

For catheterization procedures usage stainless steel, reusable, non-sterile

BAR600404 Zipser Model

DALE



Dale Foley Catheter Holder

Latex-Free. Put an end to the problems associated with tape, and reduce the potential for infection and trauma. Easy to put on and take off, the Dale foley holder features a patented locking device that securely anchors the catheter in place it helps minimize patient tampering, facilitates catheter traction. encourages mobility, and can be moved from leg to leg with total ease the legband is comfortable, due to a unique stretch material that distributes stretch material that distributes compression evenly around the thigh to avoid a tourniquet effect a Velcro fastener engages the band at any point for a custom fit

DAL316 Leg Strap