### Wound Care

## TAPES

### **3M Health Care**

### **3M Blenderm Surgical Tapes**

Occlusive, transparent plastic tape transparent to allow monitoring waterproof flexible for hard-to-tape areas hypoallergenic 5-yard rolls 10 boxes/case

MMM1525-0 <sup>1</sup>/2", 24 rolls/bx MMM1525-1 1", 12 rolls/bx MMM1525-2 2", 6rolls/bx

#### **3M Durapore Surgical Tapes**

The "Silk-like" cloth tape with strong adhesion for securing dressings or devices 10-yard rolls 10 boxes/cs

MMM1538-0 1/2", 24rls/bx MMM1538-1 1", 12rls/bx MMM1538-2 2", 6rls/bx MMM1538-3 3", 4rls/bx

### **3M Medipore Soft Cloth Surgical Tapes**

The soft, stretchy, conformable tape products those are easy to use and easy on skin 10yard rolls

MMM2961 1", 24rls/bx MMM2962 2", 12rls/bx MMM2963 3", 12rls/bx MMM2964 4", 12rls/bx MMM2966 6", 12rls/bx

#### 3M Microfoam Surgical Tapes and 3M Microfoam Surgical Tape Patch

The conformable elastic foam tape or patch for compression applications, or securing dressings on challenging areas **3M Microfoam Surgical** Tape, 5-1/2 yd rolls (stretched) 6bx/cs MMM1528-1 1", 12rls/bx MMM1528-2 2", 6rls/bx MMM1528-3 3", 4rls/bx MMM1528-4 4", 3rls/bx



### **3M Micropore Surgical Tapes**

The dependable choice for a gentle, general purpose paper tape, gentle to the skin, conformable very porous hypoallergenic costeffective for many applications suggested applications: general dressing, fragile skin, e.g. facial dressings, repeated tape dressings, securing ostomy appliances

MMM1530-0 <sup>1</sup>/<sub>2</sub>" x 10yd, 24rls/bx, 10bxs/cs MMM1530-1 1"x 10yd, 12rls/bx, 10bxs/cs MMM1530-2 2"x 10yd, 6rls/bx, 10bxs/cs MMM1530-3 3"x 10yd, 4rls/bx, 10bxs/cs MMM1533-1 Tan, 1"x 10yd, 12rls/bx, 10bxs/cs MMM1533-2 Tan, 2"x 10yd, 6rls/bx, 10bxs/cs MMM1533-0 Tan, <sup>1</sup>/2"x 10yd, 24rls/bx, 10bxs/cs MMM1535-1 Disp. Pack, 1"x 10yd, 12rls/bx, 10bxs/cs MMM1535-2 Disp. Pack, 2"x 10yds, 6rls/bx, 10bxs/cs MMM1535-0 Disp. Pack, 1/2"x 10yd, 24rls/bx, 10bxs/cs

#### **3M Micropore Plus Surgical Tape**

Breathable, gentle and comfortable the same great performance as Micropore tapes, but now with greater adhesion for more secure performance

MMM1532-0 <sup>1</sup>/2"x 10yd, 24rls/bx, 10bxs/cs MMM1532-1 1" x 10yd, 12rls/bx, 10bxs/cs MMM1532-2 2"x 10yd, 6rls/bx, 10bxs/cs MMM1532-3 3"x 10yd, 4rls/bx, 10bxs/cs MMM1536-0 Disp. Pack, 1/2"x 10yd, 24rls/bx, 10bxs/cs MMM1536-1 Disp. Pack, 1"x 10yd, 12rls/bx, 10bxs/cs MMM1536-2 Disp. Pack, 2"x 10yds, 6rls/bx, 10bxs/cs MMM1536-3 Disp. Pack, 3"x 10yd, 4rls/bx, 10bxs/cs

### BSN-JOBST Inc.



### Elastoplast Heavyweight Stretch Tape

Helps control and reduce post injury edema by providing excellent compression and support strong adhesive eliminates the need for clips or extra tape adhesive contains a dry natural rubber 100% cotton backing promotes comfort; highly twisted crepe yarn provides elasticity molds to body contours while stretching and allows flexibility of joints without loosening 1/each, 36/cs

Tan, 1"x 5yds JOB02599002 Tan, 2"x 5yds JOB02600002 Tan, 3"x 5yds JOB02601002 Tan, 4"x 5yds JOB02593002 White, 1"x 5yds JOB02594002 White, 2"x 5yds JOB02595002 White, 3'x 5yds JOB02596002 White, 4"x 5yds



### Leukotape P Extremely strong zinc oxide adhesive and high tensile strength makes it ideal for patellofemoral taping Leukotape P and Cover-Roll Stretch are materials specified for the patellofemoral (McConnell) taping technique for correcting patellar position ideal for managing chronic shoulder problems by realigning shoulder related structures, taping and retraining muscles adhesive contains a dry natural rubber 1/bx, 30/cs

JOB76168 Tape, Porous Leukotape, 1-1/2"x 15yds

#### JOB02598002

### Tyco

### Healthcare

### Kendall



### Tenderskin Hypoallergenic paper tape

Lightweight, hypoallergenic paper tape for general use strong, tangle free, easy to tear natural porosity allows skin to breathe ideal for use when frequent taping and untaping are necessary 10 yard rolls

KEN1914 1", 12/bx, 72/cs KEN2419 2", 6/bx, 36/cs



### Curasilk Hypoallergenic Cloth Tape

Ideal hypoallergenic tape for strapping or general dressing use constructed with rayon acetate for strength available in tenyard rolls

#### KEN7137

<sup>1</sup>/2", 24/bx, 12/cs KEN7138 1", 12/bx, 12/cs KEN7139 2", 6/bx, 12/cs



### CURITY Standard Porous Tape

Adhesive tape with linear porosity for maximum breathability designed to provide high tensile strength and aggressive adhesion pinked edges allow easy tearing provides general use economy white 10 yard rolls

KEN2531 1", 12/bx, 12/cs KEN3027 1-1/2', 12/bx, 12/cs KEN6613 2", 6/bx, 12/cs

### Johnson Johnson Medical



### DERMICEL Hypoallergenic Cloth Tapes

DERMICEL Hypoallergenic Cloth Tapes are made from white cellulose acetate cloth backing and latexfree, water-based adhesive commonly known as "silk" tape applications for use: Securement of surgical dressings, wound dressings and tubes JJ5143 <sup>1</sup>/2"x 10yds, 24rolls/bx, 12bx/cs JJ5144 1"x10yds, 12rolls/bx, 12bx/cs JJ5145 2"x10yds, 6rolls/bx, 12bx/cs JJ5146 3"x10yds, 4rolls/bx, 12bx/cs



### DERMIFORM Hypoallergenic Knitted Tapes

DERMIFORM Hypo-Allergenic Knitted Tapes are surgical tapes made with knitted cellulose acetate tricot material backing and water-based adhesive they provide the strength of cloth tape with the lightness of paper tape

JJ5180 <sup>1</sup>/2"x 10yds, 24rolls/bx, 12bx/cs JJ5181 1"x 10yds, 12rolls/bx, 12bx/cs JJ5182 2"x 10yds, 6rolls/bx, 12bx/cs



### ELASTIKON Elastic Tapes

ELASTIKON Elastic Tapes are made from a high twist, cotton elastic material with rubberbased adhesive use for dynamic/pressure dressings requiring elasticity (mastectomies), support strapping, fractures, muscle injuries, knees, elbows and shoulders

JJ5172 1"x 5yds, Stretched, 12rolls/bx, 12bx/cs JJ5174 2"x 5yds, Stretched, 6rolls/bx, 12bx/cs JJ5175 3"x 5yds, Stretched, 4rolls/bx, 12bx/cs JJ5177 4"x 5yds, Stretched, 6rolls/bx, 6bx/cs



Johnson & Johnson Waterproof Tapes Tapes are made from cotton cloth backing with vinyl and polyethylene coating and rubber-based adhesive for moisture resistance that is easy to write on applications for use: Windowed casts or bandage wraps

JJ5132 1"x10yds, 24rolls/bx, 12bx/cs



**ZONAS Porous Tapes** ZONAS Porous Tapes are made with cotton cloth backing and rubber-based

adhesive these strong tapes are also known as "athletic" tape or "sports" tape use for wounds requiring immobilization (sprains, strains, etc), support strapping, patient positioning and restraints

JJ5104 1"x 10yds, 12rolls/bx, 12bx/cs JJ5106 2"x 10yds, 6rolls/bx, 12bx/cs JJ5132 3"x 10yds, 4rolls/bx, 12bx/cs

## ADHESIVE BANDAGES

### **3M Health Care**



### 3M Nexcare Active Strips Flexible Foam Bandages

Flexible foam bandages that cushion and protect cuts, scrapes, blisters and other minor skin injuries available in tan and bright colors

MMM17020 <sup>3</sup>/<sub>4</sub>" x 3", 100/bx, 12bx/cs MMM17021 1" x 3", 100/bx, 12bx/cs MMM17022 7/8", 100/bx, 18bx/cs MMM17030 Brights, <sup>3</sup>/<sub>4</sub>" x3", 100/bx, 12bx/cs



### **3M Nexcare Opticlude Orthoptic Eye Patch** Designed for comfortable vision-occlusive therapy

MMM1537 Junior, 2-1/2" x 1-1/4", 20/bx MMM1539 Regular, 3-1/4" x 2-1/4", 20/bx

## Allegiance

### **Esmarch Bandages**

Our bandages are constructed of natural latex and designed to help force blood out of an extremity by compressing the superficial vessels disposable bandages are supplied sterile and individually wrapped 10 bandages per case

BGH23580-043 4"x 3yd



Honeycomb elastic Bandages With Self-Closure

Our premium quality bandages are knitted with a blend of synthetic and 100% cotton yarn increased absorbency bandages offer controlled, consistent compression (contains natural rubber latex) and feature a VELCRO self-closure at each end that adheres to the bandage for ease of application standard bandages are beige with color-coded ends to denote width also available in pastel and neon colors and in double length (550") bandages are individually packaged and are available sterile and Nonsterile in 2", 3", 4" and 6" widths intended for single patient use Non-Sterile, 210"(533.5cm) Strchd, Lngth 2" BGH23593-020 (5cm), Beige, 72/cs BGH23593-030 3" (7.5cm), Beige, 72/cs BGH23593-040 4" (10cm), Beige, 72/cs BGH23593-060 6" (15cm), Beige, 48/cs

## ANDOVER



**CoFlex Bandages** An economical bandage for use when lightweight compression is required adheres to itself, no clips or fasteners are needed removal is quick and easy, it can be away scissors without unwrapping it maintains its original size and does not shrink in width available in full color cases or an assorted color pack **Assorted Colors** AND3100CP 1" x 5yds, 30/cs AND3150CP 1-1/2" x 5yds, 30/cs AND3200CP 2" x 5yds, 36/cs AND3300CP 3" x 5yds, 24/cs AND3400CP 4" x 5yds, 18/cs Tan AND3100TN 1" x 5yds, 30/cs AND3200TN 2" x 5yds, 36/cs AND3300TN 3" x 5yds, 24/cs AND3400TN 4" x 5yds, 18/cs White AND3100WH 1" x 5yds, 30/cs AND3200WH 2" x 5yds, 36/cs AND3300WH 3" x5yds, 24/cs

### BSN-JOBST Inc.



### Cover-Roll Stretch Adhesive

Cross-elastic non-woven bandage excellent for areas requiring flexibility offers the convenience of single-sheet taping over dressings cut to size to secure virtually any dressing air and exudates permeable, hypoallergenic and radiotranslucent JOB45552 2" x10yds, 1/bx, 12/cs JOB45548 4"x 2yds, 1/bx, 12/cs



**Coverlet Fabric Dressings** 

3600 degree adhesive coverage around the island wound pad seals off the wound from dirt and contamination conforms to body curves and movement to stay on longer fabric flexibility allows for mild compression to the wound variety of shapes and sizes to fit almost every part of the body 12bx/cs

JOB00385 4 Wing, 3"x 3", 50/bx JOB01307 Finger, Lg, 50/bx JOB01306 Finger, Sm. 100/bx JOB01390 Knuckle, 1-1/2" x3", 100/bx JOB00303 Oval, 1-1/4", 100/bx JOB00330 Patch, 1-1/2" x 2", 100/bx JOB00340000 Patch, 2" x 3", 50/bx

JOB00301 Spot, 7/8", 100/bx JOB00231 Strip, 1"x 3', 100/bx JOB00230 Strip, <sup>3</sup>/4" x 3", 100/bx JOB00801 Toe Shield, 2-1/2" x 1-3/8", 100/bx



Coverlet Eye Occlusors Eye Occlusor, non-Sterile, 20/bx JOB46429 Junior Eye Occlusor, Non-Sterile, 20/bx JOB46430 Regular

### Becton Dickinson

### Ace Elastic Bandage

For preventive bandaging and for support with strains and sprains washable and autoclavable five yard lengths packed 10 per box 50 per case

BD7336	2'
BD7337	3'
BD7338	4'
BD7339	6'

# coNco

### Medi-Rip Self-Adherent Bandage

Medi-Rip is a selfadherent elastic support and compression bandage, woven of 99% cotton yarn, it is porous, lightweight and cool the cohesive treatment creates a bandage that sticks to itself but not to the patient's skin or hairno adhesive residue available in sterile and non-sterile, it can be torn easily to the exact length needed-no scissors required Medi-rip conforms readily to body contours for ease of application and maximum patient comfort it will not slip or bind there is better bandaging control with Medi-Rip Medi-Rip Tan 2" x 5yds, COC0520 12/b, 8bx/cs

COC0520 2" x 5yds, 12/b, 8bx/cs COC0530 3" x 5yds, 12/bx, 8bx/cs COC0540 4" x 5yds, 12/bx, 8bx/cs

# Shur-Band Self-Closure Elastic

Shur-band is a latex-free high quality reinforced elastic support and compression bandage with a self-closure "flap" it is easy to fasten the self-closure "flap" at any point on the bandage clips, tape or pins are not required to hold the bandage in place COC5952 2" x 5yds, 10rolls/bx

COC5953	3" x 5yds,
10rolls/bx	-
COC5954	4" x 5yds,
10rolls/bx	-
COC5956	6" x 5yds,
10rolls/bx	-

#### EZE-Band Elastic Bandage w/ self-Closure

EZE-Band is an excellent quality elastic bandage that provides firm consistent support controlled stretch allows for consistent dependable tension and strict quality control assures maximum uniformity **Non-Sterile** COC5912 2" x5yds, 12/bg, 6bg/csCOC5913 3" x 5yds, 12/bg, 6bg/csCOC5914 4" x 5-1/2yds, 12/bg, 6bg/cs COC5916 6" x 5-1/2yds, 12/bg, 4bg/cs

# Tyco

Healthcare

## Kendall



### CURITY Adhesive Bandages

Made of perforated vinyl film all bandages feature "ouchless" TELFA pad; promote healing by not sticking to wound Hypoallergenic adhesive reduces potential for skin irritation porous adhesive allows skin to breathe KEN6947 <sup>3</sup>⁄<sub>4</sub>" x3". plastic, 50/bx, 3600cs **KEN7028** 1"x 3", plastic, 50/bx, 3600cs <sup>3</sup>⁄<sub>4</sub>" x 3", KEN3065 sheer, 50/bx, 3600cs KEN3270 1" x 3", sheer, 50/bx, 3600cs KEN6403 7/8" sheer spots, 100/bx, 3600cs KEN6416 1-1/2" x1-1/2", patch, 50/bx, 1200/csKEN613000 Assorted novelty bandage, 30/bx, 720/cs



### **CURITY Flexible Adhesive Bandages** Made of flexible fabric that stretches and conforms with movement ideal for fingers, toes. elbows, and other bending places strong, hypoallergenic adhesive formula sticks even when wet "Ouchless" TELFA pad promotes healing by not sticking to wound individually wrapped in easy-to-open, pull-tab, sterile sleeves

KEN5024 <sup>3</sup>/<sub>4</sub>", Flexible 50/pkg, 1200/cs KEN6136 2" x 3-1/2", Flexible 50/pkg, 600/cs KEN6131 <sup>3</sup>/4", flexible, Neon 50/pkg, 1200/cs KEN6133 1", flexible, Neon 50/pkg, 1200/cs KEN6417 3/8" x 1-1/2", junior 100/pkg, 3600/cs



# FLEX-WRAP Cohesive Bandage

Maintains consistent compression with nonslip support applied properly, provides and excellent pressure bandage without cutting off circulation provides excellent porosity in a lightweight, comfortable bandage adhesive to itself, no clips or fasteners are needed removal is quick and easy, can be cut away with scissors without unwrapping priced less than the competition maintains original size and does not shrink

KEN4584 Tan, 4" x 5yds, Bulk, 18rolls/cs

Johnson Johnson Medical



### **BAND-AID Brand**

Adhesive Bandages

**BAND-AID Brand** bandages give you more than just an adhesive bandage, they give you the styles and shapes to cover every need and every BAND-AID Brand Bandage offers these advantages: sterile, longlasting adhesion for continuous protection, non-adherent, absorbent pad for effective wound care and rounded corners to help prevent roll-up Flexible Fabric, sterile JJ4434 Strips, <sup>3</sup>/<sub>4</sub>" x 3", 100/bx, 12bx/cs JJ4444 Strips, 1" x 3", 100/bx, 12bx/cs JJ4436 Fingertip, 100/bx, 12bx/cs JJ4438 Knuckle. 100/bx, 12bx/cs Sheer, Sterile JJ4634 Strips, <sup>3</sup>/<sub>4</sub>" x 3", 100/bx, 12bx/cs JJ4644 Strips, 1" x 3", 100/bx, 12bx/cs JJ4650 Strips, <sup>3</sup>/<sub>4</sub>" x 3", Tray Pack, 150/tray, 10 trays/cs JJ4651 Strips, 1" x 3", Tray Pack, 150/tray, 10trays/cs JJ4828 Junior Strips, 5/8" x 2-1/4", 120/bx, 12bx/cs JJ4930 Spots, 7/8", 100/bx, 36bx/cs JJ5716 Extra Large, 2" x4-1/2", 50/bx, 18bx/cs

#### Butterfly Closures, Sterile, Waterproof (contains latex) JJ4331 medium, 1-3/4" x

3/8", 100/bx, 12bx/cs



Decorated, Sterile JJ4753 Sesame Street, <sup>3</sup>/<sub>4</sub>" x 3", 100/bx, 12bx/cs JJ4755 Arthur, <sup>3</sup>/<sub>4</sub>" 100/bx, 12bx/cs

## NutraMax



**Adhesive Bandages** For quality wound protection offers full pad coverage with an absorbent, non-stick pad long lasting adhesive breathable and comfortable latex free. Available in many different types: Clear, Sheer, Plastic, Flexible fabric (DR's offices, nursing homes), Soft flexible fabric (hospital), Heavy duty flexible fabric (industrial), Butterfly and Action

NUT15215 Stat Strip, Flex, 1" x 3", 100/bx, 12bx/csNUT15210 Stat Strip, Flex, <sup>3</sup>/<sub>4</sub>"x 3", 100/bx, 12bx/cs NUT15205 Stat Strip, Sheer, 1" x 3", 100/bx, 12bx/csNUT15200 Stat Strip. Sheer, <sup>3</sup>/<sub>4</sub>" x 3", 100/bx, 12bx/cs



NUT1075413 Stat Strip Glitter, <sup>3</sup>/<sub>4</sub>" x 3", 100/bx, 12bx/cs NUT1076413 Stat Strip Neon, <sup>3</sup>/<sub>4</sub>" x 3", 100/bx, 12bx/cs

#### **Character Strips**

Classic Characters in unique assortments all bandages are 100% perforated, have a nonstick pad and are latex free

NUT1075737 Bugs Bunny & Daffy Duck, ¾" x 3", 100/bx, 12bx/cs NUT1073737 Bugs Bunny & Tasmanian Devil, ¾" x 3", 100/bx, 12bx/cs NUT1076737 Coyote & Roadrunner, ¾" x 3", 100/bx, 12bx/cs NUT15600 Circus Characters, <sup>3</sup>/<sub>4</sub>" x 3", 100/bx, 12bx/csNUT15601 Herbie the Dinosaur,  $\frac{3}{4}$ " x 3", 100/bx, 12bx/cs NUT1075538 Charlie Brown & Lucy, <sup>3</sup>/<sub>4</sub>" x 3", 100/bx, 12bx/csNUT1307538 Sport Snoopy, 7/8", 100/bx, 12bx/cs NUT1076538 Peanuts-Snoopy, <sup>3</sup>/<sub>4</sub>" x 3", 100/bx, 12bx/cs

### GAUZE

## **CONCO**



Flexicon Elastic Gauze Bandage

Flexicon is an elastic gauze bandage that holds dressings securely in place without slipping the unique filament structure and nap ensure that Flexicon clings to itself the uniform weave of cotton and polymer yarns create a gauze with stretchability that conforms even to the most difficult areas woven edges minimize lint and loose threads Flexicon is durable with great strength, yet it is soft and comfortable to wear

Non-Sterile Flexicon COC2230 3"x 4-1/2yds, 20/bg

### Cypress Medical Products



### Conforming stretch Gauze Roll

Self-adhering stretch gauze that easily molds to body contours open weave construction allows fluids passing from the primary dressings to be easily absorbed non-sterile 12/bg, 8bg/cs

CYP41-02 2"x 4.1yds CYP41-03 3"x 4.1yds



Gauze Sponge High quality long staple cotton yarn for high absorbency and purity 8ply peel back bandage sterile 2/pk, 25pk/tray, 24trays/bx

CYP42-48 4"x 4"

### Johnson Johnson Medical

### SOF-BAND Bulky Bandages SOF-BAND Bulky Bandages are 6-ply,

100% cotton gauze bandages that offer cushioning and support ideal for dressing burns, plastic surgery, or Orthopaedic wounds conformable and easy to apply

JJ6964 4"x 84", relaxed 9flexible pouch)



### **KLING Fluff Rolls**

KLING Fluff Rolls are effective as both a primary or secondary dressing, helping to provide absorbency and protection the rolls are made from textured 100% cotton gauze, which produces loft and bulk for cushioning

JJ6930 4.5" x 4.1yds Stretched; Peel Pouch –6-ply, 100rolls/cs



KLING Conforming Gauze Bandages KLING conforming Gauze bandages are 100% cotton bandages that conform, stretch and "cling" they are ideal for securing dressings, Ivs or splints and for providing support and cushioning they are easy to apply and require minimal securement **Bulk, Non-Sterile** JJ6901 1"x 131", 24/bg, 4bg/cs JJ6902 2"x 131', 12/bg, 8bg/cs JJ6903 3"x 131", 12/bg, 8bg/cs JJ6904 4"x 131", 12/bg, 8bg/cs JJ6906 6"x 131", 6/bg, 8bg/cs Sterile JJ6922 2"x 131", 12/bg, 8bg/cs JJ6923 3"x 131", 12/bg, 6bg/cs JJ6924 4"x 131", 12/bg, 6bg/cs JJ6926 6'x131", 12/bg, 4bg/cs



#### **SOF-KLING Conforming Bandages** SOF-KLING Conforming Bandages are rayon/polyester blend bandages that conform, stretch and "cling" these low-linting bandages are highly absorbent they are soft, comfortable and gentle on the skin Sterile/Patient-Ready, Relaxed JJ6991 1"x60", 12/bg, 12bgs/cs JJ6992 2"x 65", 12/bg, 8bgs/cs JJ6993 3"x75", 12/bg, 6bgs/cs

JJ6994 4"x75", 12/bg, 6bgs/cs JJ6996 6"x 85", 12/bg, 4bgs/cs **Bulk, Non-Sterile Relaxed** JJ6981 1"x60", 24/bg, 4bgs/cs JJ6982 2"x65", 12/bg, 8bgs/cs JJ6983 3"x75", 12/bg, 8bgs/cs JJ6984 4"x75", 12/bg, 8bgs/cs JJ6986 6"x85", 6/bg, 8bgs/cs

### STERI-PAD Gauze Pads

12-ply, sterile, 100% cotton pads are constructed of fine mesh, USP Type VII gauze

JJ8513 2"x2", 100/bx, 24bx/cs JJ8516 3"x3", 100/bx, 24bx/cs JJ8519 4"x4", 100/bx, 12bx/cs

> **Tyco** Healthcare

## Kendall



Conform Stretch Bandages One-ply cotton/polyester blend crocheted bandages provides softness, conformability, low lint high absorbency ideal for securing dressings, IV's and splints or for providing mild compression or support holds securely to any body contour; allows for movement and some soft tissue swelling stays in place with minimal taping **Nonsterile, Bulk** KEN2239 1"x4.1yd,

NEIN2239	1 X4.1yu,
96rll/cs	
KEN2242	2"x4.1yd,
96rll/cs	
KEN2244	3"x4.1yd,
96rll/cs	
KEN2247	4"x4.1yd,
96rll/cs	
KEN2249	6"x.5yd,
48rlls/cs	
Sterile in Soft Pouch	
	ouch
KEN2230	1"x4.1yd,
KEN2230	
KEN2230 96rll/cs	1"x4.1yd,
KEN2230 96rll/cs KEN2231	1"x4.1yd,
KEN2230 96rll/cs KEN2231 96rll/cs	1"x4.1yd, 2"x4.1yd,
KEN2230 96rll/cs KEN2231 96rll/cs KEN2232	1"x4.1yd, 2"x4.1yd,
KEN2230 96rll/cs KEN2231 96rll/cs KEN2232 96rll/cs	1"x4.1yd, 2"x4.1yd, 3"x4.1yd,
KEN2230 96rll/cs KEN2231 96rll/cs KEN2232 96rll/cs KEN2236	1"x4.1yd, 2"x4.1yd, 3"x4.1yd,
KEN2230 96rll/cs KEN2231 96rll/cs KEN2232 96rll/cs KEN2236 96rll/cs	1"x4.1yd, 2"x4.1yd, 3"x4.1yd, 4"x4.1yd,



### Vaseline Petrolatum Gauze

A sterile, non-adhering, occlusive dressing consisting of the mesh, absorbent gauze impregnated with approximately three times its weight of white petrolatum

KEN8884417601

1"x8", Cison Dressing, 200/cs KEN8884413605 3"x9" in Peelable Foil Pack, 200/cs Vaseline Non-Adhering Sterile Petrolatum Gauze

Reduces pain and trauma during dressing change by conforming to all body contours without sticking to wound a special blend of petrolatum on fine mesh gauze helps protect wound from environmental contamination packed 12 per box, 6 boxes per case Peelable foil pack

KEN8884411600 <sup>1</sup>/<sub>2</sub>"x 72" Tube KEN8884412600 1"x 36" Strips KEN8884414600 3"x 18" Strips KEN8884415600 3"x 36" Strips KEN8884416600 6 "x 36" Strips



### Petrolatum Gauze

3% Bismuth Tribromophenate in a special petrolatum blend on fine mesh gauze clings and conforms to all body contours on-adherent and sterile **Overwrap, Peelable Foil** KEN8884431302 1"x 8" Strips, 50/bx, 200/cs KEN8884431605 5"x9" Strips, 50/bx, 200/cs **Patch, Peelable Foil Packs** KEN8884433301 1"x8", 50/bx, 200/cs KEN8884433400 2"x2", 25/bx, 150/cs KEN8884433500 5"x4", 25/bx, 150/cs **Strip, Peelable Foil Packs** KEN8884433605 5"x9", 50/bx, 200/cs

## **UNNA BOOTS**

## BSN-JOBST Inc.



### Gelocast Original Unna's Boot

Impregnated with zinc oxide and calamine to soothe irritated skin nonraveling cotton gauze, uniformly coated, dries quickly providing a semirigid cast available in 3" and 4" widths 1/bx, 12bx/cs

JOB01052 3"x 10yds JOB01053 4"x 10yds

## ConvaTec



### UNNA-FLEX elastic Unna Boot

**UNNA-FLEX Unna Boot** is a stretchable, zinc oxide impregnated bandage that provides an ideal form of compression for venous leg ulcer management it delivers compression that can improve venous return and reduce and control edema it provides the benefit of an UNNA **BOOT** without being rigid and does not contain commonly used preservatives which can cause allergic reactions UNNA-FLEX is easy to wrap around the leg without gaps or pleats it stretches easily to conform closely to the patient's leg while remaining flexible and comfortable enough to encourage patient compliance and ambulation packaged individually

CVT650941 4"x10yds, 12/bx

### **Graham-Field**

**Unna's Boot Dressing** Zinc impregnated gauze bandage 10 yards per roll GFUNNA4 4", 1/dz GFUNNA-C3 3", 1/dz GFUNNA-C4 4", 1/dz

#### **Medicopaste Bandage**

Bandage dressing aids in the treatment for dislocations, sprains, eczema and as a protective covering for skin grafts of the lower extremities

GF1565-3 3"x10"yds, 1/dz GF1565-4 4"x10"yds, 1/dz GF1565C-3 3"x10"yds, 1/dz GF1565C-4 4"x10"yds, 1/dz

### Dressings

## **3M Health Care**



### Reston Self-Adhering Foam

Self-adhering foam designed to distribute body to help prevent skin breakdown and pressure ulcers, and to prevent skin damage from splints, casts and prostheses it can also be used as a padding for beds, chairs and operating room tables MMM1560M Medium Support Pad, 7-7/8" x 11-3/4" MMM1563L Light Support Roll, 4' x 196"



**3M Tegaderm Transparent Dressings** Sterile "breathable" barrier lets oxygen in and moisture vapor out to encourage healing and reduce skin maceration conformable, hypoallergenic, and versatile allows continual site observation without removing dressing Easy application frame Style w/ Label MMM9505W 2-3/8" x 2-3/4", 20/bx, 10bxs/cs **MMM1624W** 2-3/8"x 2-3/4", 100/bx, 4bxs/cs MMM9506W 4"x4-3/4", 10/bx, 10bxs/cs MMM1626W 4"x4-3/4", 50/bx, 4bxs/cs **Original Frame Style** MMM1626 4"x4-3/4", 50/bx, 4bxs/cs MMM1627 4"x10", 20/bx, 4bxs/csMMM1628 6"x8", 10/bx, 8bxs/csMMM1629 8"x12", 10/bx, 8bxs/cs First Aid Style

GFUNNA3 3", 1/dz

2-3/8" x2-MMM1620 3/4", 100/bx, 6bxs/cs MMM16214`"x5-1/2", 50/bx, 4bxs/cs**Tegaderm HP Transparent Dressing** (Higher MVTR) MMM9534HP 2-3/8"x 2-3/4", 100/bx, 4bxs/cs **MMM9536HP** 4"x4-3/4", 50/bx, 4bxs/cs MMM9545HP Oval, 1-7/8" x 2-3/8", 50/bx, 6bxs/cs MMM9546HP Oval, 3-1/2" x4", 50/bx, 4bxs/cs MMM9548HP Oval, 5-1/2"x 7", 10/bx, 8bxs/cs **Tegaderm + Pad Transparent Dressing** with Absorbent Pad Easy applications Frame Style MMM3582 2"x2-3/4", 50/bx, 4bxs/cs MMM3584 2-3/8" x 4", 50/bx, 4bxs/cs MMM3586 3-1/2" x 4", 25/bx, 4bxs/cs MMM3587 Oval, 3-1/2"x 4-1/8", 25/bx, 4bxs/cs



### Coban Self-adherent Wrap

Can be used to secure, protect, compress and support sticks to itself without adhesive; needs no pins or clips available sterile 5 yard rolls MMM1581 1", 30rls/bx MMM1582 2", 36rls/bx MMM1583 3", 24rls/bx MMM1584 4", 18rls/bx MMM1586 6", 12rls/bx **MMM1583S** 3", 24rls/bx MMM1584S 4", 18rls/bx MMM1586S 6", 12rls/bx MMM1581HP 1", Hot Pink, 30rl/bx MMM1582B 2". Blue 36rls/bx MMM1582HP 2", Hot Pink, 36rl/bx MMM1583B 3", blue, 24rls/bx MMM1583G 3", green, 24rls/bx **MMM1583R** 3", red, 24rls/bx **MMM1583W** 3", white, 24rls/bx MMM1584B 4", blue, 18rls/bx MMM1584G 4", green, 18rls/bx MMM1584R 4", red, 18rls/bx MMM1583A 3", rainbow pack, 12rls/bxMMM1583N

## Kimberly-Clark Corporation

3", neon pack, 12rls/bx

Kimberly-Clark Wound Care Products Kimberly-Clark offers a broad line of competitively priced wound care products for

your needs Self-Adherent Wrap, non-Sterile, Latex-free KC19002 3", 24rolls/cs Coverall 5'n 1', Multi-Layered wound Dressing, Latex-free KC77045 4"x4". 10/bx, 10bx/cs **Xeroform Petrolatum** Gauze KC77034 1"x 8", 12/bx Shur-Conform oil **Emulsion Dressing** KC77041 3"x3", 50/bx, 12bx/cs

### Kimberly-Clark Packing Strips

 Plain, 12 bottles/cs

 KC77051
 ¼"x5yds

 KC77052
 ½'x 5yds

 KC77053
 1"x 5yds

 Iodoform, 12bottles/cs
 KC77061
 ¼" x 5yds

 KC77062
 ¼" x 5yds

## DeRoyal



#### COVADERM Surgical Site Dressing

Covaderm adhesive wound dressings replace gauze, tape and scissors with an easy to apply sterile dressing low adherent, highly absorbent pad with a protective, air permeable adhesive backing ideal for surgical incisors, superficial lacerations and abrasions 25/bx, 8bx/cs

#### STA46-001

4"x 4"overall, 2 ½" x 2 ½" pad



**COVADERM Plus Composite Dressing** A multi-layered composite dressing, Covaderm Plus provides the ultimate wound protection, layer by layer first, a non-adherent wound contact layer prevents trauma to fragile tissue layer 2 is an absorbent pad which wicks drainage away from the skin Layer 3 is a breathable polyurethane barrier film that provides moisture vapor transmission and creates a barrier to external contaminants layer 4 is a

comfortable, adhesivecoated fabric tape

STA46-401 4"x4" overall, 2 <sup>1</sup>/<sub>2</sub>" x 4" pad, 50/cs STA46-405-1 4" x 4" VAD Vascular Access, 25/cs

### ConvaTec



### DuoDERM CGF Dressings

DuoDERM CGF (control gel formula) dressings consist of a highly flexible, polyurethane outer foam layer and an adhesive skin contact layer that contains a unique hydrocolloid composition DuoDERM CGF dressings are conformable and can be easily molded into place to enhance patient comfort and minimize leakage DuoDERM CGF dressings are indicated for use on lightly to moderately exuding acute and chronic wounds **DuoDERM CGF Sterile** Dressings CVT187660 4"x4", 5/bxCVT187661 6"x6", 5/bx



### **DuoDERM Extra Thin** CGF Sterile Dressings DuoDERM CGF Extra

Thin dressings are flexible and conformable to fit difficult anatomical locations, and offer excellent protection for sites at risk for skin breakdown they can be used as primary dressings and are highly effective secondary dressings for use over wound filters

### CVT187955 4"x4", 10/bx



DuoDERM CGF control Gel Formula Border Sterile Dressings, Squares

An adhesive hydrocolloid foam border provides extra security and a secure fit in hard-to-dress (eg. Sacral area) they can be used as a primary or secondary dressing

CVT187971 4"x4", plus ¾" Adhesive Border, 5/bx



### AQUACEL Hydrofiber Wound Dressing

A unique Hydrofiber technology wound dressing that is comprised of a nonwoven pad or ribbon dressing made of hydrocolloid fibers AOUACEL is nonadherent, conforms easily to various wound shapes, supports autolytic debridement, and can easily be removed in one piece with minimal damage to newly forming granulating tissues AQUACEL is a primary dressings that retains 3 to 4 times more exudates than gauze and twice the amount of a leading alginate **AQUACEL Hydrofiber** Wound Dressing CVT177902 4"x4". 10/bx

AQUACEL Hydrofiber wound Dressing Rope CVT177904 <sup>3</sup>/<sub>4</sub>"x 18", 5/bx



KALTOSTAT Wound Dressing

KALTOSTAT wound Dressing is a highly absorbent calciumsodium alginate wound dressing derived from natural seaweed it absorbs exudates and is designed to control minor bleeding it is indicated for use on moderately to heavily exuding wounds KALTOSTAT creates a most wound healing environment which is supportive to the healing process it maintains its integrity when fully saturated and can easily be removed in one piece without causing trauma to newly formed granulating tissue KALTOSTAT requires the use of a secondary dressings such as DuoDERM CGF, or CombiDERM to maintain contact with the wound bed

KALTOSTAT Dressing VCT168242 3"x4-3/4, 10/bx, 5bx/cs KALTOSTAT Rope CVT168117 2gm, 5ropes/bx



SAF-Gel Hydrating Dermal Wound Dressing SAF-Gel provides an

optimal moist woundhealing environment and aids in autolytic debridement it is easy to apply and is conveniently packaged in a multipleuse tube it provides a moist environment that hydrates dry or necrotic tissue and reduces wound dead space

CVT145730 3oz Tube, 12/cs



**Sterile Dressing** CombiDERM ACD (absorbent cover dressing) is a sterile dressing with an absorbent pad that combines super-absorbent technology with the proven benefits of hydrocolloids CombiDERM ACD is a primary or secondary dressing that is indicated for use on heavily exuding wounds and may be used alone or in combination with absorbent wound filter dressings, such as AQUACEL or KALTOSTAT, to manage wound exudates the hydrocolloid adhesive of CombiDERM ACD has a thin profile that is designed to minimize edge roll-up CombiDERM is indicated for use on heavily exuding wounds

CVT651031 4"x4', 10/bx



### Kendall



**TELFA "Ouchless**' **Non-Adherent dressings** Absorbent cotton pad won't disrupt healing tissue by sticking to wound each sterile dressings is individually wrapped in peel-open envelopes ideal as a primary dressing for lightly draining wounds bonded on both sides with perforated non-adherent film; can be cut to fit any shape without separating easy-to-use; an ideal take home dressing

Perforated carton for easy dispensing sterile 1's in Peel-Back Package

KEN1961 2"x3", 2400/csKEN2132 3"x4". 2400/cs KEN1050 3"x4". 900/cs KEN1169 3"x6", 750/cs KEN1238 3"x8". 600/cs Nonsterile bulk, Mylar film on one side only

## KEN2891 3"x8", 1000/cs

Nonsterile bulk, Mylar film on both sides, perforated carton for easy dispensing sterile 1's in Peel-Back Package KEN1961 2"x3", 2400/cs KEN2132 3"x4", 2400/cs



Adhesive Dressings Made of highly absorbent cotton fabric bonded on both sides with perforated non-adherent film won't disrupt healing tissue by sticking to wound adhesive strips keep dressing intact perforated carton for easy dispensing sterile 1's in Peel-back package

KEN6017 2"x3", 100/bx, 2400/cs KEN7643 3"x4", 100/bx, 1200/cs



**Curity Eye Pads** Constructed of long fiber, absorbent cotton filler covered with fine mesh gauze on both sides sterile, ready to use Oval, contoured shape one per pack in convenient dispensing carton

KEN2841 Oval, 1-5/8"x 2-5/8", 50/bx, 600/cs



**TENDERSORB WET PRUF Abdominal Pads** Consists of three key layers: soft outer nonwoven layer, fluff filler to absorb and disperse fluid and wetproof moisture barrier to retard fluid strike through unique wet-proof moisture barrier helps minimize exposure to blood-borne pathogens absorbs more fluid than ordinary ABD pads for few dressing changes all four edges sealed to prevent lint residue and leaking sterile 1's in Peelback package

KEN9190 5"x9", 16/tray, 256/cs KEN9192 7-1/2"x8", 20/tray, 200/cs KEN9194 8"x10", 18/tray, 144/cs

#### CURITY Non-Adhering Dressings

Oil emulsion blend impregnated into open mesh, knitted fabric dressing non-adhering to wound site allows free drainage of exudates away from wound ideal for donor sites, delicate suture lines, graft sites, draining wounds

KEN6112 3"x3", 600/cs

### CURITY Lakeside Cotton Roll

Soft and absorbent, no surgical bleached cotton, non-sterile

KEN6030 12-1/2" x 56', 1/bg, 25rolls/cs



**CURITY Packing Strips** All natural 100% cotton fine mesh gauze ideal for wet-to-dry packing plain or Iodoform, guaranteed sterility **Plain Packing Strips** KEN7631 <sup>1</sup>/<sub>4</sub>"x 5yds, 12/csKEN7632  $\frac{1}{2}$ "x 5yds, 12/cs**Iodoform Packing Strips** KEN7831 <sup>1</sup>/<sub>4</sub>"x 5yds, 12/csKEN7832  $\frac{1}{2}$ "x 5yds,

Johnson Johnson Medical

2/cs



**BAND-AID Brand Surgical Dressings BABD-AID** brand Surgical Dressings consist of a non-adherent facing, a highly absorbent filler and tape backing these all-in-one dressings are more absorbent than standard gauze dressings which can help reduce dressing changes applications for use: Surgical incisions, lacerations, abrasions and invasive lines – scant to medium draining wounds or surgical incisions **BABD-AID Brand Adhesive Surgical** Dressings 4"x6" JJ4311 (3"x4" pad), 25/tray, 8trays/cs 8"x 6" JJ4313 (8"x3" pad), 25/tray, 4

Enter Control of Contr

trays/cs

#### **BIOCLUSIVE Transparent Dressings BIOCLUSIVE**

Transparent Dressings are made of thin, transparent polyurethane film and are used on invasive access devices and wounds they are easy to apply and help minimize skin irritation and protect fragile tissue

JJ2461 2"x3", 100/bx, 4bx/cs JJ2463 4"x5", 50/bx, 4bx/cs



RELEASE Non-Adhering Dressings RELEASE Non-Adhering dressings are sterile dressings that are used on surgical wounds and help provide atraumatic removal during changes these dressings may be used with or without a secondary dressing depending on the need for additional cushioning or absorption

JJ2052 2"x3", 50/bx, 24bx/cs JJ2054 4"x5", 50/bx, 24bx/cs



SURGIPAD Combine Dressings SURGIPAD Combine Dressings are multilayered abdominal dressings that help provide protection and

padding for major surgical wounds, while maximizing absorbency, wicking and dispersal of wound exudates sterile

JJ2144 8"x7-1/2", 20/bx, 12bx/cs JJ2145 5"x9", 25/bx, 16bx/cs JJ2148 8"x10", 20/bx, 16bx/cs



### JOHNSON & JOHNSON Eye Pads

Johnson & Johnson Eye pads are sterile eye dressings made with a soft non-woven cover, enclosing absorbent cotton filler they are used to help provide protection to the eye and absorb secretions

JJ8771 1-5/8" x 2-5/8", Oval, 50/bx, 12bx/cs JJ8773 2-1/8" x 2-5/8", regular, 50/bx, 12bx/cs



### NU GAUZE Packing Strips

NU GAUZE Packing Strips are ready to use sterile gauze strips constructed from specially woven, ravelresistant, 100% cotton gauze applications for use: Packing wound space, dehisced amputated leg, wet-to-dry packing to aid in nasal packing and removing necrotic tissue **NU GAUZE Packing** Strips, Sterile, Plain JJ8750 <sup>1</sup>/<sub>4</sub>" 5 yds, 12jars/cs JJ8751  $\frac{1}{2}$ ' x 5yds, 12jars/cs JJ8752 1" x 5yds, 12jars/cs Iodoform JJ8755 <sup>1</sup>/<sub>4</sub>" x5yds, 12jars/cs JJ8756  $\frac{1}{2}$ " x 5yds, 12jars/cs JJ8757 1"x 5yds, 12jars/cs

#### **RED CROSS Brand** First Aid Rolled Cotton

Johnson & Johnson + First aid rolled cotton is excellent for cleansing wounds and applying medication this product may also be used for infant and beauty care where softness, sterility and absorbency are important JOHNSON & JOHNSON + First Aid Rolled Cotton is surgical grade cotton which meets all USP requirements for absorbency purity and fiber length

JJ6026 16oz, sterile, 12bx/cs

### NutraMax



**Non-Stick Pads** For cleaning and protecting wounds nonstick pads-speeds healing and won't stick to wound individually sterile wrapped soft and comfortable absorbent latex free

NUT7565033 2"x3", 10/bx NUT7575033 3"x 4", 10/bx

## **SPONGES**

## ConvaTec

### **OPTIPORE Sponge** CVT125199 25/cs

Johnson Johnson Medical



#### **SOF-WICK Sponges**

Wound dressing sponges made from a unique, nonwoven 100% rayon fabric they are used primarily as a post-operative dressing to absorb wound exudates SOF-WICK Drain Sponges and SOF-WICK IV Sponges are fenestrated, special-use sponges made from a non-woven 100% rayon fabric they are specifically designed to fit around most drains, tubes or catheters **SOF-WICK Dressings** sponges, sterile/Patient-Ready JJ2375 4"x4". 6-Ply-2's, 50/bx, 12bx/cs **SOF-WICK Drain** Sponges, Sterile/Patient-Ready JJ2391 4"x4", 6-Ply-2's, 50/bx, 12bx/cs **SOF-WICK IV Sponges, Sterile/Patient-Ready** 2"x2", 6-JJ2392 Ply-2's, 70/bx, 20bx/cs



### JOHNSON & JOHNSON Gauze Sponges

JOHNSON & JOHNSON Gauze Sponges are 100% cotton, absorbent sponges constructed of fine mesh, 100% cotton gauze, available in multiple sizes and plies

#### Sterile/Patient-ready -**RSP** Type VII Gauze JJ2315 4"x4", 8ply-2's, 100/bx, 12bx/cs JJ2317 4"x4", 12-Ply-2's, 50/bx, 24bx/cs JJ2318 2"x2", 8-Ply-2's, 100/bx, 30bx/cs JJ2321 4'x4". 8ply-10's Tray, 96trays/cs 4"x4", 12-JJ2322 Ply-10's Tray, 96 trays/cs 4"x4', 16-JJ2329 Ply-10's tray, 72 trays/cs



MIRASORB Sponges MIRASORB Sponges are cotton blend sponges that can be used for cleaning, prepping, scrubbing and wound dressing they combine attributes of woven and non-woven sponges sterile

JJ602317 4"x4", 4-Ply-2's, 50/bg, 24bg/cs JJ602319 3"x3", 4-Ply-2's, 80/bx, 30bx/cs JJ602326 4"x4", 4-Ply-10's Tray, 96trays/cs



NU GAUZE General-Use Sponges NU GAUZE General-Use Sponges are unique rayon and polyester, non-woven sponges that provide improved strength, absorbency, and scrubbing texture compared to traditional gauze it is a general-use sponge that has a lowerend cost use than gauze Sterile/Patient-Readv JJ2337 4'x4", 4-Ply-2's, 50/bg, 24bg/cs JJ2338 2"x2", 4-Ply-2's, 50/bg, 60bg/cs JJ2339 3"x3", 4-Ply-2's, 50/bg, 48bg/cs JJ2344 4"x4", 4-Ply-6's Tray, 96/cs **Bulk, Non-Sterile** 2"x2". 4-JJ27632 ply, 200/sleeve, 40 sleeves/cs JJ27633 3"x3". 4-Ply, 200/sleeve, 20 sleeve/cs 4"x4". 4-JJ27634 Ply, 200/sleeve, 10 sleeves/cs JJ27635 2"x2", 3-Ply, 200/sleeve, 25sleeves/cs JJ27636 4"x4", 3-Ply, 200/sleeve, 20sleeves/cs

### RAY-TEC X-Ray

### **Detectable Sponges**

RAY-TEC X-Ray detectable sponges are absorbent sponges constructed of fine mesh 100% cotton gauze, available in multiple sizes and plies attached to the sponge are two blue barium sulfate x-ray detectable monofilaments Sterile/patient-Ready JJ2515 4"x4", 16-Ply, 10/Tray, 72 trays/cs



**TOPPER Dressing Sponges TOPPER** dressing Sponges are multi-layer absorbent sponges covered with a fabric facing they are specifically designed for use on top of wounds or incisions Sterile/Patient-Ready JJ2435 4"x3"2's (nonwoven facing). 50/bx. 24bx/cs 4"x4"2's JJ2436 (nonwoven facing), 50/bx. 24bx/cs **Bulk, Non-Sterile** 4"x4"2's JJ6415 (nonwoven facing), 100/sleeve, 20 sleeves/cs

> **Tyco** Healthcare

## Kendall



**CURITY Gauze Sponges** Made from high-grade cotton of USP Type VII gauze economical dressing in range of sizes and plys well suited for variety of applications; wound dressing, general cleaning, minor prepping, wound packing and debriding Nonsterile sponges conveniently indexed in 10's for faster, easier handling Nonsterile KEN2146 2"x2", 8 ply, 200/bg, 5000/cs KEN2252 2"x2". 12 ply, 200/bg, 8000/cs 3"x3", 12 KEN2346 ply, 200/bg, 4000/cs 4"x4", 8 KEN2556 ply, 200/bg, 4000/cs KEN2634 4"x4", 12 ply, 200/bg, 2000/cs KEN2733 4"x4". 12 ply, 200/bg, 2000/cs Sterile 2's in Peel-Back Package KEN1806 2"x2", 8 ply, 100/bx, 3000/cs KEN1903 3x3", 12 ply, 50/bx, 2400/cs 4"x4", 8 KEN2187 ply, 50/bx, 1200/cs **KEN3033** 4"x4". 12 ply, 50/bx, 1200/cs **Sterile 5's in Soft pouch KEN3967** 4"x4", 12 ply, 50/bx, 800/cs **Sterile 10's in Soft Pouch KEN3971** 4"x4", 12 ply, 100/bx, 1600/cs **Sterile 10's in Plastic Tray** 4"x4", 8 KEN6318 ply, 1280/cs **KEN6939** 4"x4", 12 ply, 1280/cs **KEN7605** 4"x4", 16 ply, 1280/cs

### **Curity Cover Sponges**

Filled with layer of cellulose wadding covered with nonwoven fabric economical, allpurpose sponge with good absorbent capacity Nonsterile sponges conveniently indexed in 10's for faster, easier handling. **Sterile 2's in Peel-back Package** KEN2913 4"x4", 50/pkg, 1200/cs

# EXCILON Drain and I.V. Sponges

Superior absorbency and fast wicking action help reduce infections and decrease the frequency of dressing changes unique pre-cut T-slit conforms snugly around drains, catheters, chest tubes, IV's and tracheotomies no loose threads to fall into wound eliminates need for sterile scissors, conserves nursing time sterile **Sterile 2's in peel-Back** Package **KEN7086** drain Sponge 4"x4", 6 ply,

# SORB-IT Drain and IV Sponges

50/pkg, and 600/cs

Lower cost alternative to premium EXCILON Drain and IV Sponges pre-cut notch conforms well around drains, catheters, chest tubes, IV's and tracheotomies eliminates the need for sterile scissors, conserves nursing time virtually lint-free sterile **Sterile 2's in Peel-Back Package** KEN6242P Drain Sponge, 4"x4", 6 ply, 50/ray, 600/cs

KEN6241 IV Sponge, 2"x2", 6 Ply, 70/tray, 1400/cs



### VERSALON Nonwoven All-Purpose Sponges

Ideal for applying ointments, prepping, wiping needles, cleaning slides and more rayon/polyester blend construction costeffective, highly absorbent and virtually lint-free **Sterile 2's in Peel-Back** Package KEN8042 2"x2", 4 ply, 50/pkg, 3000/cs 3"x3", 4 KEN8043 ply, 50pkg, 2400/cs KEN8044 4'x4', 4 ply, 50/pkg, 1200/cs Sterile 10's in Soft Pouch 4'x4".4 KEN8047 Ply, 1600/cs Nonsterile. Bulk KEN9022 2"x2", 4 ply, 200/pkg, 8000/cs KEN9023 3"x3", 4 ply, 200/pkg, 4000/cs KEN9024 4"x4", 4 Ply, 200/pkg, 2000/cs

KEN9132 2"x2", 3 ply, 200/pkg, 5000/cs KEN9134 4"x4", 3 ply, 200/pkg, 4000/cs



### **KERLIX Super** Sponges

Prewashed, fluff-dried woven gauze, with crinkle-weave pattern and diamond-fold construction offers superior loft, bulk and high absorbency superior primary dressing designed to transfer exudates away from wound into secondary dressing absorbency, strength and virtual lintfree performance ideal for pre-op prepping and cleansing, debriding and packing **Nonsterile Bulk KEN1167** Medium, 100/pkg, 600/cs **Sterile 5's in Soft Pouch** Package KEN3085 Medium. 50/pkg, 600/cs **Sterile 2's in Soft Pouch** Package **KEN2585** Medium, 40/pkg, 480/cs

### **KERLIX Disposable** Laparotomy Sponges

Gauze is prewashed and fluff-dried prior to cutting and sewing for maximum size, absorbency and softness vacuum delinted to minimize lint and particulate matter dualcolored Radiopaque handle remains visible when saturated with blood packed in sturdy polystyrene tray for easy stacking, storing Made in U.S.A with 100% U.S. cotton sterile, ready for use Tray Package, Sterile 5's Ken4031 17"x26', 100/cs



**CURITY Disposable** Laparotomy Sponges Prewashed, packaged in paper, banded 5's fourply with all cut edges folded in contains radiopaque element for xray detectability available in soft pouch sterile package or sturdy polystyrene tray with peel back lid also available in non-sterile bulk put-up **Prewashed tray Sterile**, **Banded 5's** 18"x18', Ken6022 100/cs**Prewashed soft Pouch** Sterile. Banded 5's KEN6522 18"x18", 200/cs



### VISTEC X-Ray

**Detectable Sponges** USP Type VII gauze unique dual-colored element combines blue Radiopaque strand with white strand which is visible when sponge is saturated with blood element is latex free and runs through both sides of sponge to facilitate counting precise, double vacuum delinting process assures low lint packaged sterile in sturdy polystyrene tray with peel-back lid or Nonsterile bags of 100, either unbanded or banded in 10's Sterile, banded 10's-Trays KEN7317 4"x4"(16 ply), 1350/cs

and traumatic wounds, superficial, partialthickness, and fullthickness wounds, dry or necrotic wounds, lightly, moderately, and highly exuding wounds

CVT121253 100ml, 24bottles/cs CVT121222 20ml, 100bottles/cs

## Wound Cleansers

## ConvaTec

#### Shur-Clens Wound Cleanser

Shur-Cleans Wound Cleanser is designed to be a nonirritating cleanser with a low toxicity index, making it an ideal product for cleansing and/or debridement of wounds containing debris, foreign particles, or devitalized tissue Shur-Clens is indicated for use on: Intact skin or stage I pressure ulcers, sutured